Blood Pressure Arm Monitor – Plus

Fully Automatic
Comfort Inflation Technology
Hypertension Risk Indicator
Irregular Heartbeat Detector
60 x 2 users
Average Reading for last 3 Memories
Date/Time Function
Blood Pressure Measurement Reminder Function

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Contents:

Introduction	3
About Blood Pressure	3-5
Blood Pressure Standard	5
Blood Pressure Fluctuation	6
Features of the device	
Parts Identification.....	6
LCD Display & Symbols	7
Hypertension Risk Indicator.....	8
Irregular Heartbeat (IHB) Detector	8
Comfort Inflation Technology	9
How to Use:	
Installing Batteries.....	9
Applying the Cuff.....	10-12
Measurement Procedures.....	13-15
Recalling Results from Memory.....	16
Clearing Results from Memory	17
How to Adjust the Date and Time.....	17
Blood Pressure Measurement Reminder	18-19
Technical Data:	
Possible faults, causes and rectification.....	20
Cautionary Notes.....	21
Guarantee.....	22
Specifications	22
Blood Pressure Log	23-24
EMC guidance and manufacturer's declaration.....	25-27

Introduction

What is blood pressure?

Your heart acts as a pump to circulate blood around your body and help supply it with oxygen. Blood pressure is the force needed for the heart to push blood through the arteries. The highest pressure in this cycle is when the heart contracts, this is called the SYSTOLIC BLOOD PRESSURE. Between contractions, the heart relaxes and blood flows into it since it is at its lowest pressure, which is called DIASTOLIC BLOOD PRESSURE.

Both blood pressure readings, the systolic and diastolic, are necessary to enable a doctor to evaluate the status of a patient's blood pressure.

Many factors such as physical activity, anxiety, or simply the time of day, can influence your blood pressure.

Drinking caffeine (in tea or coffee) can also temporarily raise your blood pressure, as can the nicotine in cigarettes.

Blood pressure can also follow a daily pattern, varying from minute to minute and typically being at its lowest while we are asleep.

These variations are even more pronounced in patients with high blood pressure.

Blood pressure is measured in millimetres of mercury (mmHg) and measurements are written with the systolic pressure before the diastolic e.g. a blood pressure written as 120/80 is referred to as 120 over 80.

High Blood Pressure

The discovery that you have high blood pressure is more often than not a chance finding.

People with high blood pressure usually feel well and there may be no external signs or symptoms of the condition – unless the blood pressure has been high for some time or complications have occurred. Such complications may affect the heart, kidney, brain and other important organs. They may bring on ill health and affect your enjoyment of life. It is because of these complications, rather than the high blood pressure itself, that this condition is taken so seriously today.

It is important to remember that the higher your blood pressure the higher the risk of heart disease or stroke.

High blood pressure is sometimes referred to as raised blood pressure, elevated blood pressure or hypertension.

Some people with this condition may require long-term drug treatment, usually in the form of tablets. However, others can lower their blood pressure by other means, such as changing their diet and/or lifestyle.

What causes high blood pressure?

Blood pressure is only classed as high if it does not go down when you rest.

Permanently raised blood pressure can be caused by several factors and is linked to other medical conditions such as diabetes, kidney disease and heart disease. Your blood pressure is affected by your lifestyle and important contributing factors are:

- Being overweight
- Having a high cholesterol level
- Drinking too much alcohol
- Eating too much salt
- Not eating enough fruit and vegetables
- Not exercising enough

Smoking also increases the risks associated with high blood pressure.

Pulse Rate

This blood pressure monitor also measures your pulse rate. Your pulse reflects your heart rate and is measured in terms of the numbers of beats per minute. Pulse rate varies from minute to minute, and is affected by many things, including exercise, stress, anxiety, certain medicines and some foods. Pulse rate is normally between 60 and 100 beats per minute.

Why is it beneficial to measure your blood pressure at home?

Monitoring your blood pressure at home gives you the advantage of taking blood pressure measurements at fixed times of the day, in familiar surroundings, without external influences. As a variety of conditions affect blood pressure, a single measurement is not sufficient for an accurate diagnosis. Home monitoring enables readings to be taken over the course of weeks and helps to identify an ongoing trend.

IMPORTANT INFORMATION BEFORE USE

Blood pressure measurements should be interpreted by a doctor or trained health professional who is familiar with your medical history.

When measuring your own blood pressure there are some basic rules to follow. These are easy to learn and so applying them at home will quickly become a routine.

As much as possible, conditions before and during measurement should always be the same. Therefore, choose the same time of day to measure your blood pressure.


It is also important that your bladder is empty before taking a measurement, as a full bladder may increase the reading.

Although such cases are rare, for people having an extremely weak or irregular pulse, error readings may result which prevent proper measurement. If extreme differences are found in blood pressure readings over a period, consult your doctor or pharmacist.

Blood Pressure Standard

The United States National High Blood Pressure Education Program Coordinating Committee has developed a blood pressure standard, according to which areas of low and high-risk blood pressure are identified. This standard is a general guideline as individual's blood pressure varies among different people and different age groups.

This blood pressure classification is based on historical data, and may not be directly applicable to any particular patient.



Blood Pressure Fluctuation


Blood pressure changes over the month and even throughout the day. It is also influenced by season and temperature.

For reliable monitoring and reference of blood pressure, keeping long-term records is recommended. Please use the blood pressure log in the back of this manual for recording and tracking crucial information: date, time, blood pressure, pulse rate and body condition etc.

Parts Identification

LCD Display

Symbols

Hypertension Risk Indicator: This compares the blood pressure measurement against preset guidelines to give general guidance on the risk of hypertension.

Irregular Heartbeat (IHB) Detector: This symbol signifies that an irregular heartbeat was detected.

Memory Zone 1: Appears when the device is operated in memory zone 1.

Memory Zone 2: Appears when the device is operated in memory zone 2.

Pulse Symbol: This shows the pulse rate per minute.

Weak Battery Symbol: Appears when the batteries should be replaced.

Blood Pressure Measurement Reminder Symbol:
This appears on the LCD display once the reminder function has been activated.

Memory Average: Displays average of the last 3 readings.

Memory Symbol: Shows the memory reading number.

Hypertension Risk Indicator

The United States National High Blood Pressure Education Program Coordinating Committee has developed a blood pressure standard, which classifies blood pressure ranges into 4 risk stages:

Normal / Suspected Pre-Hypertension / Suspected Stage 1 Hypertension / Suspected Stage 2 Hypertension. This device is equipped with an indicator, which gives a guide to the assumed risk level after each blood pressure measurement. It is important that you consult with your doctor regularly. Your doctor will advise you of your personal blood pressure range as well as the point at which you will be considered at risk.

Suspected Stage 2 Hypertension

Suspected Pre-Hypertension


Suspected Stage 1 Hypertension

Normal (No Symbol)

Irregular Heartbeat (IHB) Detector

This device is equipped with an Irregular Heartbeat (IHB) Detector, which detects an irregular heartbeat. If you see this symbol frequently, you should discuss with your doctor.


Comfort Inflation Technology

This device uses the oscillometric method to detect your blood pressure. It will determine the appropriate inflation level based on pressure oscillations and will only inflate as high as necessary, providing greater comfort.

Installing batteries

1. The device requires 4 x AA (LR6) alkaline 1.5V batteries
2. When replacing batteries, ensure all 4 are replaced at the same time.
3. Ensure that the batteries are all fitted the right way around, observing the + and - signs on the battery compartment.
4. Remove the batteries if the device is not intended to be used for extended periods.


Replace the batteries when:



1. The low battery symbol appears on the display.
2. The ON/OFF/START key is pressed and nothing appears on the display.

Applying the Cuff


1. Unwrap the arm cuff, leaving the end of the cuff with the arrow printed on it through the D-ring.


2. Put your left arm through the cuff loop, with the air tube pointing down your arm towards your wrist (Fig.1). Turn your left palm upwards and place the edge of the arm cuff approximately 1.5cm to 2.5 cm above the inside of your elbow. (Fig.2). Tighten the cuff by pulling the end of the cuff.


3. Position the tube over the middle of the arm. Fasten the cuff with the fabric fastener. Allow space for 2 fingers to fit between the cuff and your arm. Position the artery mark (Φ) over the main artery (on the inside of your arm Fig.3).


Note: Locate the main artery by pressing 2 fingers approximately 2cm above the bend of your elbow on the inside of your left arm. Identify where the pulse can be felt the strongest. This is your main artery (Fig.4).


4. Plug in the air tube into the device (Fig.5).

5. Lay your arm on a table (palm upwards) so the cuff is at the same height as your heart. Make sure the tube is not bent or trapped (Fig.6).

6. This cuff is suitable for use if the arrow falls within the OK range, as shown in Fig.7. If the arrow falls outside this area, you will need to purchase a different sized cuff. Contact your local retailer.


Important Notes:

Helpful tips to help you to obtain a more accurate reading:


- Blood pressure changes with every heartbeat and is in constant fluctuation throughout the day.
- Blood pressure recording can be affected by the position of the user, his or her physiologic condition and other factors. For greatest accuracy, wait one hour after exercising, bathing, eating, drinking beverages with alcohol or caffeine, or smoking to measure blood pressure.
- Your physical state should also be the same for each measurement, in a relaxed, sitting position with your legs uncrossed – at least half an hour after your last cup of coffee or tea, your last cigarette, physical exertion or a stressful situation.
- Do not take a measurement if you are under stress or tension.
- During measurement, do not talk or move your arm or hand muscles.
- Take your blood pressure at normal body temperature. If you are feeling cold or hot, wait a while before taking a measurement.
- If the device is stored at very low temperature (near freezing), place it in a warm location for at least one hour before using.
- Wait about 5 minutes before taking the next measurement.

Measurement Procedures


1. Press the User-Switching Key to select Memory Zone 1 or Memory Zone 2. After a memory zone has been selected, press the ON/OFF/STOP Key to reset the device so it can start the measurement in the chosen memory zone.


2. Press the ON/OFF/STOP Key to turn the device on. All digits will light up, checking the display functions. The checking procedure will be completed in 2 seconds.


3. After all of the symbols have appeared, the LCD screen will show a flashing “0.” The device is ready to work and will automatically inflate the cuff slowly to start the measurement.

4. When the measurement is complete, the cuff will deflate entirely. Systolic pressure, diastolic pressure and pulse will be shown simultaneously on the LCD screen. The memory sequence number will flash and depending on the result, the Hypertension Risk Indicator symbol will appear on the display, along with the time and date. The measurement is then automatically stored in the memory.


Note: The cuff will automatically re-inflate to approximately 220 mmHg if the device detects that your body needs more pressure to measure your blood pressure.


1. The device automatically switches off approximately 1 minute after last key operation.
2. To interrupt the measurement, simply press the ON/OFF/STOP or Memory key; the cuff will deflate immediately.
3. During the measurement, do not talk or move your arm or hand muscles.

Recalling Results from the Memory

1. The device has two memory zones (1 and 2). Each zone can store up to 60 measurements.
2. To read the memory values from a selected memory zone, use the User-Switching Key to select a memory zone (1 or 2) from which you want to recall the values. Press the Memory Key. The first reading displayed is the average of the last 3 measurements stored in memory.


3. Continue to press the Memory key to view previously stored measurements. Every measurement comes with an assigned memory sequence number.

Note: The memory bank can store up to 60 readings per memory zone. When the number of readings exceeds 60, the oldest data will be replaced with the new record.


Clearing Results from the Memory

1. Press the User-Switching Key to select memory zone 1 or memory zone 2.
2. Continue to press and hold the Memory Key for approximately 5 seconds, the data in that memory zone will be erased automatically.

How to adjust the Date and Time

1. To adjust the date/time within the device, press the Mode-Switching key once. The display will show a flashing number, this is the month.
2. Change the month by pressing the Memory key. Each press will increase the number by one. Press the Enter key to confirm the entry and the screen will show a flashing number representing the day.
3. Change the date and time as described in Step 2 above, using the Memory key to change and the Enter key to confirm the entries.
4. A "0" will reappear, as the Blood Pressure Device is ready for measurement again.




Blood Pressure Measurement Reminder

This device is equipped with a reminder function. 3 short beeps sound at the selected time. Four sets of reminders can be programmed.

Press the Mode-Switching key twice to enter Blood Pressure Measurement Reminder Setting. The Enter key allows you to switch from program AL1 to AL4. Press the Memory key to select the hour and press the Enter key to confirm the entry, repeat this procedure to select the minutes. A symbol will appear on the LCD display once the reminder is set successfully. Press the ON/OFF/START key to turn the device off.

Note: This function can also be used as an alarm clock.

Note: Once you have set up the reminder function, it will only be activated when the device is in the OFF mode. The reminder function will not be activated when the device is working, including taking measurement, checking previous data and adjusting date and time.

To remove the reminder / alarm times that have been set, press the Mode-Switching key twice to enter the Blood Pressure Measurement Reminder Setting. Press the Memory key to select the hour, continue to press this key until you see the flashing line (-) and press the Enter key to confirm the entry. The minutes will now flash, press the Memory key until you see the two flashing lines (- -) press the Enter key to confirm the entry and the  symbol will disappear.

Possible faults, causes and rectification

If any abnormality should arise during use, please check the following points.

No display when the ON/OFF/START key is pressed: Have the batteries run down or have the polarities (+ -) been positioned correctly? Replace the batteries with new ones, ensuring they are in the correct position.

EE Measurement Error Make sure the Air tube connector is securely connected to the air socket on the device and the cuff is fitted correctly. Keep your arm steady and measure again. If the error keeps occurring, return the device to your local retailer.

E1 Air Circuit Abnormality: Make sure the Air tube connector is securely connected to the air socket on the side of the device and measure again. If the errors still occur, return the device to your local retailer.

E2 Pressure Exceeding 300 mmHg: Switch the device off and start again. If the error keeps occurring, return the device to your local retailer.

E3 Data Error: Remove the batteries, wait for 60 seconds, and then reload. If the error keeps occurring return the device to your local retailer.

Er Exceeding Measurement Range: Measure again quietly. If error keeps occurring, return the device to the local retailer.

Note: If the device still does not work, return it to your retailer. Under no circumstance should you dismantle and repair the device by yourself.

Cautionary Notes

1. The device contains high-precision components. Therefore, avoid extreme temperatures, humidity, and direct sunlight. Avoid dropping the main device and protect it from dust.
2. Clean the blood pressure device body and the cuff carefully with a slightly damp, soft cloth. Do not press hard. Do not wash the cuff or use chemical cleaner on it. Never use thinner, alcohol or petrol (gasoline) as cleaner.
3. Leaky batteries can damage the device. Remove the batteries when the device is not used for a long period of time.
4. The device should not be operated by children, to avoid hazardous situations.
5. If the device is stored near freezing temperatures, allow it to acclimatise at room temperature before use.
6. This device is not serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have problems, please contact the retailer.
7. For users diagnosed with common arrhythmia (atrial or ventricular premature beats or atrial fibrillation), diabetes, poor circulation of the blood, kidney problems, or those who have suffered from a stroke, or for unconscious users, the device may have difficulty in determining the correct blood pressure measurement
8. To stop a measurement at any time, press the ON/OFF/START key, and the air in the cuff will deflate immediately.
9. Once the inflation reaches 300 mmHg, the device will start deflating rapidly for safety reasons.
10. Please note that this is a home healthcare product only and it is not intended to serve as a substitute for the advice of a physician or medical professional.
11. Do not use this device for diagnosis or treatment of any health problem or disease. Measurement results are for reference only. Consult a healthcare professional for interpretation of pressure measurements. Contact your physician if you have or suspect any medical problem. Do not change your medication without the advice of your physician or healthcare professional.
12. Electromagnetic interference: The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These may lead to temporary impairment of measurement accuracy.
13. Dispose of the device, batteries, components and accessories according to the local regulations. Batteries should be disposed of at a battery collection facility to minimise hazards to health and the environment and to enable materials to be recycled. This symbol is marked on the batteries and the packaging to remind you 
14. This device may not meet its performance specification if stored or used outside temperature and humidity ranges specified in Specifications.

Guarantee:

In addition to your statutory rights, Alvita agrees that if any defect in materials or workmanship appears in the product within two years after the original date of consumer purchase, it will repair, or at its option, replace the product free of charge. This applies only if the product has been used for domestic purposes and has not been damaged through misuse, accident or neglect and has not been modified or repaired by anyone other than Alvita or its authorised agents.

If a defect appears, please check that the article is being used in accordance with the instructions. If so, please return it with this warranty and proof of the date of purchase to the store you purchased it from, as soon as possible.

Specifications

Measurement Method :	Oscillometric
Measurement Range :	Pressure: 40~250mmHg; Pulse: 40~199 beats/minute
Pressure Sensor :	Semi conductor
Accuracy :	Pressure: ± 3 mmHg; Pulse: $\pm 5\%$ of reading
Inflation :	Pump Driven
Deflation :	Automatic Pressure Release Valve
Memory capacity :	60 memories for each zone x 2 zones
Auto-shut-off :	1 minute after last key operation
Operation Environment :	10°C~40°C (50°F~104°F); 40%~85% RH
Storage Environment :	-10°C~60°C (14°F~140°F); 10%~90% RH
DC Power Source :	DC 6V four R06 (AA) Batteries
Dimensions :	156 (L) X 110 (W) X 71 (H) mm
Weight :	470g (G.W.) (w/o Batteries)
Arm circumference :	24~36 cm (9.4"~14.2")
Limited users :	Adult users
 :	Type BF Device and cuff are designed to provide special protection against electrical shocks.

*Specifications are subject to change without notice.

Blood Pressure Log

Store the reading in memory zone:

Name:

Age:

Weight:

Date																				
Time																				
mmHg																				
	220																			
	200																			
	180																			
	160																			
	140																			
	120																			
Pulse																				
	100																			
	80																			
	60																			
Body Condition																				

Name:

1

Age:

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Weight:

[illegible]

EMC guidance and manufacturer's declaration

Guidance and manufacturer's declaration – electromagnetic emissions		
The MV801f is intended for use in the electromagnetic environment specified below. The customer or the user of the MV801f should ensure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The MV801f uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The MV801f is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Compliance	

EMC guidance and manufacturer's declaration

Guidance and manufacturer's declaration – electromagnetic immunity. The MV801f is intended for use in the electromagnetic environment specified below. The customer or the user of the MV801f should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2kV for power supply lines ± 1kV for input/output lines	± 2kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1kV line(s) to line(s) ± 2kV line(s) to earth	± 1kV differential mode Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT(>95% dip in UT) for 0,5 cycle 40% UT(60% dip in UT) for 5 cycles 70% UT(30% dip in UT) for 25 cycles .5% UT(95% dip in UT) for 5 s	<5% UT(>95% dip in UT) for 0,5 cycle 40% UT(60% dip in UT) for 5 cycles 70% UT(30% dip in UT) for 25 cycles .5% UT(95% dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the MV801f requires continued operation during power mains interruptions, it is recommended that the MV801f be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

EMC guidance and manufacturer's declaration

Guidance and manufacturer's declaration – electromagnetic immunity			
The MV801f is intended for use in the electromagnetic environment specified below. The customer or the user of the MV801f should ensure that is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the MV801f including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: d= 1,2 vP d= 2.2 vP 80Hz to 800MHz d= 2.3 vP 800Hz to 2.5 GHz Where P is the maximum output power rating of the transmitter manufacturer and (d) is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <i>a</i> should be less than the compliance level in each frequency range. B Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2,5 GHz	3 V/m	
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a: Fields strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MV801f is used exceeds the applicable RF compliance level above, the MV801f should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the MV801f. b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

EMC guidance and manufacturer's declaration

Recommended separation distance between portable and mobile RF communications equipment and the MV801f			
The MV801f is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the MV801f can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MV801f as recommended below, according to the maximum output power of the communication equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmission rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Preliminary Remarks

This Blood Pressure Monitor complies with the European regulations and bears the CE mark "CE 0473". The quality of the device has been verified and conforms to the provisions of the EC Council directive 93/42/EEC (Medical Device Directive), Annex 1 essential requirements and applied harmonized standards.

EN 1060-1: 1995/A1: 2002 Non-invasive sphygmomanometers

- Part 1 – General requirements

EN 1063-3: 1997/A2: 2009 Non-invasive sphygmomanometers

- Part 2 – Supplementary requirements for electro-mechanical blood pressure measuring systems.

EN 1064-4: 2004 Non-invasive sphygmomanometers

- Part 4 – Test Procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.

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