TM/  
MC

## BM 66/U

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**READ THIS MANUAL COMPLETELY AND CAREFULLY  
BEFORE USING THIS PRODUCT**

Keep this manual in a safe location for future reference

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## Contents

1. Getting to know your instrument.....	3
2. Important information.....	4
3. Package Contents.....	5
4. Device description.....	6
5. Getting started.....	7
5.1 Inserting/Replacing batteries.....	8
5.2 Cuff information.....	8
5.3 Setting the hour format, date and time.....	8
5.4 Operation with AC Adapter.....	9
6. Measuring blood pressure.....	10
6.1 Positioning cuff.....	12
6.2 Correct posture.....	13
6.3 General Notes.....	14
6.4 Taking measurement.....	15
7. Evaluating results.....	16
7.1 Irregular heartbeat symbol.....	18
7.2 World Health Organization (WHO) classification.....	18
8. Displaying and deleting measurements.....	20
9. Care, maintenance and storing the instrument.....	22
10. Error messages/trouble-shooting.....	23
11. Technical Specifications.....	24
12. AC Adapter.....	26
13. Glossary.....	26
14. Warranty and Return Policy.....	27

## IMPORTANT SAFETY NOTES

### Signs and symbols

The following signs appear in the Safety Section (page 2,3) and in this manual on pages 7, 8, 10, 11, 12, 16, 17, 18, 22, 23.

**READ THIS ENTIRE MANUAL, THE SAFETY SECTION AND ALL INSTRUCTIONS AND WARNINGS COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT. FOLLOW ALL SAFETY INSTRUCTIONS AND WARNINGS TO AVOID HAZARDOUS SITUATIONS AND TO MAKE CORRECT USE OF THIS PRODUCT.**

**WARNING** WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.

**CAUTION** CAUTION indicates a hazardous situation which, if not avoided, may result in minor or moderate injury.

### NOTICE

NOTICE addresses practices not related to personal injury, such as product and/or property damage.

## WARNING:

- Please do not share the cuff with other people, to avoid cross infection.
- This product is not a toy. Keep it out of reach of children, toddlers, and infants.
- Keep the product out of the reach of pets.
- Packaging materials are a deadly hazard for children and can cause suffocation. Remove all packaging materials immediately and keep them away from children at all times.
- This product contains small parts that may present a choking hazard to children. Keep the unit and all parts out of reach of children. **NEVER LEAVE CHILDREN OR THOSE WHO REQUIRE CLOSE SUPERVISION UNATTENDED WITH THIS DEVICE.**
- The tubing presents a strangulation hazard. Keep this product away from children and those who require close supervision, e.g. people with mental disorders. **NEVER LEAVE CHILDREN OR THOSE WHO REQUIRE CLOSE SUPERVISION UNATTENDED WITH THIS DEVICE.**

## SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE

If the instruction manual is damaged or if you no longer have the instruction manual in your possession, please contact customer service at 1-800-536-0366 or at [info@beurer.com](mailto:info@beurer.com)

### Dear Customer,

Congratulations on your purchase of this Life Brand Blood Pressure Monitor.

Please read and understand these instructions completely and carefully before using this product and keep them for later reference. Be sure to make them accessible to other users and observe the information they contain.

## 1. Getting to know your instrument

The upper arm blood pressure monitor is used for non-invasive measurement and monitoring of adults' arterial blood pressure. This allows you to quickly and easily measure your blood pressure and to display the last recorded measurement. The values determined are classified and graphically evaluated according to World Health Organization's (WHO) guidelines.

Blood pressure measurements determined by this monitor are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method, within

the limits prescribed by the American National Standard Institute, Electronic or automated sphygmomanometers.

It is recommended the performance should be checked every 2 years or after repair. Please contact customer service.


## 2. Important information

### Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and accessories:

	Caution
	Note Note on important information
	Follow instructions for use
	Type BF applied parts The cuff is type BF applied part
	Direct current
	Manufacturer

<b>Storage</b> 	Permissible storage temperature and humidity
<b>Operating</b> 	Permissible operating temperature and humidity
	Keep dry
SN	Serial number
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products.

### Advice on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.
- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes after eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.

- The measurements taken by you are for your information only – they are not a substitute for a medical examination! Discuss the measurements with your doctor, and never base any medical decisions on them (e.g. medicines and their administration)!
- Do not use the blood pressure monitor on newborns, pregnant women or patients with preeclampsia.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the device on people who have the specified upper arm measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions, remove the cuff from the arm.
- Avoid any mechanical restriction, compression or bending of the cuff line.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Ensure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g.


intravascular access or therapy, or an arteriovenous (AV) shunt.

- Do not use the cuff on people who have undergone a mastectomy.
- Do not place the cuff over wounds as this may cause further injury.
- You can either use the blood pressure monitor with batteries or with an electricity wall outlet. Please note that data transfer and data storage is only possible when your blood pressure monitor is supplied with power. As soon as the batteries are empty or the electricity cable is disconnected from the power supply, the blood pressure monitor loses the date and time.
- To conserve the batteries, the monitor switches off automatically if no buttons are pressed for one minute.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.

### 3. Package Contents

- Blood Pressure Monitor Model BM 66/U
- Arm cuff (arm circumference 22–42 cm / 8.7" to 16.5")
- Instruction manual
- Log Book
- Storage bag
- 4 x 1.5 V AAA alkaline batteries
- AC Adapter


## 4. Device description

1. Cuff
2. Cuff line
3. Cuff connector
4. Display
5. Memory buttons **M1/M2**
6. Start/stop button **ⓘ**
7. Resting indicator display
8. WHO scale
9. Connection for cuff connector (left-hand side)
10. Connection for AC Adapter
11. Battery compartment (on back side)

## Information on the display:

1. Date/time
2. Systolic pressure
3. Diastolic pressure
4. Pulse value
5. Pulse symbol
6. Release air (arrow)
7. Number of memory space/memory display average value (F), morning (P<sup>M</sup>), evening (P<sup>M</sup>)
8. Battery replacement symbol
9. User memory
10. WHO classification
11. Cardiac arrhythmia symbol

For product information, questions or customer service call toll free 1-800-536-0366 or email: [info@beurer.com](mailto:info@beurer.com)

## 5. Getting started

### WARNING:

- This product does not and is not intended to provide a medical diagnosis. Measurements results are for reference only. Self-diagnosis and treatment, e.g. regarding medication, using measured results represent a risk for your health. Always consult with a licensed physician for determination of appropriate medication and dosage thereof. Follow the instructions of your physician or licensed healthcare provider. If you have or suspect that you have a medical problem, promptly consult your physician. If you have an emergency please call 911 immediately.
- Only a physician or a trained health care professional who is familiar with your medical history is able to accurately interpret your blood pressure measurements. Consult your physician before starting blood pressure monitoring.

### Battery Handling Safety Precautions

- Use only the size and type of batteries specified.
- Be sure to follow the correct polarity when installing the batteries. Reversed batteries may cause damage to the device.


- Do not mix different types of batteries together (e.g. Alkaline and Carbon-zinc or rechargeable batteries) or old batteries with fresh ones. Always replace batteries as a simultaneous set.
- If the batteries in the device are depleted or the device will not be used for a long period of time, remove the batteries to prevent damage or injury from possible battery leakage.
- Do not try to recharge batteries not intended to be recharged; they can overheat and rupture (follow battery manufacturer's directions.)
- Do not dispose of batteries in fire, batteries may explode or leak.
- Clean the battery contacts and also those of the device prior to battery installation.
- Remove discharged batteries from the product and dispose/recycle in compliance with all applicable laws.
- Keep batteries away from children and pets. Batteries may be harmful if swallowed. Should a child or pet swallow a battery, seek medical assistance immediately.

### NOTE:

To conserve the batteries, the monitor switches off automatically if you do not press any buttons for 1 minute.

## 5.1 Inserting/Replacing batteries

- Remove the battery cover from the back of the monitor.
- Clean the battery contacts and those of the device prior to installation with a soft dry cloth.

- Insert or replace with only four AAA 1.5 V alkaline batteries. Make absolutely sure that you insert them with the correct polarity (+/-) as indicated (see the diagram inside the battery compartment). Do not use other types of batteries.
- Replace the battery cover carefully.

If the battery change icon  is continuously illuminated, measurement is no longer possible and you must replace all the batteries.


## 5.2 Cuff information

- Proper cuff size is important for accurate measurements. Only use the device on adults who have the correct upper arm circumference for this instrument.

This unit includes 1 size of arm cuff:

- Standard size fits 22–42 cm / 8.7" to 16.5"

## 5.3 Setting the hour format, date and time

This menu allows you to set the following functions, one after another.

It is essential to set the date and time. Otherwise, you will not be able to save your measured values correctly with a date and time and access them again later.

-  If you press and hold the **M1** or **M2** memory button, you can set the values more quickly.

- Press and hold the Start/stop button ① for 5 seconds.
- Choose the desired hour format with the **M1/M2** memory buttons and confirm with the Start/stop button ①.

The year flashes on the display.

- Choose the desired year with the **M1/M2** memory buttons and confirm with the Start/stop button ①.

The month flashes on the display.

- Choose the desired month with the **M1/M2** memory buttons and confirm with the Start/stop button ①.

The day flashes on the display.

- Choose the desired day with the **M1/M2** memory buttons and confirm with the Start/stop button ①.

① If the hour format is set as 12h, the day/month display sequence is reversed.

The hour flashes on the display.

- Choose the desired hour with the **M1/M2** memory buttons and confirm with the Start/stop button ①.

The minute flashes on the display.

- Choose the desired minute with the **M1/M2** memory buttons and confirm with the Start/stop button ①.

#### 5.4 Operation with AC Adapter

You can also operate this device with an AC Adapter. When doing so, there must not be any batteries in the battery compartment.

- To prevent possible damage to the device, the blood pressure monitor must only be used with the AC Adapter included with the unit.
- Insert the connector into the device. The AC Adapter must only be connected to the power outlet with the voltage specified on the rating plate of the adapter.
- Then insert the AC Adapter in the power outlet.

- After using the blood pressure monitor, unplug the AC Adapter from the power outlet first and then disconnect it from the blood pressure monitor. As soon as you unplug the AC Adapter, the blood pressure monitor loses the date and time setting but the saved measurements are retained.

Detailed instructions concerning the use of the AC Adapter can be obtained from the manual which is included in the AC Adapter box. Please read and understand all instructions and warnings completely and carefully before using the AC Adapter.

Make sure to take the batteries out of the battery compartment before using the adapter.

## 6. Measuring blood pressure

### **WARNING:**

- People with disabilities, activity limitations, or who are physically frail, should be assisted by another person when using this unit.
- Consult your physician or pharmacist before starting to measure your blood pressure.
- Read and understand all instructions and warnings before using this device.

- Like any oscillometric blood pressure measurement devices, certain medical conditions can affect the measurement accuracy, among others:
  - disorder of the cardiac rhythm
  - very low blood pressure
  - patients in shock
  - patients with very low blood perfusion
  - diabetes
  - blood vessel anomalies
  - people with electrical implants, such as a cardiac pacemaker
  - women who are pregnant.

Due to their condition the oscillometric measurement method can produce incorrect readings. This represents a risk for your health, since values may be interpreted incorrectly. Always consult your physician to determine what will be suitable for you.

- This product is not intended to serve as a substitute for the advice of a physician or medical professional. This product is not intended to substitute for regular medical checkups. Contact your physician for specific information about your blood pressure.
- Women who underwent a breast or axillar lymph node removal operation should consult a physician or licensed healthcare provider before starting blood pressure measurements.

- Prior to use, you have to ensure that the blood pressure monitor is free of damage. If in doubt, do not use the unit and contact customer service at 1-800-536-0366.
- This device is intended only for personal, non-commercial, adult use in measuring blood pressure and pulse rate. Do not use the product for any other purpose. This product is not intended for use in a hospital, physician's office, or any other health care facility.
- The unit must be used in accordance to the specified ambient conditions, otherwise the accuracy of readings might be affected. See section "9. Care, maintenance and storing the instrument" and section "11. Technical Specifications" for details.
- Do not wrap the cuff around body parts other than your upper left arm. Misuse represents a risk to your health.
- Use this monitor only on humans.
- This product is not intended for use by or on children, toddlers and infants or on persons who cannot express their consent, e.g. persons with mental disorders or the like. Consult your physician for alternative methods of measuring a child's blood pressure.
- Do not drape tube around neck. This presents a strangulation hazard.
- People with disabilities, activity limitations, or who are physically frail, should be assisted by another person when using this unit.
- Proper cuff size is important for accurate measurements. Only use the instrument on adult persons whose

upper arm has the right upper arm circumference for the instrument. See section "11. Technical Specifications" for suitable arm circumferences.

- Do not use any cuffs, tubes, parts and accessories other than those explicitly recommended by the manufacturer for use with this product. Parts and accessories not approved for use with the device may cause damage to your health and to the product.

### **WARNING:**

- **Electromagnetic interference:** Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens) while it is in operation, as inaccurate measurements may result. To prevent such interference, use the unit at a sufficient distance from such devices or turn the devices off.
- Do not measure your blood pressure while operating a vehicle or in any situation which requires your full attention.
- Any cuff related blood pressure measurement in high repetition rates, can lead to severe measurement side effects, e.g.
  - any nerve compression with temporary arm/hand paralysis
  - the release of an arterial or venous thrombus, which can cause a life threatening situation.

Please contact your physician about the specific risks of cuff pressure in your specific case.

## CAUTION:

- The device should not be used when your arm has been wounded/injured or when a catheter has been inserted. Such use may result in injury.
- Remove any kind of upper arm jewelry or the like before taking a measurement. This could cause bruises.
- Do not place the arm cuff over heavy clothing (e.g. a jacket or sweater sleeve) as the blood pressure monitor will not be able to take a proper measurement and there is an elevated danger of acquiring hematoma or skin marks during the course of the measurement.
- In case the cuff does not stop inflating, interrupt the measurement by pressing the start/stop button  and open the cuff at once.
- When applying the arm cuff, make sure there are no wrinkles in the cuff. This could cause bruises.
- Do not exert any kind of pressure on the hose during measurement, e.g. laying your arms or any other object on the hose. This could cause incorrect measurements.
- Blood pressure measurements can lead to temporary marks on the skin at the site of the cuff placement. This is especially the case in high repetition rates, in hypertonic patients and in patients with weak pulses. In rare cases a mark may persist for a couple of days. Please

contact your physician about these specific risks of cuff pressure in your specific case.


- For hygienic reasons the cuff is intended only for use by one person.
- It is recommended that the cuff be cleaned after every 200th use.
- It is recommended that the cuff be disinfected twice a week using a cloth that is lightly moistened with ethyl alcohol (75–90 %); then allowed to air dry.

### 6.1 Positioning cuff


#### NOTE:

Always properly position the Arm Cuff on the arm at heart level, according to the following instructions.


Fit the cuff round your bare left upper arm. Do not place the arm cuff over heavy clothing (e.g. a jacket or sweater sleeve) as the blood pressure monitor will not be able to take a proper measurement. Blood circulation in the arm should not be restricted by tight clothing or other objects.


The cuff should be placed on the upper arm so that the lower edge is 2 to 3 cm (0.8" to 1.2") above the bend of the elbow and above the artery. The tube should be in line with the centre of the palm.

Now place the free end of the cuff snugly, but not too tightly, around the arm, and fix it with the Velcro fastener.




Insert the cuff tubing into the socket for the cuff attachment on the left side of the unit. Make sure it is inserted securely in the main unit and that there are no kinks in the cuff tubing.

**Important:** This device should only be operated with the original cuff. The cuff included with the device is suitable for an upper arm circumference of 22–42 cm (8.7" to 16.5"). Replacement cuffs from other models cannot be used with this model. Doing so, will result in inaccurate readings or error messages.

## 6.2 Correct posture

- Rest at least for 5 minutes before each measurement. Otherwise there may be erroneous results.
- You can perform the measurement either sitting or lying down. Always make sure that the cuff is on a level with your heart.
- In order not to distort the result, it is important to keep still during the measurement and not talk. The most precise measurement results are obtained when taken in a relaxed state. You may close your eyes and breathe calmly during the measurement.
- Please note that any muscle movement during inflation or deflation can cause measurement error.

### 6.3 General Notes

- For reliable monitoring and reference of blood pressure, keeping long-term records is recommended.
- To minimize measurement variations due to physical activity, rest at least 5 minutes before measuring your blood pressure (or 15 minutes after strenuous activity). You should not be physically tired or exhausted while taking measurement.
- In order to obtain a resting condition blood pressure, avoid eating, drinking alcohol and caffeinated beverages, smoking, exercising, and bathing for at least 30 minutes before taking a measurement.
- Stress raises blood pressure.
- Perform measurements in a quiet and relaxed environment at room temperature.
- Always wait at least 5 minutes between measurements to allow the blood circulation in your arm to return to normal. You may need to increase the wait time depending on your individual physiological characteristics. Because of normal physiologic changes and the many external factors influencing blood pressure, it is unusual to obtain identical blood pressure measurements, even when taken a few minutes apart. Speak with your physician to determine what is affecting your blood pressure and which variation may be seen as normal in your case.
- The accuracy of any blood pressure measurement with this device can be affected by a multitude of causes.

Some can be avoided some have to be accepted at the time of measurement. Speak with your physician about possible causes affecting your blood pressure.

- Do not inflate the cuff unless wrapped around arm.
- If no button is pressed for approx. 1 minute the automatic shutoff function shuts off the unit in order to preserve the batteries.
- Do not touch the unit, the tube and cuff and/or press any buttons when measurement is in progress, except for stopping the measurement. Doing so may cause incorrect readings.

## 6.4 Taking measurement

### Performing the blood pressure measurement

As described above, attach the cuff and adopt the posture in which you want to perform the measurement.

- To start the blood pressure monitor, press the Start/stop button . All display elements are briefly displayed.

The blood pressure monitor will begin the measurement automatically after 3 seconds.

The cuff automatically inflates.

-  Measuring can be cancelled at any time by pressing the Start/stop button .


The cuff's air pressure is slowly released. If you already recognize having a tendency for high blood pressure, you should reinflate the cuff and increase the cuff's pressure again. As soon as a pulse is found, the pulse symbol  is displayed.

Measurement

- Systolic pressure, diastolic pressure and pulse measurements are displayed.
- A symbol at the bottom of the display also lights up to indicate whether you were sufficiently relaxed during the blood pressure measurement (green symbol = sufficiently at rest; red symbol = not at rest).

Observe the chapter on interpreting results/measuring the resting indicator in these instructions for use.

- $E_{r\_}$  appears if the measurement could not be performed properly. Observe the chapter on error messages/troubleshooting in these instructions for use and repeat the measurement.
- Now select the desired user memory by pressing the **M1** or **M2** memory buttons. If you do not select a user memory, the measurement is stored in the most recently used user memory. The relevant symbol  or  appears on the display.

Measurement

- Using the Start/stop button **①**, switch off the blood pressure monitor. The measurement is then stored in the selected user memory.

If you forget to turn off the device, it will switch off automatically after approx. 1 minute. In this case too, the value is stored in the selected or most recently used user memory.

- Wait at least 5 minutes before taking another measurement!

### NOTE:

This monitor will re-inflate automatically if the system detects that your body requires more pressure for measurement. As soon as a heart rate is detected, the heart rate symbol is displayed.

- When the measurement is completed the cuff will deflate entirely. Heart rate, systolic and diastolic blood pressure are displayed.
- You can interrupt measurement at any time by pressing the start/stop button **①**.
- $E_r$  appears if it has not been possible to perform the measurement properly. Observe the section on error messages/troubleshooting and repeat the measurement.
- The test result is saved automatically.

- To switch off and release the pressure, press the start/stop button **①** again. If you forget to switch off the device, it switches off automatically after approx. 1 minute.

Wait at least 5 minutes before taking another measurement.

## 7. Evaluating results

### WARNING:

- This product does not and is not intended to provide a medical diagnosis. Measurements results are for reference only. Self-diagnosis and treatment, e.g. regarding medication, using measured results represent a risk for your health. Always consult with a licensed physician for determination of appropriate medication and dosage thereof. Follow the instructions of your physician or licensed healthcare provider. If you have or suspect that you have a medical problem, promptly consult your physician. If you have an emergency please call 911 immediately.
- Only a physician or a trained health care professional who is familiar with your medical history is able to accurately interpret your blood pressure measurements. Consult your physician before starting blood pressure monitoring.
- Please note that technically related measuring tolerances are possible. Please see section “11. Technical Specifications” for details.

- If you notice abnormal or suspicious variations in blood pressure measurements, consult your physician immediately.

### **Cardiac arrhythmia:**

This instrument can identify possible cardiac arrhythmia disorders during measurement and if necessary indicates the measurement with the flashing icon .

This may be an indicator for arrhythmia. Arrhythmia is a condition where the heart rhythm is abnormal as a result of defects in the bioelectrical system controlling the heart beat. The symptoms (omitted or premature heart beats, slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress or lack of sleep. Arrhythmia can only be ascertained through examination by your doctor.

Repeat the measurement if the flashing icon  is displayed after the measurement. Please note that you should rest for 5 minutes between measurements and not talk or move during the measurement. If the icon  appears often, please contact your doctor. Any self-diagnosis and treatment based on the test results may be dangerous. It is vital to follow your doctor's instructions.

### **WHO classification:**

In accordance with the guidelines/definitions of the World Health Organization and the latest findings, the measurements can be classified and assessed according to the following table.

However, these standard values serve only as a general guideline, as the individual blood pressure varies in different people and different age groups etc.

<b>Range of blood pressure values</b>	<b>Systolic (in mmHg)</b>	<b>Diastolic (in mmHg)</b>	<b>Measure</b>
Grade 3: Severe hypertension	>= 180	>= 110	Seek medical advice
Grade 2: Moderate hypertension	160–179	100–109	Seek medical advice
Grade 1: Mild hypertension	140–159	90–99	Have it checked regularly by doctor
High-normal	130–139	85–89	Have it checked regularly by doctor
Normal	120–129	80–84	Check it yourself
Optimal	< 120	< 80	Check it yourself

Source: WHO, 1999

**NOTE:**

A single measurement does not provide an accurate indication of your true blood pressure. You need to take and record several readings over a period of time. Try to measure your blood pressure at the same time each day for consistency.

**7.1 Irregular heartbeat symbol****WARNING:**

- The 'irregular heartbeat' function does not replace a cardiac examination, but may help to detect potential pulse irregularities at an early stage. Always consult your physician to determine what will be suitable for you.
- The 'irregular heartbeat' function is not designed for diagnosing or treating an arrhythmic disorder. Arrhythmia can only be ascertained by a licensed physician.

This symbol (♥), indicates that certain heartbeat/pulse irregularities were detected during the measurement. In general, an irregular heartbeat rhythm is defined as a rhythm that is more than 25 % slower or 25 % faster from the average rhythm detected while the monitor is measuring the systolic and diastolic blood pressure.

**NOTE:**

- Talking, moving, shaking or an irregular pulse during the measurement can result in the appearance of this symbol.

- Repeat the measurement if the flashing icon is displayed after the measurement. Please note that you should rest at least 5 minutes between measurements and not talk or move during the measurement. If the icon continues to be displayed in subsequent readings, we recommend you to consult your physician in order to determine if it is cause for concern. Do not try to interpret readings or attempt to treat any condition yourself. Always follow the guidance of a health care professional. If you have an emergency call 911 immediately.

**7.2 World Health Organization (WHO) classification****WARNING:**

- The WHO chart is not intended to replace a medical diagnosis. This chart is only a reference for different classifications of blood pressure.

According to WHO Guidelines/Definitions and the latest findings, the test results of adults can be classified and evaluated according to the following chart:

**NOTE:**

This standard, however, is a general guideline as an individual's blood pressure varies among different people, age groups, etc. Please consult your physician for proper diagnosis.

The bar graph in the display and the scale on the unit indicate the range of the blood pressure which has been recorded.

#### **NOTE:**

If the values for systolic and diastolic pressure are in two different WHO ranges (e.g. systolic in the high-normal range and diastolic pressure in the normal range) the graphic WHO classification on the unit indicates the higher range (high-normal in the example described).

#### **Measuring the resting indicator (using the HSD diagnosis)**

The most frequent error made when measuring blood pressure is taking the measurement when not at rest (hemodynamic stability), which means that both the systolic and the diastolic blood pressures are distorted.

While measuring the blood pressure, the device automatically determines whether you are at rest or not.

If there is no indication that your circulatory system is not sufficiently at rest, the symbol  (hemodynamic stability) lights up green and the measurement can be recorded as a reliable resting blood pressure value.

#### **GREEN: Hemodynamically stable**

The systolic and diastolic pressure measurements are taken when the circulatory system is at rest and are therefore reliable indicators of a resting blood pressure.

However, if there is an indication that the circulatory system is not sufficiently at rest (hemodynamic instability), the symbol  lights up red.

In this case, the measurement should be repeated after a period of physical and mental rest. The blood pressure measurement must be taken when the patient is physically and mentally rested, as it will be the basis for diagnosing the blood pressure level and regulating the patient's medical treatment.

If the symbol  does not light up green or red, then it could not be determined whether the circulatory system was sufficiently at rest or not. In this case, the measurement should be repeated after a period of physical and mental rest.

#### **RED: Hemodynamically unstable**

It is very probable that the systolic and diastolic blood pressures have not been measured while the patient was at rest and the resting blood pressure measurement has therefore been distorted.

Repeat the measurement after a minimum period of five minutes' rest and relaxation. Go to a sufficiently quiet and

comfortable spot and remain there calmly; close your eyes, breathe deeply and evenly and try to relax.

If the next measurement also shows insufficient stability, you can repeat the measurement after another resting period. If the measurements continue to show some instability, identify these blood pressure measurements as having been taken when the circulatory system had not been sufficiently rested.

In this case, nervousness or inner anxiety may be the cause and this cannot be cured by brief periods of rest. Existing cardiac arrhythmias (irregular heart beat) may also prevent a stable blood pressure measurement from being taken.

A lack of resting blood pressure can have various causes, such as physical or mental strain or distraction, speaking or experiencing cardiac arrhythmias during the measurement. In an overwhelming number of cases, the HSD diagnosis will give a very good guide as to whether the circulatory system is rested when taking the measurement. Certain patients suffering from cardiac arrhythmia or chronic mental conditions can remain hemodynamically unstable in the long-term, something which persists even after repeated periods of rest. The accuracy of the resting blood pressure results is reduced in these users. Like any medical measurement method, the precision of the HSD diagnosis is limited and can lead to incorrect results in some cases. The blood pressure measurements taken when the circulatory system was at rest represent particularly reliable results.

## 8. Displaying and deleting measurements

The results of every successful measurement are stored together with the date and time. If there are more than 60 measurements, the oldest measurements are deleted.

- To access memory recall mode, the blood pressure monitor must first be started. To do this press the Start/stop button **1**.
- Within 3 seconds of the full-screen display appearing, select the desired user memory (**M1** **M2**) with the **M1** or **M2** memory button.
- To view the measurements for user memory **M1**, press the **M1** memory button.
- To view the measurements for user memory **M2**, press the **M2** memory button.

Your last measurement will appear on the display.

User memory

User memory

- Press the relevant memory button (**M1** or **M2**).

**i** If you have selected user memory1, the **M1** memory button must be pressed. If you have selected user memory2, the **M2** memory button must be pressed.

**A** flashes on the display.

The average value of all saved measured values in this user memory is displayed.

- Press the relevant memory button (**M1** or **M2**).

**M** flashes on the display.

The average value of the morning measurements for the last 7 days is displayed (morning: 5.00 a.m. – 9.00 a.m.).

- Press the relevant memory button (**M1** or **M2**).

**P** flashes on the display.

The average value of the evening measurements for the last 7 days is displayed (evening: 6.00 p.m. – 8.00 p.m.).

- When the relevant memory button (**M1** or **M2**) is pressed again, the last individual measurement is displayed (in this example, measurement 03).

- When the relevant memory button (**M1** or **M2**) is pressed again, you can view your individual measurements.
- To switch the device off again, press the Start/stop button **⏻**.

**i** You can exit the menu at any time by pressing the Start/stop button **⏻**.

- To clear the memory of the relevant user memory, you must first select a user memory.
- Start individual measurement access.
- Press and hold the **M1/M2** memory buttons for 5 seconds.

All the values in the current user memory are deleted.

## 9. Care, maintenance and storing the instrument

### WARNING:

- Do not use any cuffs, tubes, parts and accessories other than those explicitly recommended by the manufacturer for use with this product. Parts and accessories not approved for use with the device may cause damage to your health and to the product.

### CAUTION:

- Do not disassemble the unit, it may result in injuries.

### NOTICE:

- The blood pressure monitor is made up of precision electronic components. Accuracy of readings and the instrument's service life depend on careful handling. Protect the unit against hard knocks (e.g. dropping the unit), moisture, water, dirt, dust, chemicals, extreme hot or cold temperatures, major temperature fluctuations, direct exposure to sunlight and heat sources which are too close (e.g. stoves, heating radiators). This may damage the unit. The device must be stored in the specified ambient conditions. Please see section "9. Care,

maintenance and storing the instrument" and section "11. Technical Specifications" for details.

- Do not use any aggressive solvents, cleaning agents, detergents or any other strong chemicals to clean the device.
- Never immerse and or spill water or any other liquid onto the monitor or any components, otherwise liquid will enter it and cause damage.
- Do not wash the unit in a washing machine, dish washer or tumble dryer.
- Observe the local regulations for material disposal. Dispose of the device, components and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.
- Batteries can contain toxins that are harmful to the environment. Always dispose of batteries in accordance with applicable local regulations.
- Never attempt to repair, open and/or disassemble the unit (including arm cuff and optional accessory) or adjust it yourself. This may damage the unit and impair the functions. If you need to have the unit repaired, please contact our customer service. Please see warranty for service contact. Before submitting any complaint, first check the batteries and replace them if necessary.
- No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibration instructions, or other information which will assist the user's appropriately qualified

technical personnel to repair those parts of equipment which are designated repairably can be supplied.

- Changes or modifications to the device will nullify the user warranty.
- Do not drop or insert any object into any opening or hose. This may damage the unit.
- Do not press the buttons with excessive force or with pointed objects.
- The calibration should be checked every 2 years or after repairs. Please contact customer service.

To keep your digital blood pressure monitor in the best condition and protect the unit from damage follow the directions listed below:

- Clean your device and cuff carefully only with a slightly moistened soft cloth and dry it immediately with a soft dry cloth. Do not press.
- When storing the device, make sure that no heavy objects are placed on top of it.
- Do not fold the cuff and tubing tightly. The cuff tube should not have any sharp kinks and keep it away from sharp edges.
- If the unit is stored near freezing, allow it to acclimate to room temperature before use.
- Always store the unit in the pouch after use.

- This unit may not meet its performance specifications if stored or used outside of the temperature and humidity ranges (see section “11. Technical Specifications”).
- Use the unit consistent with the instruction provided in this manual.

## 10. Error messages/trouble-shooting

In case of faults, the  $E_r$  message appears in the display. Error messages may appear if:

- systolic or diastolic pressure could not be measured ( $E_r 1$  or  $E_r 2$  appears on the display),
- systolic or diastolic pressure was outside the measurement range ( $H_i$  or  $L_0$  appears on the display),
- the cuff is fastened too tightly or loosely ( $E_r 3$  or  $E_r 4$  appears on the display),
- the pump pressure is higher than 300 mmHg ( $E_r 5$  appears on the display),
- pumping up takes longer than 160 seconds ( $E_r 6$  appears on the display),
- there is a system or device error ( $E_r A$ ,  $E_r B$ ,  $E_r 7$  or  $E_r 8$  appears on the display),
- the batteries are almost empty .

In such cases, repeat the measurement. Ensure you do not move or speak. If necessary, re-insert or replace the batteries.

## **Technical alarm – description**

Should the recorded blood pressure (systolic or diastolic) lie outside the limits specified in the section “Technical specifications”, the technical alarm will appear on the display indicating either “**Hi**” or “**Lo**”. In such cases, you should seek medical assistance and check the accuracy of your procedure.

The limit values for the technical alarm are factory set and cannot be adjusted or deactivated. These alarm limit values are accorded second priority under the standard IEC 60601-1-8.


The technical alarm is a non-locking alarm and must not be reset. The signal shown on the display will disappear automatically after about 8 seconds.

### **NOTE:**

If the unit still does not work, contact customer service at 1-800-536-0366.

## **11. Technical Specifications**

Model no.	BM 66/U
Type	BM 66
Measurement method	Oscillometric, non-invasive blood pressure measurement on the upper arm

Measurement range	Cuff pressure 0–300 mmHg, systolic 60–260 mmHg, diastolic 40–199 mmHg, Pulse 40–180 beats/minute
Display accuracy	Systolic $\pm 3$ mmHg, diastolic $\pm 3$ mmHg, pulse $\pm 5\%$ of the value shown
Measurement inaccuracy	Max. permissible standard deviation according to clinical testing: systolic 8 mmHg/diastolic 8 mmHg
Memory	2 x 60 memory spaces
Dimensions	L 186 mm x W 95 mm x H 56 mm
Weight	Approx. 330 g (without batteries)
Cuff size	22–42 cm / 8.7” to 16.5”
Permissible operating conditions	+10 °C to +40 °C (+50 °F to +104 °F), $\leq 90\%$ relative air humidity (non-condensing)
Permissible storage and transport conditions	-20 °C to +55 °C (-5 °F to +130 °F), $\leq 90\%$ relative air humidity, 800–1050 hPa ambient pressure
Power supply	4 x 1.5V  AAA batteries

Battery life	For approx. 100 measurements, depending on the blood pressure level and/or pump pressure
Accessories	Cuff, instructions for use, 4 x 1.5V AAA batteries, AC Adapter, storage pouch
Classification	Class II, internal supply, IPX0, no AP or APG, continuous operation, application part type BF
Auto-shut-off	1 minute after measurement, 30 sec. after memory checking and time adjusting

Technical information is subject to change without notification to allow for updates.

- This unit is in line with European Standard EN 60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this unit. More details can be requested from the stated Customer Service address or found at the end of the instructions for use.
- This device is in line with the EU Medical Devices Directive 93/42/EEC, the “Medizinproduktegesetz” (German Medical Devices Act) and the standards EN 1060-1 (non-invasive sphygmomanometers, Part 1: General require-

ments), EN 1060-3 (non-invasive sphygmomanometers, Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems) and IEC 80601-2-30 (Medical electrical equipment – Part 2–30: Particular requirements for the safety and essential performance of automated non-invasive blood pressure monitors).

- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life. If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means. Precise instructions for checking accuracy may be requested from the service address.

## 12. AC Adapter

Model no.	LXCP12-006060DEH
Input	100–240V, 50–60 Hz
Output	6V DC, 600 mA, only in connection with BM66/U blood pressure monitor.
Supplier	Globalcare Co. Ltd.
Protection	This device is double insulated and protected against short circuit and overload by a primary thermal fuse. Make sure to take the batteries out of the compartment before using the AC Adapter.
	Polarity of the the DC voltage connection
	Double insulated/equipment class 2
Enclosures and Protective Covers	Equipment enclosed to protect against contact with live parts, and with parts which can become live (finger, pin, hook test). The operator shall not contact the patient and the output plug of AC mains part simultaneously.

## 13. Glossary

**Blood pressure** – The pressure exerted by circulating blood on the walls of blood vessels; constitutes one of the principal vital signs.

Arterial blood pressure is constantly fluctuating during the course of the cardiac cycle. The highest pressure in the cycle is called the systolic blood pressure, and represents the pressure in the artery when the heart is beating. The lowest pressure is the diastolic blood pressure, and represents the pressure in the artery when the heart is at rest. Both the systolic and the diastolic pressure are necessary for a physician to evaluate the status of a patient's blood pressure.

Many factors such as physical activity, anxiety, stress, eating, smoking, medication or the time of day, can influence your blood pressure. Blood pressure is typically low in the mornings and increases from the afternoon to the evening. It is on average lower in the summer and higher in the winter.

**Electromagnetic Disturbance** – An interruption of the physical field produced by electrically charged objects; affects the behavior of charged objects in the vicinity of the field.

**mmHg** – Millimeter of mercury, the standard unit of measure used for blood pressure measurements.

**Oscillometric Methodology** – A method of pressure measurement using an electronic pressure sensor.

**Pulse Rate** – The frequency of the heartbeat, measured in beats per minute.

**World Health Organization (WHO)** – A specialized agency of the United Nations (UN) that acts as a coordinating authority on international public health.

**Resting Indicator** – This shows whether you, and consequently your circulatory system, are sufficiently at rest when the blood pressure measurement is being taken and is therefore a more precise indicator of your resting blood pressure.

## Life Brand Blood Pressure Monitors

### 14. Warranty and Return Policy

Your Life Brand Blood Pressure Monitor, Model BM 66/U, excluding the cuff (incl. tube), batteries, AC Adapter and any optional accessory, is warranted to be free from defects in materials and workmanship for the life of the product under normal conditions of intended use and service. The cuff is warranted to be free from defects in materials and workmanship appearing within two years from date of purchase when the monitor is used in accordance with the instructions provided with the unit. This warranty extends only to the original retail purchaser and does not extend to retailers or subsequent owners.

The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or three years.

If you have questions regarding the operation of your monitor call the Life Brand Blood Pressure Help Line:

**1-800-536-0366**

**or send an email to: [info@beurer.com](mailto:info@beurer.com)**

Should repair be necessary, return the unit with all component pieces, to your nearest Shoppers Drug Mart location. Please include your name, return address, phone number and email address.

Manufactured for:  
**SHOPPERS DRUG MART®**  
**PHARMAPRIXMD**  
**TORONTO, M2J 4W8**  
®/MD 911979 Alberta Ltd.

Distributed by:  
Beurer North America LP  
900 N Federal Hwy, Ste 300  
Hallandale Beach, FL 33009 USA

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Table 1

## For all ME EQUIPMENT and ME SYSTEMS

### *Guidance and manufacture's declaration – electromagnetic emission*

The BM 66/U is intended for use in the electromagnetic environment specified below. The customer or the user of the BM 66/U should assure that it is used in such an environment.

<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Group 1	The BM 66/U uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The BM 66/U is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

**Table 2**  
**For all ME EQUIPMENT and ME SYSTEMS**

<i>Guidance and manufacture's declaration – electromagnetic immunity</i>			
The BM 66/U is intended for use in the electromagnetic environment specified below. The customer or the user of the BM 66/U should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 s	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the BM 66/U requires continued operation during power mains interruptions, it is recommended that the BM 66/U be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<b>NOTE:</b> $U_T$ is the a.c. mains voltage prior to application of the test level.			

**Table 3**  
**For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

*Guidance and manufacture's declaration – electromagnetic immunity*

The BM 66/U is intended for use in the electromagnetic environment specified below. The customer or the user of the BM 66/U should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 V<sub>rms</sub> 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V<sub>rms</sub></p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the BM 66/U, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance:</b></p> $d = 1.2 \sqrt{P}$ <p><math>d = 1.2 \sqrt{P}</math> 80 MHz to 800 MHz</p> <p><math>d = 2.3 \sqrt{P}</math> 800 MHz to 2.5 GHz</p> <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

*Guidance and manufacture's declaration – electromagnetic immunity*

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BM 66/U is used exceeds the applicable RF compliance level above, the BM 66/U should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BM 66/U.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

**Table 4**  
**For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

*Recommended separation distances between portable and mobile RF communications equipment and the BM 66/U*

The BM 66/U is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BM 66/U can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BM 66/U as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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®TM/MC trademarks of/