beurer herzfrequenz pulsuhr heart rate monitor

german|engineering

PM 25

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Table of contents

1. Scope of delivery	27
2. Important notes	28
3. General information for training	30
4. Functions of the hr monitor	31
4.1 HR functions	31
4.2 Settings	31
4.3 Time functions	31
5. Transmission of signal and methods of devices measurement	31
5.1 Receiving range of the heart rate monitor	31
5.2 Sensors on the chest strap	31
6. Getting started	32
6.1 Putting on the heart rate monitor	32
6.2 Activating the heart rate monitor from the power saving mode	32
6.3 Putting on the chest strap	32
7. General operation of the hr monitor	33
7.1 Buttons on the HR watch	33
7.2 Display	33
7.3 Menus	34
7.4 Power saving mode	34
7.5 Main display	34
7.6 Overview of the menus	34
8. Basic settings	35
8.1 Overview	35
8.2 Enter personal data	35
8.3 Set training zone	36
8.4 Set units	37
9. Watch settings	37
9.1 Overview	37
9.2 Set time of day and date	38
9.3 Set alarm clock	39
10. Recording of training	40
10.1 Overview	40
10.2 Things to know about the recordings	40
10.3 Stop run times	41
10.4 Reset recording	41
11. Results	42
11.1 Overview	42
11.2 Display training data	42
12. Technical specifications and batteries	44
12.1 Replacing the batteries	44
12.2 Distortion and interference	44
12.3 How can I recognize interference? ..	45
12.4 What kind of causes for interference are there?	45
12.5 Interference from other devices for measuring heart rates	45
13. List of abbreviations	45
14. Index	48

1. Scope of delivery

Heart rate monitor with wrist band	
Chest strap	
Elastic tension strap (adjustable)	

In addition:

- Pouch
- Operating instructions and separate short instructions

2. Important notes

Please carefully read through these operating instructions, keep it for future reference and also make it available to others.

Training

- This product is not a therapeutic device. It is a training instrument, which was developed for the measurement and representation of the human heart rate.
- With high-risk sports, please note that the use of the heart rate monitor can represent an additional source of injury.
- If in doubt or for illnesses, please consult your doctor concerning the values of your upper and lower target heart rate as well as the duration and frequency of training. In such a way, you can achieve optimal results when training.
- **WARNING:** People suffering from heart and circulatory illnesses or those with pacemakers should only use this heart rate monitor after having consulted their doctor.

Application

- The device is only intended for private use.
- This device may only be used for the purpose it was designed for and in the manner outlined in this operating instructions. Any form of improper use can be dangerous. The manufacturer is not liable for damages that are caused by improper or incorrect use.
- The heart rate monitor is waterproof. (see chapter "Technical Specifications").
- Electromagnetic disturbances may occur. (see chapter "Technical Specifications").

Cleaning and care

- From time to time, carefully clean your chest strap, elastic chest strap and heart rate monitor with a soap and water solution. Then rinse off all components with clean water. Carefully dry them with a soft towel. You can wash the elastic chest strap in the washing machine at 30° without fabric softener. The chest strap is not made for the dryer!
- Store the chest strap in a clean and dry place. Dirt impairs the elasticity and functioning of the transmitter. Sweat and moisture can keep the electrodes wet and consequently activate the transmitter, which shortens battery life.

Notes on handling batteries

- If your skin or eyes come into contact with battery fluid, flush out the affected areas with water and seek medical assistance.
-  **Choking hazard** Small children may swallow and choke on batteries. Store the batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect the batteries from excessive heat.
-  **Risk of explosion!** Never throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, split or crush the batteries.

Repair, Accessories and Disposal

- Prior to use, it should be ensured that the device and accessories are free of damage. If in doubt, do not use the device and consult your sales representative or the customer service address provided.
- Repairs may only be conducted by customer service or authorized dealers. Please do not in any case attempt to repair the device!
- Use the device only in combination with the accessories provided with it.
- Avoid contact with sunscreen or the like since these can cause damage to the imprint or plastic components.
- For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

3. General information for training

This heart rate monitor serves to measure the human pulse. By means of various settings, you can support your individual training program and monitor your pulse. You should monitor your pulse during exercise because you should prevent your heart from excessive strain on the one hand and on the other in order to achieve the optimal training results. The table provides you with some information on the selection of your training zone. The maximum heart rate reduces with advancing age. Information on the training heart rate must therefore always be in relation to the maximum heart rate. The following rule of thumb is used for the identification of the maximum heart rate:

220 - age = maximum heart rate

The following example applies to a 40-year-old person: $220 - 40 = 180$

	Heart/health programme	Fat burning zone	Fitness programme	Endurance training programme	Anaerobic training
Max. heart rate range	50-60 %	60-70 %	70-80 %	80-90 %	90-100 %
Effect	Strengthens the heart and the circulatory system	The body burns the highest percentage of calories from fat. Strengthens the heart and the circulatory system, enhances fitness	Improves the respiratory and circulatory system. Ideal to promote the basic endurance	Improves speed maintenance and increases basic speed	Overloads specific muscles. High injury risk for athletes at amateur level, risk related to the heart in case of disease
Suitable for	Ideal for beginners	Weight control and loss	Athletes at amateur level	Ambitious athletes at amateur level, professional athletes	Only professional athletes
Training	Regenerative training		Fitness training	Endurance training	Fitness training addressing specific zones

See the chapter "Basic Settings", section "Setting your training zone" on how to individually set your training zone.

4. Functions of the hr monitor

4.1 HR functions

- Exact ECG heart rate measurement
- Transmission: analogue
- Individual training zone adjustable
- Acoustic and visual alarm for exit from the training zone
- Average heart rate (average)
- Maximum training heart rate
- Calorie consumption in Kcal (basis jogging)
- Fat burning in g/oz

4.2 Settings

- Gender
- Weight
- Unit of weight kg/lb
- Height
- Unit of height cm/inch
- 12/24-hour format

4.3 Time functions

- Time of day (12- or 24-hour format)
- Calendar and weekday
- Date
- Alarm / alarm clock
- Stop watch

5. Transmission of signal and methods of devices measurement

5.1 Receiving range of the heart rate monitor

Your heart rate monitor picks up the heart rate signals sent by the chest strap transmitter within a range of 70 centimeters.

5.2 Sensors on the chest strap

The entire chest strap consists of two components: the actual chest strap and an elastic tension belt. On the inside of the chest strap contacting the body, there are two rectangular, corrugated sensors. Both sensors determine your heart rate ECG exactly and transmit this to the heart rate monitor.

6. Getting started

6.1 Putting on the heart rate monitor

You can wear your heart rate monitor like a wrist watch.

6.2 Activating the heart rate monitor from the power saving mode

You can activate your heart rate monitor by pressing any key. If the heart rate monitor has not received any signal for approx. five minutes, it switches to the power saving mode. At that point, only the current time will be displayed. Press any key in order to reactivate the watch.

6.3 Putting on the chest strap

- Tighten the chest strap with the elastic tension belt. Adjust the length of the belt such that the belt is snug, but neither too loose nor too tight. Position the belt around the chest such that the logo is facing outwards and in the correct position directly over the breast bone. For men, the chest belt should be located directly below the pectoral muscles, for women, directly below the breast.
- Since optimal contact between the skin and the heart rate sensors can not be created immediately, it may take some time until the heart rate is measured and displayed. If necessary, alter the position of the transmitter in order to achieve optimal contact. Most of the time, sufficient contact is achieved through the accumulation of sweat beneath the chest strap. You can, however, also wet the contact points on the inside of the chest strap. For this, lift the chest strap up from the skin slightly and wet the two sensors with saliva, water or ECG gel (available at the chemistry). The contact between the skin and chest strap may not be interrupted also with greater movements of the thorax such as during deep breathing. Significant chest chair can interfere with and even prevent contact.
- Put on the chest strap a few minutes before starting in order to warm it up to body temperature and establish optimal contact.

7. General operation of the hr monitor

7.1 Buttons on the HR watch

- | | | |
|---|--|---|
| 1 | | Starts and stops functions and adjusts values. Press and holding accelerates the entry. |
| 2 | | Switches between the functions within one menu and confirms settings. |
| 3 | | Switches between the different menus. |
| 4 | | Pressing quickly illuminates the display for five seconds.
Pressing and holding activates/deactivates the acoustic alarm when exiting the training zone. |

7.2 Display

- | | |
|---|-------------|
| 1 | Top line |
| 2 | Middle line |
| 3 | Bottom line |
| 4 | Symbols |

- ▶ Alarm clock is activated.
- ♥ Heart rate is being received.
- ▲ Training zone has been exceeded.
- ▼ Training zone has been underrun.
- 📶 Acoustic alarm for exit from the training zone is activated.

An explanation of the display abbreviations can be found in the list of abbreviates at the end of the operating instructions.


7.3 Menus

You can switch menus with the  Button.

7.4 Power saving mode

If your heart rate monitor has not received any signal for approx. five minutes, it switches to the power saving mode. The current time and date only are then displayed.


Press any key to reactivate the HR monitor.

7.5 Main display

When changing menus, the name of the menu is first shown, then your HR monitor will automatically switch to the main display of the active menu.

From there, you have to first choose a submenu or you can directly access functions.

7.6 Overview of the menus

Menu	Meaning	Description: see chapter...
Time	Time	"Watch Settings"
Training	Training	"Recording of Training"
Result	Result	"Results"
Setting	Setting	"Basic settings"

8. Basic settings

8.1 Overview

In the menu Setting you can: • Enter your personal data (user). Based on this data, your HR monitor determines your calorie consumption and your fat burning during training! • Set your training zone (limits). See chapter "General Information on Training" for explanations on how to set the training zones. Your HR monitor issues an acoustic and visual alarm during your training when the training zone is exited. • Change the units of measurement (units) for height and weight.	
Once you are in the menu Setting, your HR monitor automatically switches to the main display. This displays User.	

8.2 Enter personal data

MENU	Switch to the menu Setting.	
START/STOP	Gender is displayed. START/STOP Select M, if you are male. Select F, if you are female.	
OPTION/SET	Weight is displayed. START/STOP Set your height (80-227 cm or 2'66"-7'6" feet).	

OPTION/SET	Height is displayed. START/STOP Set your height (80-227 cm or 2'66"-7'6" feet).	 Height
OPTION/SET	Return to the display User .	

8.3 Set training zone

You can set your individual upper and lower heart rate limit as described below.

MENU	Switch to the menu Setting .	 User
OPTION/SET	Limits is displayed.	
START/STOP	Lo limit is displayed. START/STOP Set the lower limit for your training zone (40-238).	 Lo limit
OPTION/SET	Hi limit is displayed. START/STOP Set the upper limit (41-239).	 Hi limit
OPTION/SET	Return to the display Limits .	

8.4 Set units

MENU	Switch to the menu Setting .	138 User
OPTION/SET	Limits is displayed.	
OPTION/SET	Units is displayed.	
START/STOP	kg/lb is displayed. The weight unit blinks. START/STOP Switch between the units.	49 kg/lb
OPTION/SET	cm/inch is displayed. The unit of measurement for height blinks. START/STOP Switch the units.	cm cm/inch
OPTION/SET	Return to the display Units .	

9. Watch settings

9.1 Overview

In the Time menu you can: • Set the time and the date. • Set the alarm clock.	
When you are in the Time menu the heart rate monitor automatically activates the main display. In this mode both the weekday and the current date are displayed. On the display: Fri 12.04 (Friday 12th April).	

9.2 Set time of day and date

MENU	Switch to the menu Time.	
START/STOP 5 sec	Set Time is displayed briefly. 24hr is displayed and flashes. START/STOP Changing the hour format (12 or 24 hours). Notes: • Display of the date for 24h format: “Weekday day.month”. • Display of the date for 12h format: “Weekday month.day”. Display AM (morning), e.g. 2:00 AM means 2:00 . Display PM (afternoon), e.g. 2:00 PM means 14:00.	
OPTION/SET	Hour is displayed. The hour setting flashes. START/STOP Set the hour: • 0 - 24 for 24h format • 0 - 12 for 12h format	
OPTION/SET	Minute is displayed. The hour setting blinks. START/STOP Set the minute (0-59).	
OPTION/SET	Year is displayed. The year setting blinks. START/STOP Set the year. The calendar goes up to 2099.	

OPTION/SET	Month is displayed. The monthly setting blinks. START/STOP Set the month (1-12).	
OPTION/SET	Day is displayed. The day setting blinks. START/STOP Set the day (1-31).	
OPTION/SET	Back to the main display.	

9.3 Set alarm clock

You can turn off the alarm clock by pressing any button. On the following day, we will be woken up at the same time.

MENU	Switch to the menu Time .	
OPTION/SET	Alarm is displayed.	
START/STOP 5 sec	Set Alarm is shown briefly. On Alarm is displayed. START/STOP Turns the alarm clock off (Off) and on (On).	
OPTION/SET	Hour is displayed. The hour setting blinks. START/STOP Set the hour.	
OPTION/SET	Minute is displayed. The hour setting blinks. START/STOP Set the minute.	
OPTION/SET	Back to the main display.	

10. Recording of training

10.1 Overview

In the menu Training you can: • Stop your run times and thereby record your training session. • Reset recording to zero.	 138  00:00.00 Training
Once you are in the menu Training, your HR monitor automatically switches to the main display. This displays the already recorded training time on the middle line (after reset: 00:00:00) and on the line below Start.	 138  00:00.00 Start

10.2 Things to know about the recordings

- As soon as you have started the recording, your HR monitor starts recording your training data. You can have your training data displayed in the menu Result.
- Your HR monitor continually records your training data until you have reset this. This also applies when you first start training on the following day. The running time also allows you to stop your training session at any time.
- Reset your recording each time you would like to begin a new recording. The memory of your heart rate monitor can record 23:59:59 hours

10.3 Stop run times

MENU	Switch to the menu Training . If you have already stopped a run time, the recorded training time is displayed on the middle line.	
START/STOP	Start the recording. Run is displayed. Even if you change to another menu during the recording, the time continues to run in the background.	
START/STOP	Stop the recording. Stop is displayed. You can interrupt and resume your recordings at any time by pressing the button START/STOP .	

10.4 Reset recording

	You are in the menu Training . No training is currently being recorded. Stop is displayed.	
START/STOP 5 sec	Reset the recording. Hold Reset is displayed.	
	The display automatically switches back to the main display.	

11. Results

11.1 Overview

In the menu Result you can: Display your training data: mean and maximum heart rate, times within and outside of the training zone as well as calorie consumption and fat burning.	138 Result
Once you are in the menu Result, your HR monitor automatically switches to the main display. This displays bPM AVG.	120 bPM AVG

11.2 Display training data

MENU	Switch to the menu Result .	
	BPM AVG (average training heart rate) is displayed.	120 bPM AVG
OPTION/SET	HRmax (maximum training heart rate) is displayed.	180 bPM HRmax
OPTION/SET	In [min] (time within the training zone) is displayed.	12 1:15 [min]

OPTION/SET	Lo [min] (time below the training zone) is displayed.	Lo 0:13 [min]
OPTION/SET	Hi [min] (time above the training zone) is displayed.	Hi 0:09 [min]
OPTION/SET	[kcal] (the total calorie consumption during the recording of the training session) is displayed. The basis for the display of the calorie consumption and the fat burning are: • Personal data (basic settings) • the time and heart rate measured during the training session. The calculation is conducted based on moderately fast jogging.	128 [Kcal]
OPTION/SET	Fat[g] (Fat burning during the entire training recording) is displayed.	8 Fat[g]
OPTION/SET	Back to the main display.	

12. Technical specifications and batteries

Heart rate monitor	Water-resistant up to 30 m (suitable for swimming).  The water resistance of the device cannot be guaranteed if you jump into the water (high pressure loads) or press the buttons in the rain.
Chest strap	Water-resistant (suitable for swimming).
Battery heart rate monitor	3V Lithium battery, type CR2032 (Recommendation: Energizer batteries). battery life: approx. 20 months (daily use of approx. 1 hour)
Battery chest strap	3V Lithium battery, type CR2032. battery life: approx. 28 months (daily use of approx. 1 hour)

12.1 Replacing the batteries

The battery of the HR monitor may only be replaced by a specialized dealer (any lock-maker's shop) and at the service address provided. Otherwise, the claims to warranty become void.

You can replace the battery in the chest strap on your own. The battery compartment is located on the inner side. Use a coin or the battery to be used to open the compartment. The positive pole of the battery inserted must face upwards. Carefully close the battery compartment again.

- The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.
- The codes below are printed on batteries containing harmful substances:
Pb = Battery contains lead,
Cd = Battery contains cadmium,
Hg = Battery contains mercury.

12.2 Distortion and interference

In the HR monitor, a receiver for signals from the chest strap transmitter is integrated. During use, it may arise that the HR monitor receives signals other than those transmitted by the chest strap.

12.3 How can I recognize interference?

Interference can be recognized in that the values displayed are unrealistic. If your normal heart rate amounts to 130 bpm and it is suddenly shown to be 200, then a short-term disturbance is present. Once you have removed yourself from the vicinity of the source causing the disturbance, the values will automatically resume being accurate. Please note that distorted values can mildly influence the training result displays and overall values.

12.4 What kind of causes for interference are there?

In some situations, short-range electromagnetic interference can occur. In order to avoid these, keep well away from these sources.

- Nearby power lines
- Walking through security gates in airports.
- Nearby overhead lines (tram)
- Nearby traffic lights, television sets, car motors, bike computers
- Nearby exercise equipment
- Nearby mobile phones +

12.5 Interference from other devices for measuring heart rates

The HR monitor as a receiver can react to interfering stray signals such as other heart rate monitors within a vicinity of one metre. This can result in incorrect measurements particularly when training in groups. If you e.g. job next to one another, it can be beneficial to wear the HR monitor on the other wrist or to increase the distance to one's training partner.

13. List of abbreviations

Abbreviations on the display in alphabetical order:

Abkürzung	Deutsch
12 hr	12-hour format AM/PM
24 hr	24-hour format
Alarm	Alarm
AM	Ante meridiem: time between midnight and midday
AVG	Average heart rate during training
BPM	Heart beats per minute
bPM AVG	Average heart rate during training, beats per minute

Abkürzung	Deutsch
bPM Hrmax	Maximal heart rate during training, beats per minute
cm/inch	Centimeter/Inch
Day	Day
f	female
Fat[g]	Fat in gram
Fri	Friday
Gender	Gender
Height	Height
Hi	Time above the training zone
Hi limit	Highest training limit
Hour	Hour
HR	Heart Rate
HRmax	Max. heart rate during training
In	Time within the training zone
[Kcal]	Calorie consumption in Kcal
kg	Kilogram
kg/lb	Kilogram/pound
Limits	Training limits
Lo	Time below the training zone
Lo limit	Lowest training limit
m	Lowest training limit
Minute	Minute
[min]	Minute
Mon	Monday
Month	Month

Abkürzung	Deutsch
PM	Post meridiem: time between midday and midnight
Reset	Reset values to zero
Result	Result
Run	Run time is measured
Sat	Saturday
Setting	Settings (Menu)
SEt	Setting
Start	Start
Stop	Stop
Sun	Sunday
Time	Time (Menu)
Thu	Thursday
Training	Training (Menu)
Tue	Tuesday
Units	Units
User	User
Wed	Wednesday
Weight	Weight
Year	Year

14. Index

A		R	
Accessories	29	Recording	40
Alarm clock	39	Recording of training	40
Application	28	Repair	29
B		Replacing the batteries.....	44
Basic settings	35	Reset	41
Buttons	33	Reset recording	41
C		Reset stopwatch	41
Care	28	Result (Menu).....	34
Cleaning	28	Results.....	42
D		S	
Display	33	Scope of delivery.....	27
Display training data	42	Set alarm clock.....	39
Disposal.....	29	Set date	38
Distortion and interferences	44	Set time of day and date	38
G		Set training zone.....	36
Getting started.....	32	Set units	37
L		Setting (Menu)	34
List of abbreviations	45	Stop run times	41
M		Stopwatch	41
Main display	34	Symbol	33
Menu	34	T	
Methods of devices measurement ...	31	Technical specifications	44
N		Time (Menu).....	34
Notes	28	Training	30
P		Training (Menu)	34
Pacemakers.....	28	Transmission of Signal from the	
Power saving modus	34	Devices.....	31
		W	
		Watch settings.....	37

beurer herzfrequenz pulsuhr

Beurer GmbH • Söflinger Straße 218 • 89077 Ulm, Germany • www.beurer.com
www.beurer-gesundheitsratgeber.com • www.beurer-healthguide.com

754.601-0117