Women's Wrist Digital Blood Pressure Monitor

Item # 04-825-001

Please read this guidebook
completely before operating this unit.

**DETAILED
GUIDEBOOK**

English • Spanish

Limited Five-Year Warranty

The warrantor guarantees that its Digital Blood Pressure Monitor will be free from manufacturing defects under normal use for five years from the original purchase date. This warranty covers only normal use, and does not apply to use in any clinical or commercial applications. This warranty does not cover batteries or other power sources that may be provided with or used with the Digital Blood Pressure Monitor. This warranty is voided if the Digital Blood Pressure Monitor product is misused or abused in any manner.

If the Digital Blood Pressure Monitor fails to operate during the first five years from the time originally purchased, please call our Customer Care Help Line at 800-622-4714 and we will try to resolve the issue as soon as possible. If the warrantor determines the unit failed to operate due to a manufacturing defect, the unit will be replaced at the option of the warrantor. Replacement is the sole remedy under this limited warranty. This warranty gives you specific legal rights, which vary from state to state. As a condition of this warranty, the enclosed warranty registration card must be completed and sent to us within 10 days of the purchase date or can be filled out online at www.mabisdmi.com.

This warranty constitutes the warrantor's only responsibility and obligation to replace materials or components. We make no other express or implied warranties, arising by operation of law or otherwise, or any warranty of merchantability or fitness for a particular use or purpose whether or not the use or purpose has been disclosed to the warrantor in specifications, drawings or otherwise, and whether or not the warrantor's products are specifically designed and/or manufactured by the warrantor for the buyer's use or purposes, except for the limited warranty stated above. The warrantor will not be responsible for any indirect, incidental, special, consequential, or punitive damages or other loss, including, but not limited to, damage to or loss of other property or equipment and personal injuries, whether to purchaser or others. The warrantor shall in no event be liable to the purchaser for any amount in excess of the cost of repair and/or replacement of the unit.

Index

Introduction	2
Product Identification and Caution	3
General Blood Pressure Information	4-6
Important Information Before Use	7
World Health Organization (WHO) Indicator	8
Irregular Heartbeat Detection	9
Battery Installation/Replacement	10
Setting the Date/Time	11
Applying Your Blood Pressure Cuff	12-13
Taking Your Blood Pressure Reading	14-15
Recalling Measurements in Memory	16
Description of Display Symbols	17-18
Troubleshooting	19
Care and Maintenance	19
Product Specifications	20
Sample Blood Pressure Log	21
Blood Pressure Log	22

Introduction

To achieve the maximum benefit from your blood pressure monitor, we recommend that you first consult with your physician or trained health care professional.

Thank you for purchasing an Automatic Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.


The method of measurement that your Automatic Monitor uses is called the oscillometric method. The monitor detects your blood's movement through the artery in your wrist and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor easy-to-use. The preformed cuff requires little effort to apply and is comfortable to wear.

Clinical research has proven a direct relationship between blood pressure in the wrist and blood pressure in the arm. Because the arteries in the wrist and arm are connected, changes in wrist blood pressure reflect changes in arm blood pressure. Frequent wrist blood pressure measurements will provide you and your physician with an accurate indication of change in your true blood pressure.

People with severe hypertension, severe arteriosclerosis or severe diabetes should consult their physician about monitoring blood pressure at the wrist.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained observer using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Product Identification and Caution

CAUTION

While taking a measurement, you can stop the inflation or deflation process of the cuff at any time by pressing the POWER button.

General Blood Pressure Information

What is Blood Pressure?

Blood pressure is the pressure that is exerted by blood flowing against the walls of the blood vessels throughout your body.

Your heart, which is the center of the circulatory system, provides the force for the blood to flow or circulate. When your heart contracts, or “beats”, the blood is forced through the blood vessels increasing the pressure. This is the highest pressure in the cycle or what is referred to as **SYSTOLIC BLOOD PRESSURE**. In between beats, your heart relaxes and your blood pressure decreases. This is referred to as **DIASTOLIC BLOOD PRESSURE**.

This complete series of events which occurs in a single heartbeat is known as the **CARDIAC CYCLE**.

Your Digital Blood Pressure Monitor will automatically read your blood pressure and display both systolic and diastolic readings on the screen. Your systolic will be positioned as the upper number and the diastolic reading will be the lower number.

(systolic) **120/80** (diastolic)

Blood pressure is measured in millimeters (mm) of Mercury (Hg) and is generally recorded with the systolic pressure (120) listed first and the diastolic pressure (80) listed second. The numbers are typically separated by a slash mark (/) as shown above.

*Both pressure readings, the **SYSTOLIC** and the **DIASTOLIC**, are necessary for a physician to evaluate the status of a patient's blood pressure.*

Please contact your physician for specific information regarding your own blood pressure.

What Influences Blood Pressure?

Many factors such as genetics, age, sex, altitude, physical activity, anxiety, muscular development, certain medications or even the time of day can influence blood pressure. Influences such as sleep or relaxation decreases blood pressure, while anxiety or exercise increases blood pressure.

Why Monitor Your Blood Pressure at Home?

A visit to a physician's office can be a stressful situation for a patient. Anxiety is a known factor in raising blood pressure. This temporary occurrence of elevated blood pressure at the physician's office is commonly referred to as “white coat syndrome”.

Whether or not you experience “white coat syndrome”, home blood pressure monitoring provides you with the opportunity to supplement your physician's office measurements. These home readings, when taken over a period of time, can show an accurate indication of change. Furthermore, your records can assist your physician in evaluating your health and in making important decisions in the diagnosis and treatment of your condition. Because of this, it is important to take consistent, daily measurements of your blood pressure.

The variations in your individual readings should only be interpreted by your physician or a trained health care professional.

WHO Blood Pressure Classifications


Standards for assessment of high or low blood pressure, regardless of age, have been established by the World Health Organization (WHO) as shown in the chart on the following page.

CATEGORY	SYSTOLIC (mmHg)	DIASTOLIC (mmHg)	COLOR INDICATOR
Optimal	<120	<80	GREEN
Normal	120-129	80-84	GREEN
High - Normal	130-139	85-89	GREEN
Stage 1 - Mild Borderline Hypertension	140-159	90-99	YELLOW
Stage 2 - Moderate Hypertension	160-179	100-109	ORANGE
Stage 3 - Severe Hypertension	≥180	≥110	RED
Isolated Systolic Hypertension	≥140	<90	
	140-159	<90	YELLOW
	160-179	<90	ORANGE
	≥180	<90	RED

This chart is only a general guideline. Please contact your physician to determine your **NORMAL BLOOD PRESSURE**.

Variations in Blood Pressure

Blood pressure can be influenced by many factors. This change can occur from moment to moment. Normally, blood pressure is lowest during the sleeping period and rises during the day. The graph below represents the variations in blood pressure shown over a day with measurements taken every 5 minutes.

The dotted line represents the sleep period. The rise in blood pressure at 1 p.m. (A in the graph) corresponds to a stressful occurrence and at 7 p.m. (B in the graph) a period of exercise.


Important Information Before Use

- Blood pressure measurements should only be interpreted by a physician or a trained health care professional who is familiar with your medical history. Through regular use of this device and recording of your measurements, you can keep your physician informed of the changes in your blood pressure.
- Perform your measurement in a quiet place. You should be seated in a relaxed position.
- Avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading. If you are exhibiting signs of stress, avoid taking your measurement until the feeling subsides.
- Rest 15 minutes prior to taking a reading.
- Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
- Apply the cuff to the proper position on the left wrist. The unit should be located on the inside of your left wrist with the display panel in view.
- Keep the monitor stable during measurement to achieve an accurate reading. Remain still; do not talk during the measurement.
- Record your daily readings on the chart provided.
- Take your readings at the same time, each day (or as recommended by your physician) to get an accurate indication of change in your true blood pressure.
- Wait a minimum of 15 minutes between readings to allow for the blood vessels to return to normal. The wait time may vary depending on your individual physiological characteristics.
- Although such cases are rare, for those with an extremely weak pulse or irregular pulse, errors may result which prevent proper measurement. If abnormal variations are noticed, consult with your physician or trained healthcare professional.
- This device is intended for adult use.

While taking a measurement, you can stop the inflation or deflation process of the cuff at any time by pressing the **POWER** button.

World Health Organization (WHO) Indicator

The WHO Indicator, located on the left side of the LCD screen, will identify your reading within a range in accordance with the guidelines set by the World Health Organization (WHO). See the chart below for reference.

SYS  DIA 	Stage 3 Hypertension Systolic pressure: ≥ 180 mmHg -or- Diastolic pressure: ≥ 110 mmHg
SYS  DIA 	Stage 2 Hypertension Systolic pressure: 160-179 mmHg -or- Diastolic pressure: 100-109 mmHg
SYS  DIA 	Stage 1 Hypertension Systolic pressure: 140-159 mmHg -or- Diastolic pressure: 90-99 mmHg
SYS  DIA 	High Normal Systolic pressure: 130-139 mmHg -or- Diastolic pressure: 85-89 mmHg
SYS  DIA 	Normal Systolic pressure: 120-129 mmHg -or- Diastolic pressure: 80-84 mmHg
SYS  DIA 	Optimal Systolic pressure: <120 mmHg -and- Diastolic pressure: <80 mmHg

Irregular Heartbeat Detection

Your digital blood pressure monitor features Irregular Heartbeat Detection. This feature allows users to accurately monitor blood pressure even if an irregular heartbeat should occur.

When an irregular heartbeat is detected, the IHB icon will appear on the display, *Fig. 1*.

NOTE: Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat if this symbol appears frequently at the end of a reading.


Fig. 1


Battery Installation/Replacement

Replace the batteries when the Low Battery symbol appears on the display, *Fig. 1*, or when the display remains blank when the POWER Button is pressed.

To change the batteries:

1. Open the battery cover, *Fig. 2*.
2. Making sure the battery polarity is correct, install two AAA alkaline batteries, *Fig. 3*.
3. Close the battery cover, *Fig. 4*.

NOTE: It is recommended to remove the batteries if the unit will not be used for an extended period of time. Rechargeable batteries are not suitable for this monitor.

Setting the Date/Time

NOTE: When the unit is off, the date and time will still appear at the top of the screen.

1. When the unit is off, press both the POWER and MEMORY button at the same time for two seconds. The month will begin to flash, *Fig. 1*.
2. Press the MEMORY button to change until the desired number appears.

NOTE: To scroll through the numbers quicker, hold down the MEMORY button.

3. Press the POWER button to set the entry.
4. The date, *Fig. 2*, followed by the hour and minutes will blink. Repeat steps 2 and 3 to adjust the date/time setting.
5. After setting the minutes, the unit will automatically exit out of the date/time setting mode, *Fig. 3*.

NOTE: The time is expressed in 24hr military time. The chart below explains the time conversion.

Regular Time (12 hours)	Military Time (24 hours)	Regular Time (12 hours)	Military Time (24 hours)
12 am Midnight	0	12 pm Noon	12
1 am	1	1 pm	13
2 am	2	2 pm	14
3 am	3	3 pm	15
4 am	4	4 pm	16
5 am	5	5 pm	17
6 am	6	6 pm	18
7 am	7	7 pm	19
8 am	8	8 pm	20
9 am	9	9 pm	21
10 am	10	10 pm	22
11 am	11	11 pm	23


Fig. 1


Fig. 2


Fig. 3

Applying Your Blood Pressure Cuff

Cuff Dimensions

3-1/4" x 2-3/4" x 1-1/8"

Wrist Circumference Range

Suitable for 5-1/2" – 7-5/8" (14 – 19.5 cm)

Avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading. If you are exhibiting signs of stress, avoid taking your measurement until the feeling subsides. Rest 15 minutes prior to taking a reading.

1. Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
2. You should be seated with your feet flat on the floor.
3. Hold your left arm in front of you with your palm facing upward.
4. Apply the preformed cuff to your left wrist. The digital display panel should be positioned on the inside area of your wrist facing you, *Fig. 1*.
5. Adjust the cuff as shown in *Fig. 2* from the edge of the 'Head of the Ulna Bone.'


Fig. 1


Fig. 2

Applying Your Blood Pressure Cuff

6. Wrap and fasten the cuff securely, yet comfortably around your wrist, *Fig. 3*.

NOTE: The left wrist should be used unless you are unable to or should not use your left arm or wrist. In consultation with your physician, you may modify the following instructions and apply the cuff to your right wrist.


Fig. 3

Taking Your Blood Pressure Reading

Proceed only after reading the previous sections of this manual.


Fig. 1

1. Rest your elbow on a solid surface with your palm facing upward. Elevate your hand so that the cuff is at the same level as your heart, Fig. 1. Relax your left hand.
2. This unit will beep and run a self-test, Fig. 2, and show the last measured reading, Fig. 3, or if there is no previous measurement, "0" as shown in Fig. 4.
3. When the display appears as shown, Fig. 5, the monitor is ready. The cuff will automatically begin to inflate.
4. The display will show the increasing pressure in the cuff. Along with the measurement progress on the left side of the screen.
5. The cuff will automatically inflate to approximately 190 mmHg or 40 mmHg higher than the last systolic pressure stored in memory.
7. Once the cuff inflates to the appropriate pressure level, the measurement will begin.
8. The unit will begin to automatically deflate and the numbers appearing on the display will decrease.

Fig. 2


Fig. 3


Fig. 4

Fig. 5

Taking Your Blood Pressure Reading

9. When all the air is deflated, the measurement is complete. The unit will display your blood pressure (systolic and diastolic) measurements and pulse rate, Fig. 6.
10. Record your readings on a chart.
11. At this point, your measurements will be automatically saved in memory.
12. To conserve the batteries, press the POWER button to turn the power off. Otherwise, the unit will automatically shut off after approximately 1 minute.

Fig. 6

NOTE: This monitor inflates the cuff to approximately 190 mmHg or 40 mmHg higher than the last systolic pressure stored in memory. If the system detects that the cuff needs additional pressure, it will automatically reinflate to approximately 40mmHg higher than the last pressure level.

Recalling Measurements in Memory

You can recall up to 60 measurements in the memory. These measurements can be shared with your physician or trained healthcare professional.

1. Press and release the MEMORY button. It will show how many measurements are stored within the bank, followed by the last result, Fig. 1.

Fig. 1

2. Press the MEMORY button again to scroll through the saved measurements. After you've scrolled through all the readings, the unit will automatically turn off.

NOTE: When the number of readings exceeds 60, the oldest data will be replaced with the new record.

To erase all stored measurements:

1. Press and release the MEMORY button. The unit will display how many measurements are stored within that bank.
2. Press and hold down the Memory button until two horizontal lines flash. This indicates that all measurements within the selected memory bank have been erased, Fig.2.



Fig. 2

Description of Display Symbols

Indicates a measurement within Stage 3 Hypertension.

Indicates a measurement within Stage 2 Hypertension.

Indicates a measurement within Stage 1 Hypertension.

Appears when batteries should be replaced.

-- Description of Display Symbols Continued

Description of Display Symbols

Shows the pulse rate per minute.

Occurs when an error was made during measurement.

Appears when there is an irregular heartbeat detected during measurement.

Appears when the cuff is unstable or it has too much air left in it. Once the arrow stops flashing, the unit is ready to take a measurement.

Appears when a measurement error has occurred. Restart the measurement.

Appears when recalling a measurement from memory.

Troubleshooting

If any abnormality occurs during use, please check and correct the following:

Condition

Display is blank when power is on.

Measurement incomplete or abnormally low or high values displayed.

Measurements are different from those typically measured by physician or every measurement reading is different.

EE Symbol displays

Correction

Check and correct the polarity of the installed batteries. Reinstall or replace batteries.

Review and follow 'Applying Your Blood Pressure Cuff' and 'Taking Your Blood Pressure Reading' sections. Replace the batteries.

Remember, blood pressure readings are influenced by physical and mental conditions and/or even the time of day. Daily measurement recordings should be interpreted by your physician.

Do not talk or move during the measurement. Review "Applying Your Blood Pressure Cuff" section.

Care and Maintenance

1. Only use a soft, dry cloth to clean your blood pressure monitor.
2. Avoid using any types of liquids on the monitor or cuff.
3. Store your monitor in the provided storage case when not in use.
4. Do not store the unit where it will be exposed to direct sunlight, dust or humidity. Avoid extreme temperatures.
5. Never disassemble the monitor or cuff.
6. Dropping or subjecting your blood pressure monitor to strong shocks should be avoided.
7. Remove batteries prior to storing.

Product Specifications

Name and Model Number	HealthSmart Women's Wrist Blood Pressure Monitor, 04-825-001
Display System	Digital display/LCD
Measuring Method	Oscillometric Method Automatic Air Inflation and Measurement
Power Source	2 AAA batteries
Measuring Range	Pressure: 0 to 300 mmHg Pulse: 40-180 beats/minute
Accuracy	Pressure: ± 3 mmHg; Pulse: $\pm 5\%$ of reading
Pressurization	Automatic pressurization by pump
Deflation	Automatic pressure release valve
Memory	Built-in memory enabling display of up to 60 measurements
Automatic Power Off	Approximately 1 minute after last button operation
Battery Life	Approximately 1 month with normal usage
Operation Environment	Temperature 41°F – 104°F (5°C – 40°C) Humidity <85% RH max
Storage Environment	Temperature -4°F – 131°F (-20°C – 55°C) Humidity <95% RH max
Weight	4.3-oz. (without batteries)
Dimensions	3-5/8" x 2-7/8" x 1-1/8"
Wrist Circumference Range	Suitable for 5-1/2" – 7-5/8" (14 – 19.5 cm)
Accessories	Storage case, Instruction Guidebook, two AAA alkaline batteries, Quick Start Guide

Specifications are subject to change without notice.

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1931 Norman Drive South • Waukegan, IL 60085
Toll-Free Customer Care Help Line: 1-800-622-4714
Monday - Friday 8:00 am - 4:30 pm CST

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