FORA

**FORA D15
BLOOD GLUCOSE
PLUS BLOOD PRESSURE
MONITORING SYSTEM**

English

Thank you for choosing

FORA D15
BLOOD GLUCOSE PLUS
BLOOD PRESSURE
MONITORING SYSTEM

Dear FORA D15 System Owner:

This manual contains important information you must know about the system. Please read it thoroughly and carefully.

This product can assist you in monitoring your blood glucose level as well as your blood pressure.

For other questions regarding this system, please contact our customer service team at 1-888-307-8188, 1-866-469-2632. You can also visit www.foracare.com/usa for helpful information.

IMPORTANT SAFTY INSTRUCTIONS READ THIS BEFORE USING

The following basic safety precautions should always be taken.

1. Close supervision is necessary when the device is used by, on, or near children, handicapped persons or invalids.
2. Use the device only for the intended use described in this manual.
3. Do not use accessories which are not supplied by the manufacturer.
4. Keep the equipment and its flexible cord away from hot surfaces.
5. Do not use the equipment where aerosol sprays are being used, or where oxygen is being administered.
6. Do not use device if it is not working properly or has been damaged.
7. Before using product to test your blood glucose, read all instructions thoroughly and practice the test. Do all quality control checks as directed and consult with a diabetes healthcare professional.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

TABLE OF CONTENTS

IMPORTANT INFORMATION	06
Warnings and Precautions	06
Health Information	07
Alternative Site Testing (AST)	09
INTRODUCTION TO THE SYSTEM	11
Intended Use	11
Principle of Measurement	11
Major Features	12
Contents of the System	13
Appearance and Key Function of the Monitor	14
Monitor Display	15
Test Strip	16
BEFORE USE	17
Battery Replacement	17
Setting the Monitor	19
BLOOD GLUCOSE MEASUREMENT	21
Important Information	21
Coding the Monitor	22
Control Test	23
Important Control Solution Information	24
Performing a Control Solution Test	25
Test Procedure	28

Comparing Monitor and Laboratory Results	33
Risk Message	35

BLOOD PRESSURE MEASUREMENT 36

Suggestions Before Measuring	36
Correct Application of the Pressure Cuff	36
Testing Your Blood Pressure	38

USING THE MONITOR MEMORY 40

Viewing Results on the Monitor	40
Deleting Memory	42

VIEWING RESULTS ON A PERSONAL COMPUTER 43

TAKING CARE OF YOUR MONITOR AND TEST STRIPS 49

Cleaning	49
Storage	49

TROUBLESHOOTING GUIDE 52


Error Message	53
Operating Problem	54

SPECIFICATIONS 59

FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT 61

IMPORTANT INFORMATION

Warnings and Precautions

- The FORA D15 System is designed for use on individuals **age 16 and above**. It shall **NOT** be used under any circumstances on newborns, infants or persons who cannot communicate.
- This device does **NOT** serve as a cure of any symptoms or diseases. The data measured are for reference only. Always consult your physician to have the results interpreted.
- This device is **NOT** able to measure in the presence of common arrhythmia, such as arterial or ventricular premature beats or arterial fibrillation. It may produce reading error.
- Do **NOT** use the device for purposes other than measuring blood glucose and blood pressure for human beings.
- Do **NOT** wrap the pressure cuff around anything other than your arm.

Health Information

Blood Glucose

Blood glucose monitoring plays an important role in diabetes control. A long-term study showed that **keeping blood glucose levels close to normal value** can reduce the risk of diabetes complications by up to 60%*1. The results you get with the FORAD15 system can help you and your healthcare professional monitor and adjust your treatment plan to gain better control of your diabetes.

Time of day	Plasma glucose range (mg/dL) for people without diabetes*	Your target range (mg/dL)
Fasting and before meal	Less than 100 mg/dL (5.6 mmol/L)	
2 hours after meals	Less than 140 mg/dL (7.8 mmol/L)	

*1: American Diabetes Association position statement on the Diabetes Control and Complications Trial (1993).


* Source: American Diabetes Association (2010). Standards of Medical Care in Diabetes. Diabetes Care, 33 (Supplement 1): S1-100.

NOTE

Please work with your doctor to determine a target range that works best for you.

Blood Pressure

Clinical studies show that adult diabetes is often accompanied by elevated blood pressure. People with diabetes can reduce their heart risk by managing their blood pressure along with diabetes treatment*². Knowing your routine blood pressure patterns can help determine whether or not your body is in good condition. Human blood pressure naturally increases after reaching middle age. This symptom is a result of continuous aging of the blood vessel. Other causes include obesity, lack of exercise, and cholesterol (LDL) adhering to the blood vessels. Rising of blood pressure accelerates hardening of the arteries, making the body become more susceptible to apoplexy and coronary infarction. The WHO (world health organization) published the following guideline for blood pressure range:

Source: 1999 WHO/ISH guidelines for the management of hypertension.

*²: American Diabetes Association: The Diabetes-Heart Disease Link Surveying Attitudes, Knowledge and Risk (2002).


Alternative Site Testing (AST)

Important: There are limitations for doing AST.

Please consult your healthcare professional before you begin AST.

What is AST?

Alternative site testing (AST) refers to when individuals check their blood glucose levels using areas of the body other than fingertip. This system allows AST at the palm, the forearm, the upper arm, the calf or the thigh with results as reliable as those obtained from fingertip testing.

What's the advantage?

Fingertips feel pain more readily because they are full of nerve endings (receptors). Since nerve endings are not so condensed at other body sites, tests performed there will be less painful.

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at the fingertip reflects these changes faster than capillary blood at other sites. Therefore, when testing blood glucose during or immediately after a meal or physical exercise, **results from AST and fingertip may be significantly different.**

We strongly recommend you perform AST **ONLY** during the following intervals:

- In a pre-meal or fasting state (more than 2 hours since the last meal).
- Two or more hours after taking insulin.
- Two or more hours after exercise.

Do NOT use AST if:

- You think your blood glucose is low.
- You are unaware of hypoglycemia.
- Your AST results do not correspond the way you feel.
- You are testing for hyperglycemia.
- Your routine glucose results are often fluctuating.

How to increase the accuracy?

Stimulating blood perfusion by rubbing the puncture site prior to blood extraction has a significant influence on the glucose value obtained. Blood from a site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site is rubbed prior to blood extraction, the difference is markedly reduced.

Please follow suggestions below before puncturing the skin:

- Rub the puncture site approximately 20 seconds before penetration.
- Use a clear cap (included in the kit) while setting the lancing device.

INTRODUCTION TO THE SYSTEM

Intended Use

The FORA D15 Blood Glucose plus Blood Pressure Monitoring System is intended for the quantitative measurement of glucose in fresh capillary whole blood (from the finger, the palm, the forearm, the upper-arm, the calf and the thigh) for self testing by persons with diabetes in the home or by healthcare professionals in healthcare facilities. The system is also intended to be used to measure non-invasively the systolic and diastolic blood pressure and pulse rate of an adult individual over age 16, at home.

Principle of Measurement




Blood glucose is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The monitor measures the current and displays the corresponding blood glucose level. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

Blood pressure is measured non-invasively at the arm based on the Oscillometric method.


Both functions work independent of each other (one measurement, either blood glucose or blood pressure, at a time) in order to avoid any interference problems.

Major Features

Four subtypes of FORA D15 were developed. Each uses a slightly different method of data transmission as the following table indicates.

MODEL	DATA TRANSMISSION FUNCTION
<p>FORA D15c</p>  <p>The diagram shows a handheld device labeled 'FORA D15c' connected to a desktop computer labeled 'PC' by a cable. The device has a screen and several buttons.</p>	<p>You can view results on the PC through cable connection.</p>
<p>FORA D15b</p>  <p>The diagram shows a handheld device labeled 'FORA D15b' with a wireless signal icon above it, connected to a desktop computer labeled 'PC'.</p>	<p>You can view results on the PC through a wireless connection using Bluetooth technology.</p>
<p>FORA D15z</p>  <p>The diagram shows a handheld device labeled 'FORA D15z' with a wireless signal icon above it, connected to a desktop computer labeled 'PC'.</p>	<p>You can view results on the PC through a wireless connection using Zigbee Gateway technology.</p>
<p>FORA D15f</p>  <p>The diagram shows a handheld device labeled 'FORA D15f' connected to a 'FORA Telehealth System' via a cable. The system includes a 'FORA Telehealth System' box and an 'INTERNET' cloud icon.</p>	<p>Data Transmission Function: You can transmit results to the Telehealth System through cable connection using Telehealth Gateway technology.</p>

Contents of the System


- Your system includes:
- 1 Blood glucose plus blood pressure monitor
 - 2 Pressure cuff
 - 3 Quick start user guide
 - 4 Owner's manual
 - 5 Daily logbook
 - 6 Protective Storage Case
 - 7 Warranty card
 - 8 4 x 1.5V AA alkaline batteries

These products have been designed, tested, and proven to work together as a system to produce accurate blood glucose test results. Use only FORA D15 test strips and FORA control solution with your FORA D15 Monitor.

NOTE

- If your system has been opened prior to use and/or does not contain all parts listed above, please return it to the place of purchase.
- Please note that RS232 Cable for model FORA D15c, Bluetooth Adapter for FORA D15b, Zigbee Gateway Adapter for FORA D15z, Telehealth Gateway for FORA D15f, test strips, control solutions and lancets are not included in the standard kits. Please contact Toll Free Customer Care at 1-888-307-8188, 1-866-469-2632 or place of purchase for availability.

Appearance and Key Function of the Monitor

1 RF Indicator

lights up red or green to assist in data transmission.

Only for FORA D15b & FORA D15z

2 Memory Button

is used to enter the system memory and adjust settings.

3.1 Button

is a power button for blood pressure measurement.

3.2 RF Button

assists you in wireless data transmission.

Only for FORA D15b & FORA D15z

4 Set Button

is used to set up monitor.

5 LCD Screen

displays test results are displayed with symbols and simple messages.

6 Test Slot

is where you insert the test strip.

7 Air Tube

8 Arm Cuff


9 Air Plug

10 Air Plug

(Connect to the air jack)

11 Air Jack

Monitor Display

1 Date & Time

2 Test Result Area

Test results are displayed here. (A decimal point appears when the unit of measurement is millimoles per litre.)

3 Blood Drop Symbol

Flashes when it is ready to apply the sample.

4 Ketone Warning

Appears when glucose result is above 240mg/dL (13.3mmol/L).

5 Test Strip Message

Appears when the monitor is ON or the test strip has not been inserted.

6 Memory

Appears when you review the memory.

7 Pulse Symbol

8 Pulse Rate

9 Day Average

Indicates that the displayed test result is an average.

10 Temperature Message

Appears with the ambient temperature after test strip insertion.

11 Low Battery

Appears when the battery is low or must be replaced.

12 Diastolic Pressure Symbol

13 Units of Measurement

mmol/L or mg/dL (mmHG or Kpa) will appear with the test result.

14 Systolic Pressure Symbol

15 Code

Appears with the code number of the test strips currently in use.

Test Strip

This system measures the amount of sugar (glucose) in whole blood. Blood is applied to the absorbent hole of the test strip and is automatically drawn into the reaction cell where the reaction takes place.

The test strip consists of the following parts:

Contact Bars

Insert this end of the test strip into the monitor.
Push it in firmly until it will go no further.

Test Strip Handle

Hold this part to insert the test strip into the slot.

Confirmation Window

This is where you confirm if enough blood has been applied to the absorbent hole of the strip.

Absorbent Hole

Apply a drop of blood here and the blood will be absorbed automatically.

Please see section,
“Test Procedure”, for complete instructions.


NOTE

Always hold the test strip by the Test Strip Handle when inserting it into the monitor. Be careful not to touch the Absorbent Hole as it may contaminate the sample.

BEFORE USE

Battery Replacement

Your meter comes with four 1.5V AA alkaline batteries. The meter will alert you when the power is low by displaying two different messages:


1. With battery exhausted symbol “” appears on the display: the meter is functional and the result remains accurate for another 30 measurements, but battery should be replaced.


2. With battery exhausted symbol “” flashes on the display: the batteries can not provide enough power for a test. You must replace the batteries immediately.

To replace the battery, make sure that the meter is turned off.

Step1

Press the buckle on the battery cover and lift it up to remove cover. Remove old batteries.

Step2

Insert four 1.5V AA size alkaline batteries. Make sure the + (positive) and – (negative) marks match in the symbols in the battery compartment. If inserted correctly, you will hear a “beep”.

NOTE


- Do not use different type, date or brand name batteries together. Use only new batteries of the required size and type.
- Replacing the batteries does not affect previous test results stored in the monitor; however, the settings may need to be updated.
- Batteries might leak chemicals if unused for a long period of time. Remove the batteries if you are not going to use the monitor for an extended period (i.e. 3 months or more).
- Batteries should be kept away from small children. If swallowed, promptly seek medical assistance.
- Discard batteries according to your local regulations.

Setting the Monitor

Your monitor comes with the time, date, memory deletion, voice volume and language preset. If you need to change the time or replace the batteries, you may have to enter the setting mode and reset. Please set the time, date, memory deletion, voice volume and language according to the following steps.


Step1

Press the S button to start setting.

Step3

Press the M button to make changes. Press the **S button** to advance to the next setting.

Step2

Year flashes first.

Step4

Date will flash. Follow Step3 for changes.

Step5

Time flashes. Follow Step3 for changes.

Time 10:08 PM

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Unit of blood pressure will flash. Follow Step3 for changes.

mmHg
KPa

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Unit of glucose concentration will flash. Follow Step3 for changes.

mg/dL
mmol/L

Step8 (For FORA D15c, D15b, D15z)

Temperature will flash. Follow Step3 for changes.

°C °F

NOTE

- The time, date, unit of measurement and unit of temperature can **ONLY be changed** in the setting mode. Therefore, these parameters cannot be changed while performing a glucose or blood pressure test.
- The monitor will turn off automatically if left idle in the Setting Mode for 3 minutes.
- You may also press the “” button at any time to turn off the monitor and exit the Setting Mode.
- The milligram per deciliter (mg/dL) is the standard unit in the United States. The millimole per liter (mmol/L) is the standard unit in Canada. Use of the wrong unit of measure may cause you to misinterpret your blood glucose level, and may lead to incorrect treatment.



BLOOD GLUCOSE MEASUREMENT

Important Information

1. Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
2. If your blood glucose results are lower or higher than usual, and you do not have symptoms of illness, repeat the test. If you have symptoms or continue to get results that are higher or lower than usual, follow the treatment advice from your healthcare professional.
3. Use only capillary whole blood sample to test your blood glucose. Using other substances will cause wrong results.
4. If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all instructions described in this owner's manual, call your healthcare professional.
5. Inaccurate results may occur in severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis.
6. Please refer to your test strip package insert for additional important information.

Coding the Monitor

When using the system or the first time, or starting to use a new vial of test strips, you need to check the code number. It is important to make sure that the code that appears on the LCD is the same as the code on the test strip vial before you proceed. To ensure accurate test results, only test blood glucose when the codes are identical.

Step1

Insert the code strip when the monitor is off. Wait until the code number appears on the LCD screen.

Make sure the code numbers on the LCD screen, code strip, and strip vial are all the same.

Step2

Remove the code strip. The LCD screen will display “OK”, indicating that the monitor has finished coding and is ready for blood glucose testing.

WARNING!

If the code number displayed on the monitor does not match the number printed on the vial, test results will be inaccurate.

Control Test

FORA control solutions contains a known amount of glucose that reacts with test strips. By comparing your control solution test results with the expected range printed on the test strip vial label, you can ensure that the monitor and the test strips are working together as a system and that you are performing the test correctly. It is very important that you do this simple check routinely to make sure your results are accurate.

When should the control solution test should be performed?

- When you use this system to test your blood for the first time, practice the procedure using control solution. When you can do three tests in a row that are within the expected range, you are ready to test your blood.
- To routinely check the meter and test strips, perform a single test for each level of control solution at least once a week.
- When you first get your Monitor.
- When you begin using a new vial of test strips.
- Whenever you suspect that the monitor or test strips are not working properly.
- When your blood glucose test results are inconsistent with how you feel, or when you think your results are inaccurate.
- When you want to practice running the test.
- If you drop the monitor.

Important Control Solution Test Information


- Use only **FORA** control solutions. Check the expiration date on the control solution vial. Do not use if expired.
- Control solution, monitor, and test strips should be at room temperature, 68 to 77°F (20 to 25°C), before testing.
- Shake the vial before use. Discard the first drop of control solution. After squeezing, wipe off the dispenser tip to avoid contamination. The above tips ensure you to get a clean sample and an accurate result.
- Use for only 3 months after the first opening date. Record the opening date on the control solution vial. Discard the remaining control solution 3 months after this date.
- Store the control solution tightly closed at temperatures 36-86°F (2-30°C). Do not freeze.

NOTE

The control solution range is printed on FORA test strip vial. Only use FORA control solution with FORA test strips. It is used to evaluate monitor and test strip performance. It is not the recommended range for your blood glucose level.

Performing a Control Solution Test

Take a test strip out with clean and dry hands first.

Step1 Insert Test Strip.

Fully insert the test strip, with the contact bars end first facing up into test slot. The meter will turn on automatically.

“CHK” and “” will appear first, followed by

- ▶ ambient temperature
- ▶ code number and flashing “”

Be sure the code number on the screen is the same as the code number on the test strip vial. If the code numbers do not match, please see “Coding the Monitor” section of this manual.


Step2 Press the M button.

After the “” symbol appears on display, press the M button so that “CTL” appears. With the “CTL” sign on the display, the meter will not store your test result in memory. If you decide not to perform a control solution test, press the M button again, and the “CTL” sign will disappear.


NOTE

- Contact bars must be inserted all the way into the monitor or you may receive an inaccurate test result.
- Every time you perform a control solution test, you must enter into the “CTL” test mode so that the test result will not be stored in the monitor memory. Failure to do so will cause confusion between the glucose test result and the control solution test result in system memory.


Step3 Obtain control solution.

Shake the control solution vial well. Remove the cap from the control solution bottle and place it on a flat surface. Squeeze the vial, discard the first drop, and wipe off the dispenser tip to prevent contamination. **Squeeze the vial again to produce another drop and place this drop on the top of cap.**

Step4 Apply Control Solution.


Hold the meter up to meet the absorbent hole of the test strip and the drop will be automatically absorbed. Once the confirmation window is completely filled, the meter will start counting down.

To avoid contaminating the control solution with the contents of the test strip, you will have to place a drop of control solution on a clean surface. Do not directly apply control solution to a test strip.

Step5 Result Appears in 7 Seconds.

After counting to 0, the control solution test result will appear. Compare this result with the range printed on the test strip vial. The result should fall within this range.

Out-of-range results


If test results fall outside the printed range, check the “Operating Problem” section in the troubleshooting guide and repeat the test. If your results continue to be out-of-range, the system may not be working properly. Do NOT test your blood. Please contact Customer Care Line at 1-888-307-8188, 1-866-469-2632 for help.

Test Procedure


Wash and dry your hands before starting. To test your blood glucose, you need: the monitor, test strips, the lancing device, and lancets.

Step1 Set the lancing device.

Pull off the cap of the lancing device.


Insert a lancet into the lancet holder and push down firmly until it is fully secured.


Twist the protective disk off the lancet.


Replace the cap by aligning the arrow on the cap with the release button.

Select the depth of penetration by turning the adjustable tip in either direction so that the arrow on the cap points to the desired depth.

Pull the cocking control back until it clicks. You will see a color change inside the release button when it is ready.

If it does not click, the device may have been cocked when the lancet was inserted.

Blood from sites other than the fingertip

A clear cap, together with the kit, make it easier to get a drop of blood for AST. When you want to obtain blood from sites other than the finger, replace the lancing device cap with the clear cap. Turn the clear cap until it is snug but not too tight, and then slide the ejection/cocking control back until it clicks.

Ensure the removable clear cap is changed every time a new patient is tested. Do not use the same end cap on different patients.

The lancing device is now ready for use. Set aside for later use.


WARNING!

To reduce the chance of infection:

- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and lancing device.

NOTE

Your lancing device may operate differently. See manufacture's instruction for your device.


Step2 Insert the test strip and check the code number.

Fully insert test strip face up with contact bars end first into the test slot. The monitor will turn on automatically and display the followings in sequence:

“CHK” and “”

▶ ambient temperature


▶ code number and flashing “”

Make sure the code number displayed on the screen is the same as the code number printed in the test strip vial. If the code numbers are different, please refer to the section, “Coding the Monitor”, for the coding procedure.

Step3 Get a drop of blood.

Select the puncture site at fingertips or at other parts of the body (AST). Clean the puncture site using cotton moistened with 70% alcohol and let it air dry.

► Fingertip testing

Hold the lancing device firmly against the side of your finger. Press the release button. You will hear a click, indicating that the puncture is complete.

► Sites other than fingertips




Please refer to the section, “Alternative Site Testing”, for possible puncture sites.

After penetration, discard the first drop of blood with clean tissue paper or cotton. Gently squeeze the punctured area to obtain blood. Be careful **NOT to smear the blood sample**. The volume of blood sample must be at least 0.7 microliters (≈ actual size).

NOTE

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- Before you decide to begin AST, please consult your health professional.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid and serum, which may affect the test result.


Step4 Apply blood into the test strip.

When “  ” flashes on the display, apply your blood to the absorbent hole of the test strip until the confirmation window is **fully covered** with blood. The meter will then begin to count down automatically.

Step5 Obtain an accurate result in 7 seconds.

The result of your blood glucose test will appear after the monitor counting down to 0. This reading will automatically be stored in the monitor memory.

Step6 Remove the lancet.

Always use caution when removing the lancet.

Take the lancet out carefully. Place the disk on a hard surface and push the exposed tip into the protective disk.

NOTE

- Do not push your punctured finger (with blood on it) against the test strip or try to apply a smeared blood sample to the test strip.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to restart the test procedure.
- The blood should completely fill the confirmation window before the meter begins to count down. If you find that the confirmation window is not filled with blood when the meter is counting. NEVER try to add more blood to the test strip. Discard the test strip and retest with a new one.
- If you have trouble filling the confirmation window, contact Customer Care Line at 1-888-307-8188, 1-866-469-2632 for help.

WARNING!

The used lancet and the used test strip may be biohazardous. Please discard them carefully under your healthcare provider's instructions.

Comparing Monitor and Laboratory Results

The monitor provides you with whole blood equivalent results. The result you obtain from your monitor may differ somewhat from your laboratory result due to normal variation. Monitor results can be affected by factors and conditions that do not affect laboratory results in the same way (See test strip package insert for typical accuracy and precision data, and for important information on Limitations). To make an accurate comparison between monitor and laboratory results, follow the guidelines below.

Before you go to the lab:

- Perform a control solution test to make sure that the monitor is working properly.
- Fast for at least eight hours before doing comparison tests.
- Take your monitor with you to the lab.

While at the lab:

Make sure that the samples for both tests, the monitor test and the lab test, are taken and tested within 15 minutes of each other.

- Wash your hands before obtaining a blood sample.
- Never use your monitor with blood that has been collected in a gray-top test tube.
- Use fresh capillary blood only.

You may still experience a variations from the result because blood glucose levels can change significantly over short periods, especially if you have recently eaten, exercised, taken medication or experienced stress.*³ For example, if you have just eaten, the blood glucose level from a finger stick can be up to 70 mg/dL(3.9 mmol/L) higher than blood drawn from a vein (venous sample) used for a lab test.*⁴ Therefore, it is best to fast for eight hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of body fluid (severe dehydration) may also cause a monitor result to differ from a laboratory result.

*³: Surwit, R.S., and Feinglos, M.N.: Diabetes Forecast (1988), April, 49-51.

*⁴: Sacks, D.B.: "Carbohydrates. " Burtis, C.A., and Ashwood, E.R.(ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.

Risk Message

Special symbols and messages appear together with your test result.


MESSAGE	WHAT IT MEANS	WHAT TO DO
Lo	Your result is below measurement limit, 20 mg/dL (1.1 mmol/L). This indicates hypoglycemia (low blood glucose.)	You should seek immediate medical assistance.
H i	Your result is above measurement limit, 600 mg/dL (33.3 mmol/L). This indicates severe hyperglycemia (high blood glucose).	You should seek immediate medical assistance.
KETONE	Your result is equal to or higher than 240 mg/dL (13.3 mmol/L). This indicates that there is a possibility of ketone accumulation if you have Type 1 diabetes patients.	Please seek immediate medical assistance.

BLOOD PRESSURE MEASUREMENT

Suggestions Before Measuring



1. Avoid caffeine, tea, alcohol, and smoking for at least 30 minutes before measurement.
2. Wait 30 minutes after exercising or bathing before measurement.
3. Sit down for at least 10 minutes before measuring.
4. Do not measure when feeling anxious or tense.
5. Take a 5-10 minutes break between measurements. This break can be longer if necessary, depending on your physical conditions.
6. Keep the records for your physician as reference.
7. Blood pressure naturally varies from two hands. Always measure your blood pressure using the same arm.

Correct Application of the Pressure Cuff

Step1

Connect the air plug of the tubing to the air jack at the side of the monitor.

Step2

Make sure the smooth cloth is on the inside of the cuff loop and the metal D-ring does not touch your skin.

Step3

Hold your left arm in front of you with your palm facing up. Slide the cuff onto your left arm, above your elbow. The red line on the edge of the cuff should be approximately **0.8 to 1.2 inch (2 cm to 3 cm)** above your elbow. Arrange the tubing over the inside of your arm where your main arteries are located.


Step4

When the cuff is positioned correctly, pull the end of the cuff to make it snugly around your upper arm. You should be able to fit your index finger between the cuff and your arm.

Step5

Press the hook material firmly against the pile material. The top and bottom edges of the cuff should be tightened evenly around your upper arm.

Testing Your Blood Pressure

Step1

Sit down for at least 10 minutes before measuring.

Step2


Place your elbow on a table or other object. Relax your hand with the palm facing up.


Step3

Relax and make sure the cuff is at the same height as your heart. Press the “” button. Remain still and do not talk or move during the measurement.

Step4

Measurement is in progress. After the monitor turns on, all the LCD symbols will appear with a long “beep” sound. Then the cuff will begin to inflate automatically.

After the cuff pressure is reached, the cuff will begin to deflate. You will see the number decreasing and “” flashing in the display.

Step5

Read a result. After the measurement, the monitor displays the systolic pressure, diastolic pressure, and heart rate. Press “” button to turn off. Otherwise, it will turn off automatically if left idle over 3 minutes.



NOTE

- If a higher pressure value is needed, the monitor will stop deflation and inflate again.
- If you press the “” button during the measurement, the monitor will turn off immediately.

USING THE MONITOR MEMORY

Your monitor stores 352 most recent blood glucose test results along with their respective dates and times in its memory. It also provides you with 7-, 14-, 21-, 28-, 60- and 90-day averages of your blood glucose test results. You can review these test result averages in memory through the following easy steps.

Viewing Results on the Monitor


Step 1 Read the Average of Blood Glucose Results:

When the monitor is off, press and release the M button. The 7-day average will appear, indicating that you are in the memory mode. If you continue to press the M button, the 14-, 21-, 28-, 60-, and 90-day average will appear in order.

When using the monitor for the first time, or when there are no test results in memory, you will only see flashing “M”. This means that there are no test results in memory.

The 7-day average is calculated from the blood glucose results obtained during the last 7 days. The example displayed means that the blood glucose average performed over the last 7 days was 128 mg/dL.

Step2 Recall the Stored Test Results.

After 90-day average, the other results will be displayed in the order of time and date taken. You will see the number of the result stored, followed by the result. When the memory is full, the oldest results are deleted as new ones are added.

Step3 Exit the Memory Mode.


You can also exit memory mode by pressing the “” button to turn off the monitor.

NOTE

- The control solution results are **NOT** stored in the memory. The list of past results and the result averages are for blood glucose results only.
- If no button is pressed for 3 minutes, the meter will display “OFF” and turn off automatically.

Deleting Memory

To clear the memory, turn off the monitor. Press and hold the M button for 3 seconds. "CLr/ALL" then appears to indicate that **ALL** the memory stored in the monitor will be deleted.


VIEWING RESULTS ON A PERSONAL COMPUTER

Results in memory can be transmitted to the personal computer by either cable or wireless connection. Accessories needed to use this function are:

- Health Care System Software: downloadable from the home page of Fora Care Inc.
- Interface Cable: an accessory for model FORA D15c.
- Bluetooth Adapter: an accessory for model FORA D15b.
- Zigbee Gateway: an accessory for model FORA D15z.
- Telehealth Gateway: an accessory for model FORA D15f.

To learn more about Health Care Software System or to obtain an interface cable separately, please contact Customer Care Line at 1-888-307-8188, 1-866-469-2632 for help.

Data transmission via cable (FORA D15c)

Step1 Installing the Software

Install Health Care System Software on your computer by following the instructions provided on

Fora Care's website: www.foracare.com/usa

Step2 Connecting to a Personal Computer


Connect the interface cable to a serial port of your computer. With the meter turned off, connect the cable to the data port of the meter. Then "PCL" and the time will appear on the monitor screen to indicate that the meter is ready to transmit data.

Step3 Data transmission

Follow the instructions to transmit data from the monitor's memory to your computer. Results transmitted will include date and time. Remove the cable and the monitor will automatically turn off.

Data transmission via bluetooth adapter (FORA D15b)

Please follow the steps below before transmitting data to your computer via a Bluetooth adapter.

Step1 Setting FORA D15b monitor for Bluetooth pairing:

Make sure your FORA D15b monitor is off, then hold down the RF button until “PCL” and time appear. The Bluetooth indicator on the monitor will blink quickly. This indicates that the monitor is now ready for Bluetooth wireless pairing.

Step2 Pairing FORA D15b monitor with your computer:

Select the FORA D15b monitor device code, “TaiDoc-BTM.” Enter the Bluetooth passkey (PIN), 111111, for pairing. If pairing succeeds, the Bluetooth indicator on the monitor will blink slowly.

Step3 Connecting Bluetooth Serial Port Service:

Right click on icon “TaiDoc-BTM.” Select “connect” and “Bluetooth Serial Port Service” for connection. If the connection succeeds the Bluetooth indicator on your monitor will stay lit. Now you are ready for data transmission via a Bluetooth adapter.

NOTE

If you are using the Bluetooth Adapter provided by Fora Care Inc., please follow the installation guide “Bluetooth Easy Card” in the package.

Data transmission via zigbee gateway (FORA D15z)

Step1 Installing the Software

Install Zigbee Gateway and turn on the device. (Please contact Customer Care Line at 1-888-307-8188, 1-866-469-2632 to set up the device).

Step2 Data transmission


- After you finish measuring, your results will be transmitted automatically.
- Press the RF button to transmit the latest measurement data.

The RF indicator on the monitor will light up and you will hear two beeps from the gateway. It means your result has been successfully transmitted to the server. The gateway will “beep” twice again and then the monitor will turn off.

NOTE

- While the monitor is connected to the PC or the server, it is unable to perform a blood glucose or blood pressure test.
- When the Zigbee Gateway device is normal, the RF indicator will be green. If it is red, please contact Customer Care Line at 1-888-307-8188, 1-866-469-2632 for help.
- If you hear a long “beep” from the gateway, an error has occurred. Please try the connection again. If the problem persists, please contact Customer Care service for assistance.

Data transmission via telehealth gateway (FORA D15f)

FORA D15f Blood Glucose plus Blood Pressure Monitoring System can transmit results remotely via FORA Telehealth Gateway.

NOTE

- Telehealth serves as a tool to provide delivery of health-related information to the hospital and enables the health care professionals to monitor the health condition of patients. Please check with your healthcare provider about the availability of this service.
- FORA Telehealth Gateway (9014 series) can be purchased separately.

FORA D15 Monitoring System	To use with FORA Telehealth Gateway
FORA D15f	GW 9014A (RS232)
FORA D15b	GW 9014B (Bluetooth)

Before transmitting results, turn on the power of FORA Telehealth Gateway.

Method one:

- ▶ Taking measurements while interface cable is connected.

Connect both ends of the RS232 interface cable to both of the data port of FORA D15f and RS232 port of FORA Telehealth Gateway. Take measurements according to the sections of testing procedure in this Owner's Manual.


After measurement is finished the monitor will automatically switch to the transmission mode within one minute. After transmission is completed, the monitor will automatically turn off.

Method two:

- Use the interface cable to transmit results after measuring.

When the meter is off, connect both ends of the RS232 interface cable to the data port of FORA D15f and RS232 port of FORA Telehealth respectively. Press S button, PCL will appear on the display to indicate the monitor is switched to the transmission mode. After transmission is completed, the monitor will automatically turn off.

Step1

Step2

NOTE

The S button of FORA D15f has two functions. When you press the S button for one second, PCL will appear on the display and the monitor is switched to the transmission mode. When you press and hold the S button for 5 seconds, the monitor will switch to the setting mode. In the setting mode, only date and time can be changed.

TAKING CARE OF YOUR MONITOR AND TEST STRIPS

To protect the monitor and test strips from dirt, dust and other contaminants, please wash and dry your hands thoroughly before use.

Cleaning

1. To clean the monitor exterior, wipe using a cloth moistened with tap water or a mild cleaning agent. Dry, then the device with a soft and dry cloth. Do not flush with water.
2. Do not use organic solvents to clean the monitor and cuff.
3. The cuff can be cleaned by wiping with a moistened cloth and soap. Do not immerse in water.

Storage


1. Monitor Storage


- Storage condition: -4°F to 140°F (-20°C to 60°C), below 95% relative humidity.
- Always store or transport the monitor in its original storage case.
- Avoid dropping and strong impact.
- Avoid direct sunlight and humidity.
- Do not disassemble, modify or try to repair the monitor or wrist cuff by yourself.

- Do not over-wring the cuff or turn it inside out.
- If you are not going to use the monitor for an extended period of time, please remove the batteries.

2. Strip Storage


- Storage condition: 39°F to 104°F (4°C to 40°C), below 85% relative humidity. **Do NOT freeze.**
- Store your test strips in their original vial only. Do not transfer to another container.

- Store test strip packages in a cool and dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately replace the vial cap and close it tightly.


- Handle the test strip with clean, dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the vial label when you first open it. Discard the vial along with remaining test strips 3 months after this date.

- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children as the cap and the test strip may be a choking hazards. If swallowed, promptly seek medical assistance.

3. Control Solution Storage

- Storage condition: Store the control solution tightly closed at temperatures 36 to 86°F (2 to 30°C). **Do NOT freeze.**
- Record the opening date on the control solution vial. Discard the remaining control solution 3 months after this date.

TROUBLESHOOTING GUIDE

The following is a summary of various display messages. These messages help to identify certain problems but do not always appear in all cases when a problem has occurred. Improper use may cause an inaccurate result without producing an error message. In the event of a problem, refer to the information under “WHAT TO DO.”

Never try to disassemble the monitor under any circumstance. If you encounter any error messages not listed below or if you have followed the actions recommended below but the problem persists, please call our Customer Care Line at 1-888-307-8188, 1-866-469-2632 for support.

Error Message

MESSAGE	CAUSE	WHAT TO DO
Err.00	Weak pulse.	Refit cuff tightly, relax, and repeat measurement as shown in the section, "Correct Application of the Pressure Cuff".
Err.02	Monitor can not determine systolic pressure.	
Err.04	Monitor can not determine systolic or distolic pressure or it is out of range.	
Err.07 Err.08		
Err.22	Invalid coding.	Try a new code strip and insert again.
Err.24	Used strip inserted.	Use a new strip.
Err.25	Ambient temperature is below operating temperature.	Please operate this monitor between 50°F and 104°F (10°C and 40°C).
Err.26	Ambient temperature is over operating temperature.	
Err.05 Err.06 Err.09 Err.10 Err.11 Err.15 Err.20 Err.28	Those messages indicate internal problems with the monitor. Do not attempt to fix them. Please contact Customer Care Line at 1-888-307-8188, 1-866-469-2632.	

Operating Problem

1. Blood Glucose Measurement

1. If the monitor does not display a message after you have inserted a test strip.

PROBABLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Batteries incorrectly installed or absent.	Check that the batteries are correctly installed.
Test strip inserted upside down or incompletely.	Insert the test strip correctly with the contact bars end first and facing up.
Defective monitor.	Contact Customer Care Line at 1-888-307-8188, 1-866-469-2632 for help.

2. You turned off the monitor but it makes a beeping sound.

PROBABLE CAUSE	WHAT TO DO
The strip is still inside the slot.	Remove the strip. If you need to test your blood glucose, insert an unused strip.

3. The test does not start after applying the sample.

PROBABLE CAUSE	WHAT TO DO
Insufficient blood sample.	Repeat the test with a new test strip and a larger sample.
Defective test strip.	Repeat the test with a new test strip.
“  ” was not flashing on the display when sample was applied.	Repeat the test with a new test strip. Apply sample only when “  ” appears on the display.
Defective monitor.	Contact Customer Care Line at 1-888-307-8188, 1-866-469-2632 for help.

4. The control solution test result is out of range.

PROBABLE CAUSE	WHAT TO DO
Error in performing the test.	Read the instructions thoroughly and repeat the test again.
Improper code number.	Check if the code number on the display matches the code number on the test strip vial.
Control solution vial shaken poorly.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiration date or the opening date of the control solution.
Control solution that is too warm or too cold.	Allow control solution, monitor and test strips to be at room temperature (68-77°F/20-25°C) before testing.
Test strip deterioration.	Repeat the test with a new test strip.
Monitor malfunction.	Contact Customer Care Line at 1-888-307-8188, 1-866-469-2632 for help.

2. Blood Pressure Measurement

5. No display after pushing the “” button.

PROBABLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Batteries incorrectly installed or absent.	Check that the batteries are correctly installed.

6. Heart rate is higher/lower than user’s average.

PROBABLE CAUSE	WHAT TO DO
Moving during measurement.	Repeat measurement.
Measuring right after exercise.	Rest at least 30 minutes before measurement.

7. A result is higher/lower than user's average measurement.

PROBABLE CAUSE	WHAT TO DO
May not have been in correct position while measuring.	Adjust to a correct position to measure.
Natural variation of blood pressure from time to time.	Keep in mind for next measurement.

8. Cuff inflates again while measuring.

PROBABLE CAUSE	WHAT TO DO
Cuff is not fastened.	Fasten the cuff again.
Normal action. If user's blood pressure is higher than the previous value, the monitor would automatically pump to a more suitable pressure. Relax and wait for the measurement.	

SPECIFICATIONS

Power Source: 4x AA (1.5V) alkaline batteries

Size of Monitor w/o Cuff: 137mm (L) x 90 mm (W) x 54 mm (H)

Memory: 352 measurement results with respective dates & times

Power Saving: Automatic power off if left idle for 3 minutes

Data Transmission: RS232/Bluetooth/Zigbee/Telehealth Gateway

System Operating Conditions:

50°F to 104°F (10°C to 40°C), below 85% R.H.

Monitor Storage Conditions:

-4°F to 140°F (-20°C to 60°C), below 95% R.H.

Strip Storage Conditions:

39°F to 104°F (4°C to 40°C), below 85% R.H.

KETONE Warning: glucose value over 240 mg/dL

Measurement Unit: mg/dL or mmol/L

Linear Range: 20-600 mg/dL (1.1 -33.3 mmol/L)

Precision: ±5% (CV)

Accuracy: ±15mg/dL when glucose < 75mg/dL;

±20% when glucose ≥ 75mg/dL

Pressure Range: 0-300 mmHg

Heart Rate Range: 40-199 beats per minute

Measurement Unit: mmHg or KPa

Accuracy of Pressure: ±3mmHg or ±2% of reading

Accuracy of Heart Rate: ±4% of reading

Maximum Inflation Pressure: 300 mmHg

Bluetooth Carrier Frequency: 2400MHz to 2483.5MHz
Bluetooth Modulation Method: GFSK, 1Mbps, 0.5BT Gaussian
Transmission Power: +3dBm to -20dBm; Power control 4 stage
Receiving Signal Range: -88dBm to -20 dBm
Receiver IF Frequency: 1.5MHz center frequency
Maximum Data Rate: Asynchronous:723.2kbps/57.6kbps;
Synchronous: 433.9kbps/433.9kbps

Telehealth Gateway Specification

Model no.: TD-9014A(GW 9014A) /TD-9014B(GW 9014B)
Input interface: RS232(TD-9014A) / Bluetooth(TD-9014B)
Output interface: ISP dial-up service or Ethernet
Dimension: 170mm x 120mm x 30mm
Weight: 221.8 g
Power AC power adapter: Input: 100-240V AC; Output: 6V, 1A
Operation Environment: 10°C to 40°C, below 90% R.H. (Non-condensing)
Storage Environment: -4°F to 149°F(-20°C to 65°C), below 90% R.H. (Non-condensing)

This device has been tested to meet the electrical and safety requirements of:
IEC 60601-1, EN 60601-1, IEC 61010-1, EN 61010-1, EN 61010-2-101,
EN 60601-1-2, EN 61326

FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT

15.21

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

15.105(b)

Federal Communications Commission (FCC) Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference with radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference with radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a different circuit than that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation of the device.

FCC RF Radiation Exposure Statement:

1. This Transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
2. This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 centimeters between the radiator and your body.

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