Dear **FORA D30** System Owner:

Thank you for purchasing the **FORA D30** Blood Glucose plus Blood Pressure Monitoring System. This manual provides important information to help you operate this system correctly. Before using this product, please read the following contents thoroughly and carefully.

According to the clinical studies from American Diabetes Association, elevated blood pressure often accompanies adult diabetes. These studies also suggest that diabetes patients could reduce the risk of cardiovascular diseases by managing their blood glucose levels and blood pressure. Therefore, with the compact size and easy operation of this **FORA D30** Blood Glucose plus Blood Pressure Monitoring System, you can easily monitor your blood glucose levels and blood pressure by yourself at any time or place. In addition, this system can help you and your healthcare professionals to monitor and adjust your treatment plans, and keep your diabetes and blood pressure under control.

If you have other questions regarding this product, please contact the place of purchase or call Customer Care Line.

IMPORTANT SAFETY INSTRUCTIONS

READ BEFORE USE

1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or damaged.
4. Do **NOT** use the device where aerosol sprays are being used or where oxygen is being administered.
5. Do **NOT** use under any circumstances on neonates or infants.
6. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only.
7. Before using this device to test blood glucose, read all the instructions thoroughly and practice the test. Do all quality control checks as directed.

8. Keep the device and testing supplies away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
9. Keep the equipment away from hot surfaces.
10. Do **NOT** apply the cuff to areas other than the place directed.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

TABLE OF CONTENTS

BEFORE YOU BEGIN	06
Important Information	06
Intended Use	07
Test Principle	07
Contents of System	08
Monitor Overview	09
Screen Display	10
Test Strip	11
Speaking Function	12
SETTING THE MONITOR	17
THE FOUR MEASURING MODES	22
BEFORE TESTING	23
Control Solution Testing	23
Performing a Control Solution Test	24
TESTING WITH BLOOD SAMPLE	27
Overview of the Lancing Device	27
Setting up the Lancing Device	28
Preparing the Puncture Site	30
Alternative Site Testing	32
Performing a Blood Glucose Test	34
TESTING YOUR BLOOD PRESSURE	37
Before Measurement	37
Fitting the Cuff Properly	38
Proper Measurement Position	39
Taking Measurements	41

MONITOR MEMORY	44
Reviewing Test Results	44
Reviewing Day Average Results	45
Downloading Results onto a Computer	47
MAINTENANCE	50
Battery	50
Using the AC Adapter	52
Caring for Your Monitor	53
Caring for Your Test Strips	54
Important Control Solution Information	55
SYSTEM TROUBLESHOOTING	56
Result Readings	56
Error Messages	57
Troubleshooting	59
DETAILED INFORMATION	62
Reference Values	62
Comparing Monitor and Laboratory Results	64
SPECIFICATIONS	66

BEFORE YOU BEGIN

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have any symptoms of illness, repeat the test. If you have symptoms or continue to get results higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood samples to test your blood glucose. Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all the instructions described in this owner's manual, call your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hyperglycaemic-hyperosmolar state, with or without ketosis. Please consult the healthcare professional before use.

Intended Use

This system is a 2 in 1 system designed to measure blood glucose outside the human body (*in vitro* diagnostic use) and to measure blood pressure non-invasively. It is intended for use in the home and in clinical settings as an aid to monitoring the effectiveness of diabetes control. It should not be used for the diagnosis of diabetes and hypertension, or for testing neonates.

The alternative site testing in this system can be used only during steady-state blood glucose conditions.

This monitor has some speaking function but is not intended for use by the visually impaired.

Test Principle


Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The monitor measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

The blood pressure is measured non-invasively at the arm, based on the oscillometric method.

This device is NOT able to take measurements in the presence of common arrhythmia, such as arterial or ventricular premature beats or arterial fibrillation. It may produce reading errors.

Contents of System

Your new **FORA D30** system kit includes:


- | | |
|--------------------------|----------------------------------|
| 1 Monitor | 6 4 x 1.5V AA alkaline batteries |
| 2 Protective Wallet | 7 Daily Log Book |
| 3 Quick Start User Guide | 8 Pressure Cuff (arm type) |
| 4 Owner's manual | |
| 5 Warranty Card | |

Test strips, control solutions, sterile lancets or lancing device are not included in the kit. They can be purchased separately. Please make sure you have those items needed for a blood glucose test beforehand.

NOTE


If any items are missing from your kit or opened before use, please contact your local customer service or place of purchase for assistance.

Monitor Overview

- | | |
|------------------------------------|--|
| 1 Test Slot | 10 Earphone Jack |
| 2 Display Screen | 11 DC Adapter Port |
| 3 Up Button | 12 Air Jack |
| 4 Down Button | 13 Air Tube |
| 5 On/Off Button | 14 Air Plug |
| 6 M Button | 15 BT Indication Light
(for FORA D30f) |
| 7 S Button | 16 Test Strip Ejector |
| 8 Repeat Button | 17 Battery Compartment |
| 9 Data Port (for FORA D30e) | |

Screen Display



- | | |
|--------------------------------|---|
| 1 Date and Time | 11 Unit of Pressure Measurement |
| 2 Blood Drop Symbol | 12 Memory Mode Symbol |
| 3 Test Strip Symbol | 13 Day Average |
| 4 Volume Symbol | 14 Pulse Rate |
| 5 Systolic Pressure Symbol | 15 Pulse Rate Symbol |
| 6 Unit of Glucose Measurement | 16 User Number |
| 7 Unit of Pressure Measurement | 17 Measuring Modes |
| 8 Diastolic Pressure Symbol | 18 Diastolic Pressure/Blood Glucose Value |
| 9 Battery Symbol | 19 Systolic Pressure Value |
| 10 Ketone Warning | |

Test Strip

Front Back

The front side of the test strip should face up when inserting the test strip.

1 Absorbent Hole

Apply a drop of blood here. The blood will be automatically absorbed.

2 Confirmation Window

This is where you confirm if enough blood has been applied to the absorbent hole in the strip.

3 Test Strip Handle

Hold this part to insert the test strip into the slot.

4 Contact Bars

Insert this end of the test strip into the meter. Push it in firmly until it will go no further.

ATTENTION:

Test results might be wrong if the contact bar is not **fully** inserted into the test slot.

NOTE

The **FORA D30** monitor should only be used with **FORA D30** Test Strips. Using other test strips with this meter can produce inaccurate results.

Speaking Function

FORA D30 “speaks” aloud with step by step instructions to guide you through the process of blood glucose testing. The following table tells you when and what the monitor “says”.

WHEN does the monitor say?	WHAT does the monitor say?
Setting up the monitor	
Setting the year	You are now in setting mode. The Year is (number).
Setting the month	The month is (month).
Setting the date	The date is the (number)th.
Setting the hour mode	The time is (12/24) hour mode.
Setting the hour	The hour is (number) PM/AM.
Setting the minutes	The minutes are (number).
Setting the blood glucose unit	Blood glucose unit is (milligrams per deciliter/ millimoles per liter.)
Setting the temperature unit	Temperature unit is degree (Fahrenheit/ Celsius).
Setting the blood pressure unit	Blood pressure unit is (millimeters of mercury / kilopascals).
Setting the voice volume	Voice volume is (number).
Deleting the memory	1. Delete memory record function. Press Memory key to clear all records. 2. Are you sure you want to clear ALL records?

	3. Press Memory key again to clear all records. 4. All memory records have been deleted.
Selecting language	Thank you for using FORA Tele-Health products.
At the end of monitor setting	1. Today's date is (day) (month) (day)th, (year). 2. The time is (hour minute) AM/ PM. 3. Blood glucose unit is (milligrams per deciliter/ millimoles per liter). 4. Temperature unit is degree (Fahrenheit/ Celsius). 5. Blood pressure unit is (millimeters of mercury / kilopascals) 6. Voice volume is (number) 7. Have a nice day.
Using the monitor	
Insert test strip When the monitor is ready to test. (symbol appears on display)	1. Thank you for using FORA Tele-Health products. Test strip inserted. 2. Today's date is (week day) (month) (day), (year). 3. The time is (hour minute) PM/ AM.

	4. The room temperature is (number) degree (Fahrenheit/ Celsius). 5. The user's number is one. 6. Please apply blood onto the test strip.
When selected AC for testing	This is before meal measurement.
When selected PC for testing	This is after meal measurement.
When selected QC for control solution testing	You are now in control solution mode.
After sample is applied onto the test strip	Now testing.
When the test is completed. (result appears on display)	The user's number is one. Your blood glucose is (number) (milligrams per deciliter/ millimoles per liter).
After displaying the results. (for FORA D30e only)	Connecting to your computer.
When the test result is < 20 mg/dL.	Your blood glucose is less than 20 milligrams per deciliter.
When the test result is > 600 mg/dL.	Your blood glucose is over 600 milligrams per deciliter.
When turning on the monitor	Thank you for using FORA Tele-Health products. The user's number is (number). Please select user's number.
After selecting the user number.	The user's number is (number). Please relax during measurement.

When the test is completed. (result appears on display)	Your systolic blood pressure is (number) (millimeters of mercury/ kilopascals). Your diastolic blood pressure is (number) (millimeters of mercury/ kilopascals). Your heart rate is (number) beats per minute.
Recall the Stored Test Results	1. Thank you for using FORA Tele-Health products. 2. Today's date is (day) (month)th, (year). 3. The time is (hour minute) AM/ PM. 4. You are now in memory mode. 5. The user's number is (number). 6. Please select user's number. 7. The user's number is (number). 8. Memory Record for (month) (day) th, (year), (hour minute) AM/PM. 9. Your blood glucose was (number) (milligrams per deciliter/ millimoles per liter). OR Your systolic blood pressure was (number) millimeters of mercury. Your diastolic blood pressure was (number) millimeters of mercury. Your heart rate was (number) beats per minute.

Reading the Average of Blood Glucose Results:	1. You are now in day average mode. 2. The (number) day blood glucose average is (number) (milligrams per deciliter/ millimoles per liter) for (number) records.
If no results stored in memory	No memory record available.
If no results stored in memory	The (number) day blood glucose average does not exist.
When the monitor is turned off	Have a nice day.
When room temperature is outside operating range, which is 50-104°F (10-40°C).	Room temperature is out of range, thus unable to measure.
When the power is insufficient	Battery is dead. Please replace.
Removing the test strip during the measurement.	Please ensure proper operation then try again.
When inserting a used test strip.	Test strip has been used. Please replace.

NOTE


If you would like to listen to the voice messages by earphone, you can plug a pair of earphones into the earphone port on the side of your monitor.

SETTING THE MONITOR

Before using your monitor for the first time, check and update these settings. Make sure you complete the steps below and save the desired settings.

► Entering the setting mode

Start with the monitor switched off. Keep pressing **S** to turn on the monitor.

Step 1

Set the date.


With the year flashing, press **1/7** to select the correct year. Press **S** to set.

With the month flashing, press **1/7** to select the correct month. Press **S** to set.


With the day flashing, press **▲▼** to select the correct day. Press **S** to set.

Step 2

Set the time format.



Press **▲▼** to select the desired time format — 12h or 24h. Press **S** to set.

Step 3

Set the time.



With the hour flashing, press **▲▼** to select the correct hour. Press **S** to set.

Step 4

Set the unit of glucose measurement.



Press  to select mg/dL or mmol/L. Press  to set.

Step 5

Set the unit of temperature.



Press  to select °F or °C. Press  to set.

Step 6

Set the unit of pressure measurement.

Press  to select mmHg or KPa. Press  to set.


Step 7

Set the speaking volume.

There are seven (7) speaking volume options for you to choose. Press **AV** to select the desired speaking volume appears. Press **S** to set.

Volume 0 indicates that the speaking function is turned off, and will not display “”. Volume 1 to 7 indicates speaking volume from low to high, and “” will be displayed during testing.


Step 8

Delete the memory.

While “**del**” and a blinking “**M**” symbol appear on the display, if you do not wish to delete the saved results, press **S** to skip this step.

If you wish to delete **ALL** the results, press **M** twice. “**CLr**” and “**M**” are displayed on the monitor, which indicates that all results have been deleted.

Step 9

Choosing a language.

Press  to select L1/L2. The default language for the monitor is L1, which is English. Press .

Then the monitor will turn off automatically or you can press  to turn off the monitor. The monitor repeats all your settings before shutdown.

Congratulations! You have completed all settings!

NOTE

- These parameters can **ONLY** be changed in the setting mode.
- If the monitor is idle for 3 minutes during the setting mode, it will turn off automatically.
- The milligram per deciliter (mg/dL) is the standard unit in the United States. The millimole per liter (mmol/L) is the standard unit in Canada. Use of the wrong unit of measure may cause you to misinterpret your blood glucose level, and may lead to incorrect treatment.

THE FOUR MEASURING MODES


The monitor provides you with four modes for measuring: General, AC, PC and QC.

Modes	Use when
General (do not display)	any time of day without regard to time since last meal
AC	no food intake for at least 8 hours
PC	2 hours after a meal
QC	testing with the control solution

You can switch between each mode by:


Step 1

Start when the monitor is off. Insert a test strip to turn on the monitor, the screen will display a flashing blood drop.

Step 2

Press **M** to switch between General, AC, PC and QC mode.

BEFORE TESTING

Control Solution Testing

FORA Control Solution contains a known amount of glucose that reacts with test strips and is used to ensure your monitor and test strips are working properly together.

Do a control solution test when:


- you first receive the monitor
- at least once a week to routinely check the monitor and test strips
- you begin using a new vial of test strips
- you suspect the monitor or test strips are not working properly
- your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate
- testing process or
- you dropped or think you have damaged the monitor.

Performing a Control Solution Test

Step 1

Insert the test strip to turn on the monitor.


Insert the test strip into the monitor. Wait for the monitor to display the test strip and blood drop symbols.

Step 2

Press **M** to mark this test as a control solution test.

If you press **M** again, the “QC” will disappear and this test is no longer a control solution test.

Step 3


Applying test control solution.

Shake the control solution vial thoroughly before use.

Squeeze out a drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap.

Hold the monitor and move the absorbent hole of test strip to touch the drop. Once the confirmation window fills completely, the monitor will begin counting down.


To avoid contaminating the control solution, do not apply control solution onto a strip directly.

Step 4

Read and compare the result.

After counting down to 0, the control solution test result will appear on the display. Compare this result with the range printed on the test strip vial and it should fall within this range. If not, please read the instructions again and repeat the control solution test.

Out-of-range results.


If you continue to have test results fall outside the range printed on the test strip vial, the monitor and strips may not be working properly. Do NOT test your blood. Contact the local customer service or place of purchase for help.

NOTE

- The control solution range printed on the test strip vial is for control solution use only. It is not a recommended range for your blood glucose level.
- See the **Maintenance** section for important information about your control solutions.

TESTING WITH BLOOD SAMPLE

Overview of the Lancing Device

- | | | |
|--------------------------|--------------------------|----------------------|
| 1 Cap | 5 Lancet Holder | 9 A Clear Cap |
| 2 Adjustable Tip | 6 Ejector | |
| 3 Depth Indicator | 7 Release Button | |
| 4 Arrow | 8 Cocking Control | |

If your lancing device differs from the one shown here, please refer to the manufacturer's manual to ensure proper usage.

NOTE


To reduce the chance of infection:

- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.

Setting up the Lancing Device


Step 1

Pull off the cap of the lancing device.


Step 2

Insert a lancet into the lancet holder and push down firmly until it is fully secured.


Step 3

Twist the protective disk off the lancet.


Step 4

Replace the cap by aligning the arrow on the cap with the release button.

Step 5


Select the depth of penetration by turning the adjustable tip in either direction so that the arrow on the cap points to the desired depth.

Step 6

Pull the cocking control back until it clicks. You will see a color change inside the release button when it is ready.

If it does not click, the device may have been cocked when the lancet was inserted.

Hold your fingers on the body. Not the cap.

The lancing device is now ready for use. Set it aside for later use.

Preparing the Puncture Site

Stimulating blood perfusion by rubbing the puncture site before blood extraction has a significant influence on the glucose value obtained.

Blood from a site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site was rubbed prior to blood extraction, the difference was significantly reduced.

Please follow the suggestions below before obtaining a drop of blood:

- **Wash and dry your hands before starting.**
- Select the puncture site either at fingertips or other body parts (please see section “Alternative Site Testing” (AST) on how to select the appropriate sites).
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry.**
- Rub the puncture site for about 20 seconds before penetration.
- Use a clear cap (included in the kit) while setting up the lancing device.


► Fingertip Testing

Press the lancing device tip firmly against the side of your fingertip. Press the release button to prick your finger, then a click indicates that the puncture is complete.

► Blood from Sites Other Than the Fingertip

Replace the lancing device cap with the clear cap for Alternative Site Testing (AST). Pull the cocking control back until it clicks. When lancing the forearm, upper arm, hand, thigh, or calf, avoid lancing the areas with obvious veins because of excessive bleeding.

NOTE


- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- Please consult your health care professional before you begin AST.
- It is recommended to discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Alternative Site Testing

IMPORTANT: There are limitations with AST (Alternative Site Testing). Please consult your health care professional before you perform AST.

What is AST?

Alternative site testing (AST) means that people use parts of the body other than the fingertips to check their blood glucose levels. This system allows you to test on the palm, the forearm, the upper arm, the calf or the thigh with the equivalent results to fingertip testing.

What is the advantage?

Fingertips feel pain more readily because they are full of nerve endings (receptors). At other body sites, since nerve endings are not so condensed, you will not feel as much pain as at the fingertips.

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at fingers reflects these changes faster than capillary blood at other sites. Thus, when testing blood glucose during or immediately after a meal, exercise, or any other event, **take the blood sample from your finger only.**

We strongly recommend that you perform AST **ONLY** at the following times:

- In a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do **NOT** use AST if:


- You think your blood glucose is low.
- You are unaware of hypoglycemia
- You are testing for hyperglycemia
- Your AST results do not match the way you feel.
- Your routine glucose results are often fluctuating.

Performing a Blood Glucose Test

Step 1

Insert the test strip to turn on the monitor.


Wait for the monitor to display the strip symbol and blood drop symbol.

Step 2

Obtain a blood sample.

Use the previously set lancing device to puncture your desired site. After penetration, discard the first drop of blood with a clean tissue or cotton. Gently squeeze the punctured area to obtain another drop of blood. Be careful NOT to smear the blood sample.


The volume of blood sample must be at least 0.7 microliters in volume. (* actual size).

Step 3

Apply the sample.

Hold the blood drop to the absorbent hole of the test strip. Blood will be drawn in and after the confirmation window is completely filled, the monitor begins counting down.


NOTE

- Do not press your punctured site against the test strip or try to smear the blood.
- If you do not apply a blood sample to the test strip within 3 minutes, the monitor will automatically turn off. You must remove and reinsert the test strip to start a new test.
- The confirmation window should be filled with blood before the monitor begins to count down. NEVER try to add more blood to the test strip after your drop of blood has moved away. Discard the used test strip and retest with a new one.
- If you have trouble filling the confirmation window, please contact your health care professional or the local customer service for assistance.

Step 4

Read your result.

The result of your blood glucose test will appear after the monitor counts down to 0. This blood glucose result will automatically be stored in the memory.

For **FORA D30f**, the monitor will turn on the Bluetooth connection and connect to your personal computer. If you do not want to transmit the data, press **⏻** to turn off the monitor. If you want to transmit the test result, follow the instructions provided in the software to transmit data.

Step 5

Remove the used test strip and the lancet.

Remove the test strip by hand, throw the strip into a disposal container. The monitor will turn itself off automatically after the test strip is removed.

Always follow the instructions in the lancing device insert when removing the lancet.

WARNING

The used lancet and test strip may be potentially biohazardous. Please discard them carefully according to local regulations.

TESTING YOUR BLOOD PRESSURE

This monitor provides you with two different ways to measure your blood pressure. Select from the operations below:

- **Single measurement:**
Perform an individual blood pressure measurement.
- **Averaging measurement:**
Automatic perform three (3) consecutive blood pressure measurements and display the average result in the end.


Before Measurement

- Avoid caffeine, tea, alcohol, and tobacco at least 30 minutes before measurement.
- Wait 30 minutes after exercising or bathing before measurement.
- Sit or lie down for at least 10 minutes before measurement.
- Do not measure when feeling anxious or tense.
- Take a 5-10 minute break between measurements. This break can be longer if necessary, depending on your physical condition.
- Keep the records for your doctor as reference.
- Blood pressure varies between each arm. Always measure your blood pressure on the same arm.
- To take a blood pressure measurement after performing a blood glucose test, make sure that the test strip has been removed from the monitor.

► Fitting the Cuff Correctly


Step 1

Connect the air plug of the tubing to the air jack of the monitor.


Step 2

Assemble the cuff as shown on the right. The smooth surface should be inside the cuff loop and the metal D-ring should not touch your skin.

Step 3

Stretch your left/right arm in front of you with your palm facing up. Slide and place the cuff onto your arm with the air tube and artery mark region (in red) toward the lower arm.


Wrap and tighten the cuff above your elbow. The red line on the edge of the cuff should be approximately **0.8 to 1.2 inches (2 to 3 cm)** above your elbow. Align the tube over the main arteries on the inside.

Step 4

Leave a little free space between the arm and the cuff, you should be able to fit 2 fingers between them. Clothing must not restrict the arm. Remove all clothing covering or constricting the measurement of the arm.

Step 5

Press the hook material firmly against the pile material. The top and bottom edges of the cuff should be tightened evenly around your upper arm.

► The range index of cuff should fall into this range.

► Proper Measurement Position

Step 1


Sit down for at least 10 minutes before measuring.

Step 2

Place your elbow on a flat surface. Relax your hand with the palm facing up.

Step 3

Make sure the cuff is about the same height as the location of your heart. Press  Remain still and do not talk or move during the measurement.

WARNING

If the cuff is relatively lower/higher than the heart, the obtained blood pressure value could be higher/lower than the actual value. A 15 cm difference in height may result in an error of around 10 mmHg.

Step 4

Measurement is in progress.

After the monitor is turned on, the cuff will begin to inflate automatically.


Taking Measurement

► Taking a single measurement

Always apply the pressure cuff before turning on the monitor.


Step 1

Press . All the LCD symbols will appear. Press  to select user “” (numbers from 1 to 4) whose test reading will be stored in its memory. Then the cuff will begin to inflate automatically.

Step 2

The heart symbol “” will flash when a pulse is detected during the inflation.

Step 3

After the measurement, the monitor displays the systolic pressure, diastolic pressure and pulse rate.

Step 4

Press  to switch off. If left idle for 3 minutes, it will switch off automatically.

NOTE


- If you press  during measurement, the monitor will be turned off.
- If the pulse rate symbol is shown as “” instead of “”, this indicates that the monitor has detected an irregular heart beat.

► Averaging Measurement Mode

Always apply the arm cuff before turning on the monitor.


Step 1

Press and hold  for three seconds. The monitor will turn on and enter the averaging mode. Press  to select user “” (numbers from 1 to 4) whose test reading will be stored in its memory and press  to confirm. Then the cuff will begin to inflate automatically.


Step 2

After the first measurement is completed, the monitor will start counting down before the second measurement begins. The number on the right represents the remaining countdown between each measurement. The monitor will take three (3) measurements consecutively with intervals of 20 seconds.

Step 3

After taking three measurements, the results are averaged to produce the blood pressure measurement. Press  to turn off the monitor.

MONITOR MEMORY


Your monitor stores the 400 most recent blood glucose test results along with respective dates and times in the monitor memory. To recall the memory, **start with the monitor off.**

Reviewing Test Results

Step 1

Press and release .


 will appear on the display. The user symbol  (numbers from 1 to 4) appears first, press  to select which one you used to store in the memory.

Step 2

Press .

Review all test results stored in the monitor.

Step 3

Exit the memory.

Press **⏻** and the monitor will turn off.


Reviewing Day Average Results

Step 1

Press and release **M**.


When **M** appears on the display, keep pressing **M** for 3 seconds until the flashing “^{DAY}
AVERAGE” appears. Release **M** and then your 7-day average result will appear on the display.

Step 2

Press .

Review your 14-, 21-, 28-, 60- and 90- day average in AC, PC and QC mode. For example:

Step 3


Exit the memory.

Press  and the monitor will turn off.


NOTE

- Any time you wish to exit the memory, press  or leave it without any action for 3 minutes. The monitor will turn off automatically.
- If using the monitor for the first time, “—” displays when you recall the test results or review the average result. This indicates that there is no test result in the memory.

Downloading Results onto a Computer

There are two subtypes of **FORA D30** system available. Each one of them offers a different way to transmit allowing users to choose their own desired ways to manage their test results. These feature functions are listed as below:

Model	Data Transmission Function
FORA D30e 	Test results can be viewed on the PC through the cable connection.
FORA D30f 	Test results can be viewed on the PC through a wireless connection by using Bluetooth technology.

Results in the memory can be transmitted to a personal computer. The Health Care Software System and an interface cable are required before installation. To learn more about the Health Care Software System or to obtain an interface cable separately, please contact your place of purchase.

Data transmission via cable (For FORA D30e only)

Step 1


Install Software.

Download the Health Care System Software and instruction manual provided on FORA's website: <http://www.fora.com/usa> to your computer. Follow the instruction manual to install the software on your computer.

Step 2

Connect the interface cable to a serial port on PC.

With the monitor turned off, connect the interface cable to the data port located at the side of the monitor. "PCL" will appear on the display.

Step 3

Data transmission.

Follow the instructions provided in the software to transmit data. Results will reflect the real dates and times. Remove the cable or press  and the monitor will turn off.

NOTE

The monitor is unable to perform a blood glucose or blood pressure test while connecting to the PC or the server.

Data transmission via Bluetooth (For FORA D30f only)

Step 1

Install Software.

Download the Health Care System Software and instruction manual provided on FORA's website: <http://www.fora.com/usa> to your computer. Follow the instruction manual to install the software on your computer.

Step 2

Connect the monitor to PC.

With the monitor turned off, press  to initiate the data transmission of the monitor. "PCL" will appear on the display and the Bluetooth indicator on the right-hand corner will blink, indicating the monitor is ready for transmission.

Step 3

Data transmission.

Follow the instructions provided in the software to transmit data (results and the corresponding dates and times). The blue indicator will start blinking, indicating the monitor is transmitting your results onto the PC. Press  and the monitor will turn off.

NOTE

Make sure your PC has Bluetooth before transmitting the data and the monitor is within the receiving range.

MAINTENANCE

Battery

Your monitor comes with four 1.5V AA size alkaline batteries.


► Low Battery Signal

The monitor will display the two messages below to alert you when the monitor power is getting low.

Step 1

The “” symbol appears along with display messages:


The monitor is functional and the result remains accurate, but it is time to change the batteries.

Step 2

The “” symbol appears with **E-b** and **low**:

The power is not enough to do a test. You must change the batteries immediately.

► Replacing the Battery

To replace the batteries, make sure that the monitor is turned off.

Step 1

Press the edge of the battery cover and lift it up to remove.

Step 2

Remove the old batteries and replace with four 1.5V AA size alkaline batteries.

Step 3

Close the battery cover. If the batteries are inserted correctly, you will hear a “beep” afterwards.

NOTE


- Replacing the batteries does not affect the test results stored in memory.
- As with all small batteries, these batteries should be kept away from small children. If swallowed, promptly seek medical assistance.
- Batteries may leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period (i.e. 3 months or more).
- Properly dispose of the batteries according to your local environmental regulations.

Using AC Adapter

► Connect AC adapter to the monitor


Step 1

Connect AC adapter plug to AC adapter jack of the monitor.

Step 2


Plug AC adapter power plug into an electrical outlet. Press ON/OFF to start the measurement.

► Remove AC adapter from the monitor


Step 1

When the monitor is off, remove the AC adapter power plug from the electrical outlet.

Step 2

Disconnect AC adapter plug from AC adapter jack of the monitor.

Caring for Your Monitor

To avoid the meter and test strips attracting dirt, dust or other contaminants, please wash and dry your hands thoroughly before use.

► **Cleaning**

- To clean the monitor exterior, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft dry cloth. Do **NOT** flush with water.
- Do **NOT** use organic solvents to clean the monitor.
- Do **NOT** wash the pressure cuff.
- Do **NOT** iron the pressure cuff.

► **Monitor Storage**

- Storage condition: -4 °F to 140 °F (-20 °C to 60 °C), below 95% relative humidity.
- Always store or transport the monitor in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

Caring for Your Test Strips

- Storage condition: 39.2°F to 104°F (4°C to 40°C), below 85% relative humidity. Do not freeze.
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool and dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Handle the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the date you first opened the vial on the vial label
- Discard the remaining test strips after 3 months.
- Do not use the test strips beyond their expiration date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, immediately see a doctor for help.

For further information, please refer to the test strip package insert.



Important Control Solution Information

- Use only **FORA** control solutions with your monitor.
- Do not use the control solution beyond the expiration date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test should be done at room temperature (68°F to 77°F / 20°C to 25°C). Make sure your control solution, monitor, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 36°F and 86°F (2°C and 30°C). Do NOT freeze.




SYSTEM TROUBLESHOOTING



If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service. Do not attempt to repair by yourself and never try to disassemble the monitor under any circumstances.

Result Readings

MESSAGE	WHAT IT MEANS
	Appears when your result is below the measurement limit, which is less than 20 mg/dL (1.1 mmol/L).
	Appears when your result is equal to or higher than 240 mg/dL (13.3 mmol/L). This indicates the possibility of ketone accumulation for type 1 diabetes. Please seek medical assistance immediately.
	Appears when your result is higher than the limit of measurement, which is higher than 600 mg/dL (33.3 mmol/L).

Error Messages

MESSAGE	CAUSE	WHAT TO DO
	Inflation or pressure error.	Please contact local customer service for help.
	Deflation error.	
	Blood Pressure measurement error.	Refit cuff tightly and correctly. Relax and repeat measurement. If error still exists, please contact the local agent for help.
	Battery is too low.	Replace the batteries or input DC adaptor.
 	Problem with the monitor.	Review the instructions and repeat the test. If the monitor still does not work, please contact the local customer service for help.

MESSAGE	CAUSE	WHAT TO DO
	Strip has been used.	Repeat the test with a new strip.
 	Environmental temperature is outside the system's operational range.	System operational range is 50°F to 104°F (10°C to 40°C). Repeat the test after the monitor and test strip have reached the above temperature.
	User remove strip while measuring blood glucose.	Repeat the test with a new strip.

Troubleshooting

► Blood Glucose Measurement

1. If the monitor does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective monitor or test strips.	Please contact customer service.

2. If the test does not start after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatically shutoff (2 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing “  ” appears on the display.
Defective monitor.	Please contact customer service.

3. If the control solution testing result is out of range:

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiration date of the control solution.
Control solution that is too warm or too cold.	Control solution, monitor, and test strips should be at room temperature (68° F to 77° F / 20° C to 25° C) before testing.
Defective test strip.	Repeat the test with a new test strip.
Monitor malfunction.	Please contact customer service.

► Blood Pressure Measurement

1. If there is no display after pressing  :

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Batteries incorrectly installed or absent.	Check that the batteries are correctly installed.

2. If the heart rate is higher or lower than user's average:

POSSIBLE CAUSE	WHAT TO DO
Movement during measurement.	Repeat measurement.
Measurement taken just after exercise.	Rest at least 30 minutes before repeating measurement.

3. If the result is higher or lower than user's average measurement:

POSSIBLE CAUSE	WHAT TO DO
May not be in correct position while measuring.	Adjust to the correct position to measure.
Blood pressure naturally varies from time to time.	Keep in mind for next measurement.

4. If the cuff inflates again while measuring:

POSSIBLE CAUSE	WHAT TO DO
Cuff is not fastened.	Fasten the cuff again.
If the user's blood pressure is higher than the pressure the device has inflated to, the device will automatically increase its pressure and start to inflate again. Stay relaxed and wait for the measurement.	

DETAILED INFORMATION

Reference Values

Blood Glucose

Blood glucose monitoring plays an important role in diabetes control. A long-term study showed that maintaining **blood glucose levels close to normal** can reduce the risk of diabetes complications by up to 60%*¹. The results provided by this system can help you and your healthcare professional monitor and adjust your treatment plan to gain better control of your diabetes.

Time of day	Normal plasma glucose range for people without diabetes (mg/dL)
Fasting and before meal	Less than 100 mg/dL (5.6 mmol/L)
2 hours after meals	Less than 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association (2008). Clinical Practice Recommendations. Diabetes Care, 31 (Supplement 1): S1-108.

Please consult your doctor to determine a target range that works best for you.

References:

*1: American Diabetes Association position statement on the Diabetes Control and Complications Trial (1993).

Blood Pressure

Clinical studies show that adult diabetes is often accompanied by elevated blood pressure. People with diabetes can reduce their heart risk by managing their blood pressure along with diabetes treatment^{*2}.

Knowing your routine blood pressure trend tells whether your body is in good condition or not. Human blood pressure naturally increases after reaching middle age. This symptom is a result of continuous ageing of the blood vessels. Further causes include obesity, lack of exercise and cholesterol (LDL) adhering to the blood vessels. Rising blood pressure accelerates hardening of the arteries, and the body becomes more susceptible to apoplexy and coronary infarction. The recommended blood pressure range is as below:

Classification	Systolic Pressure (mmHg)	Diastolic Pressure (mmHg)
Hypotension ^{*3}	Less than 90	Less than 60
Normal ^{*4}	Less than 120	Less than 80
Pre-hypertension ^{*4}	120–139	80–89
Stage 1 Hypertension ^{*4}	140–159	90–99
Stage 2 Hypertension ^{*4}	160 or more	100 or more

References:

^{*2}: National Heart, Lung, and Blood Institute, Diseases and Conditions.

^{*3}: The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure. NIH Publication. 2003. No. 03-5233.

^{*4}: American Diabetes Association: The Diabetes-Heart Disease Link Surveying Attitudes, Knowledge and Risk (2002)

Comparing Monitor and Laboratory Results

The monitor provides you with whole blood equivalent results. The result you obtain from your monitor may differ somewhat from your laboratory result due to normal variation. Monitor results can be affected by factors and conditions that do not affect laboratory results in the same way. To make an accurate comparison between monitor and laboratory results, follow the guidelines below.

Before going to the lab:

- Perform a control solution test to make sure that the monitor is working properly.
- Fast for at least eight hours before doing comparison tests, if possible.
- Take your monitor with you to the lab.

While staying at the lab:

Make sure that the samples for both tests are taken and tested within 15 minutes of each other.

- Wash your hands before obtaining a blood sample.
- Never use your monitor with blood that has been collected in a gray-top test tube.
- Use fresh capillary blood only.

You may still have a variation from the result because blood glucose levels can change significantly over short periods of time, especially if you have recently eaten , exercised, taken medication, or experienced stress^{*5}. In addition, if you have eaten recently, the blood glucose level from a finger prick can be up to 70 mg/dL (3.9 mmol/L) higher than blood drawn from a vein (venous sample) used for a lab test^{*6}. Therefore, it is best to fast for eight hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of body fluid (dehydration) may also cause a monitor result to be different from a laboratory result.

References:

*5: Surwit, R.S., and Feinglos, M.N.: Diabetes Forecast (1988), April, 49-51.

*6: Sacks, D.B.: "Carbohydrates. " Burtis, C.A., and Ashwood, E.R.(ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.

SPECIFICATIONS

SYSTEM PERFORMANCE	
Power Source	Four 1.5V AA alkaline batteries
Size of Monitor w/o Cuff	152 (L) x 99 mm (W) x 55 mm (H), 280g without batteries.
Cuff Size	M (medium): 24-35 cm (9.4-13.8 inches) with air tube 80 cm
Memory	400 memory records.
Power Saving	Automatic power off if system idle for 3 minutes (normal mode).
System Operating Condition	50 °F to 104 °F (10 °C to 40 °C), below 85% RH
Monitor Storage Condition	-4 °F to 140 °F (-20 °C to 60 °C), below 95% RH
Power Supply Input	DC +6V / 1A (max) via Power Plug

BLOOD GLUCOSE MEASUREMENT PERFORMANCE	
Measurement Unit	mg/dL or mmol/L
Linear Range	20 to 600 mg/dL (1.1 to 33.3mmol/L)
Precision	±5 % (CV)
Accuracy	±15 mg/dL (0.83 mmol/L) when glucose < 75mg/dL
	±20% when glucose ≥ 75mg/dL
Ketone Warning	glucose value is over 13.3mg/dL

BLOOD PRESSURE MEASUREMENT PERFORMANCE	
Measurement Unit	mmHg or KPa
Systolic Measurement Range	50 mmHg -250 mmHg
Diastolic Measurement Range	30 mmHg -180 mmHg
Pulse Rate Measurement Range	40 -199 beats / minute
Maximum Inflation Pressure	300 mmHg
Accuracy of Pressure	± 3 mmHg or $\pm 2\%$ of reading
Accuracy of Pulse Rate	$\pm 4\%$ of reading
Measurement unit	Either mmHg or KPa

This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1, IEC/EN 60601-1-2, IEC/EN 61010-1, IEC/EN 61010-2-101, IEC/EN 61326-2-6.

Reference to Standards:

- EN 1060-1 /-3, NIBP-requirements
- IEC60601-1 General requirement for safety
- IEC60601-1-2 Requirements for EMC
- EN1060-4, NIBP clinical investigation
- ANSI/AAMI SP10, NIBP requirements

FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT

15.21

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

15.105(b)

Federal Communications Commission (FCC) Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference and
- 2) This device must accept any interference received, including interference that may cause undesired operation of the device.

FCC RF Radiation Exposure Statement:

1. This Transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
2. This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 centimeters between the radiator and your body.

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