

microlife BP 3BU1-5

Wrist Watch Blood Pressure Monitor Instruction Manual (1-31)

Прибор для измерения

кровяного давления на запястье (32-63)

Automatyczny nadgarstkowy ciśnieniomierz krwi

Instrukcja obsługi (64-93)

Csuklóra helyezhető vérnyomásmérő készülék Használati útmutató (94-123)

Armenia

Alfa Pharmmedimpex 8-33 Komitas ave. AM-375033 Yerevan / Armenia Tel. +374 1 569 698 Fax +374 1 544 772

Belarus

Degreemed Av. Fabrichnaia 22, 2 floor BY-220033 Minsk / Belarus Tel. +375 172 06 5349 Fax +375 172 06 5349

Bulgaria

BULOPTICA 29 "Iskar" Str. BG-Sofia 1000 / Bulgaria Tel. +359 298 08 498 Fax +359 298 11 709

Croatia

MEDICUS MATEA Goleska 20 HR-10 020 Zagreb / Croatia Tel. +385 1 654 25 55 Fax +385 1 654 25 60

Estonia

Oriola AS Estonia Saku tn. 8 EE-11314 Tallinn / Estonia Tel. +372 65 15 100 Fax +372 65 15 111

Greece

Karabinis ALFA – GAUZE Ltd. 151 Lauriou str. Peania GR-19002 Athens Tel. +30 1 66 45 751 Fax +30 1 66 46 000

Georgia

Babutsa Ltd. 13 Gamsakhurdia av. GE-380015 Tbilisi / Georgia Tel. +995 (32) 379 873 Fax +995 (32) 530 532

Hungary

MEDTRON Ltd. Budapest 62 Pf. 920 HU-1386 Budapest / Hungary Tel. +361 412 8004 Fax +361 412 8004

Latvia

Oriola-Riga SIA Dzelzavas cela 120M LV-1082 Riga / Latvia Tel. +371 780 2461 Fax +371 780 2460

Lithuania

UAB "ORIOLA VILNIUS" Laisves pr. 75 LT-2022 Vilnius / Lithuania Tel. +370 2 688 401 Fax +370 2 688 400

Poland

ECO CLIMA Mc Gregor UI. 11-go Listopada 97 PL-41800 Zabrze / Poland Tel. +48 32 274 14 26 Fax +48 32 274 14 25

Romania

Algomed Str. V.Poloni nr. 9 Sector 5 RO-Bucuresti / Romania Tel. +40 1 411 89 08

Slovakia

Interpharm Slovakia a.s. Podzahradna ul. C. 47 SI-821 06 Bratislava / Slovakia Tel. +421 245 24 4962 Fax +421 245 25 7961

Slovenia

Simps'S d.o.o. IOC Trzin Motnica 3 SI-1236 Trzin / Slovenia Tel. + 386 1 5621 338 Fax + 386 1 5621 339

Ukraine

Microlife Ukraine ul. Serova-Naberesznaya 5a UA-49000 Dniepropetrovsk / Ukraine Tel. +38 (0562) 342 206 Fax +38 (056) 778 57 63 UA-03187 Kiev / Ukraine Tel. +38 (044) 568 20 92 Fax +38 (044) 267 73 84

Wrist Watch Blood Pressure Monitor

Instruction Manual

Table of contents

1. Introduction

- 1.1. Features of the BP 3BU1-5
- 1.2. Important information about self-measurement

Important information on the subject of blood-pressure 2. and its measurement

- 2.1. How does high/low blood-pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done, if regular high/low values are obtained?
- 2.4 MAM technology (Microlife Average Mode)

The various components of the blood-pressure monitor 3.

Putting the blood-pressure monitor into operation

- 4.1. Inserting the batteries
- 4.2. Reading the set date
- 4.3. Setting the time and date and user selection

Downloaded from w 5. Manual Incompaning out a measurement

6. Printer functions

- 6.1. How to print the last stored measurement (single data print)
- 6.2. How to print all measurements stored in memory (data print out with chart)
- 6.3. Stop printing
- 7. Error messages/malfunctions
- 8. Care and maintenance, recalibration
- 9. Guarantee
- 10. Reference to standards
- 11. Technical specifications
- 12. www.microlife.com

1. Introduction

1.1. Features of the BP 3BU1-5

The blood-pressure monitor BP 3BU1-5 (with integrated time/date display) is a fully automatic, digital blood-pressure measuring device for use on the wrist, which enables very fast and reliable measurement of the systolic and diastolic blood-pressure as well as the pulse frequency by way of the oscillometric method of measuring. The device offers a very high and clinical tested measurement accuracy and has been designed to provide a maximum of user-friendliness. Before using, please read through this instruction manual carefully and then keep it in a safe place. For further questions on the subject of blood-pressure and its measurement, please contact your doctor.

Attention!

1.2. Important information about self-measurement

- Do not forget: self-measurement means control, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.
- The pulse display is not suitable for checking the frequency of heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after con-

2. Important information on the subject of blood-pressure and its measurement

2.1. How does high/low blood-pressure arise?

The level of blood-pressure is determined in a part of the brain, the so-called circulatory centre, and adapted to the respective situation by way of feedback via the nervous system.

To adjust the blood-pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. The latter is effected by way of fine muscles in the blood-vessel walls.

The level of arterial blood-pressure changes periodically during the heart activity: During the «blood ejection» (Systole) the value is maximal (systolic blood-pressure value), at the end of the heart's «rest period» (Diastole) minimal (diastolic blood-pressure value).

The blood-pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2. Which values are normal?

Blood pressure is too high if at rest, the diastolic pressure is above 90mmHg and/or the systolic blood-pressure is over 160mmHg. In this case, please consult your doctor immediately. Long-term values at this level endanger your health due to the associated advancing damage to the blood vessels in your body.

Even with normal blood-pressure values, a regular self-check with your blood-pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately. If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. Never use the results of your measurements to alter independently the drug doses prescribed by your doctor.

Table for classifying blood-pressure values (units mmHG) according to World Health Organization:

Range	Systolic	Diastolic	Measures
	Blood-pressure	Blood-pressure	
Hypotension	lower than 100	lower than 60	Check with
			your doctor
Normal range	between 100	between 60	Self-check
	and 140	and 90	
Mild	between 140	between 90	Consult your
hypertension	and 160	and 100	doctor
Moderately	between 160	between 100	Consult your
serious	and 180	and 110	doctor
hypertension			
ı www.Maruls S.com manuals sea	higher than 180	higher than 110	Consult your

Downloaded from

- hypertension». Please consult your doctor if you suspect that this might be the case.
- Correctly measured diastolic blood-pressure values above 120mmHg require immediate medical treatment.

2.3. What can be done, if regular increased/low values are obtained?

- a) Please consult your doctor.
- b) Increased blood-pressure values (various forms of hypertension) are associated long- and medium term with considerable risks to health. This concerns the arterial blood vessels of your body, which are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can be the result. Furthermore, with long-term continuously increased blood-pressure values, the heart will become structurally damaged.
- c) There are many different causes of the appearance of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.

A) Eating habits

- Strive for a normal weight corresponding to your age. Reduce overweight!
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

B) Previous illnesses

- Follow consistently any medical instructions for treating previous illness such as:
- Diabetes (Diabetes mellitus)
- Fat metabolism disorder
- Gout

C) Habits

- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (Coffee)

D) Physical constitution:

- After a preliminary medical examination, do regular sport.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please

IB BP 3BU1-5 VarB 0702 axd 13.08.2002 8:41 Uhr

- An advanced measurement accuracy is achieved by the automatic analysis of three successive measurements.
- The new system provides reliable values for the doctor and can be used as the basis for reliable diagnostics and medication therapy for high blood pressure.

A) Why MAM?

- Human blood pressure is not stable
- Scattering of devices

B) Key advantages

Reduction of:

- Device scattering
- Insufficient rest prior to measurement
- Movement artefacts
- Cuff positioning influences


C) Medical benefits

- Improved accuracy
- Reliable patient self-measurement data for the doctor
- Safe hypertension diagnostic
- Reliable therapy control

D) Measurement sequence

3. The various components of the blood-pressure monitor

The illustration shows the blood-pressure monitor BP 3BU1-5, consisting of:



Wrist cuff Type WC1 13.5–19.5 cm, for wrist circumference of 13.5 to 19.5 cm

4. Putting the blood-pressure monitor into operation

4.1. Inserting the batteries

After you have unpacked your device, first insert the batteries. The battery compartment is located on the left lateral side of the device (see illustration).

- b) Insert the batteries (2 x size AAA 1.5V), thereby observing the indicated polarity.
- c) If the battery warning appears in the display, the batteries are empty and must be replaced by new ones.

Attention!


- After the battery warning appears, the device is blocked until the batteries have been replaced.
- Please use «AAA» Long-Life or Alkaline 1.5V Batteries. The use of 1.2V Accumulators is not recommended.
- If the blood-pressure monitor is left unused for long periods, please remove the batteries from the device.

Functional check: Press and hold the 0/I button to test all the display elements. When functioning correctly all segments appear.

4.2. Reading the set date


Please press the «Time» key.






You must then re-enter the date and current time. For this, please proceed as follows (Example: Entering 2001-06-20 Time 09:30 o'-clock):

1) Press the TIME button for at least 3 seconds. The display now indicates the set year, during which the last two characters blink.



2) The correct year can be entered by pressing the MEMORY button.(Example: 1 x press)



3) Press the TIME button again. The display now switches to the current date, during which the Downloaded from www.ManufilestnCharacters(month)





4) The corresponding month can now be entered by pressing the MEMORY button.

Example: 6 x presses)





5) Press the TIME button again. The last two characters (day) are now blinking






can now be entered by pressing the MEMORY button.


(Example: 15 x presses)






8) The corresponding hour can now be entered by pressing the MEMORY button.


(Example: 9 x presses)



 Press the TIME button again. The last two characters (Minutes) now blink.



10) The exact time can now be entered by pressing the MEMORY button. (Example: 30 x presses)




User selection:


This advanced blood pressure monitor allows you to track blood pressure readings for 2 individuals independently.

- a) Before measurement make sure you set the unit for the intended user. The unit can track results for 2 individuals. (User 1, User 2)
- b) The unit is set to User 1. Click the TIME button to change to User 2.
- c) click the MEMORY button to change into User 1
- d) We suggest the first person to take their pressure to be User 1.









Further Information

5. Carrying out a measurement

5.1. Before the measurement:

- Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quite atmosphere for about ten minutes before the measurement.
- Measure always on the same wrist (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.

5.2. Common sources of error:

Note:

Comparable blood-pressure measurements always require the same conditions! These are normally always quiet conditions.

 All efforts by the patient to support the arm can increase the blood-pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary. peated after a 5 minute pause or after the arm has been held up in order to allow the accumulated blood to flow away (after at least 3 minutes).

5.3. Fitting the cuff

- a) Remove all eventual objects and jewellery (e.g. wristwatch) from the wrist in question. Draw the cuff over the wrist.
- b) The distance between the cuff and the hand should be approx. 10 mm.
- c) Secure the cuff with the Velcro fastener, so that it lies comfortably and not too tight, whereby no space should remain between the cuff and the wrist.
- d) Lay the arm on a table, with Downloaded from www.Manual People Imp Up Wards. Support the







5.4. Measuring procedure


After the cuff has been appropriately positioned, the measurement can begin:

a) Press the 0/I-button, the pump begins to inflate the cuff. In the display, the increasing cuff-pressure is continually displayed.



b) After reaching the inflation pressure, the pump stops and the pressure slowly falls away. The cuff-pressure (large characters) is displayed during the measurement. When the device has detected the pulse, the heart symbol in the display begins to blink and a beep tone is audible for every pulse beat.






c) When the measurement has been concluded, a long beep Downloaded from www.Manualoge Sounds. The measured



The measurement results are displayed, until you switch the device off. If no button is pressed for 5 minutes, the device switches automatically off, to save the batteries.

d) When the unit is set to **Average Mode**setting generally 3 separate measurements will take place in succession and
calculates your detected blood pressure
value. There will be 35 seconds resting
time in-between each measurement. A
count down indicates the remaining time
and a beep will sound 5 seconds before the
2nd and 3rd readings will begin.



In case that the single data of each cycle differ too much from each other, a fourth measurement is performed before the result will be displayed. In rare cases the blood pressure is such unstable that even after four measurements the data vary too much. In this case "ERR 6" is shown and no result can be given. If one of measurement causes an error message it is repeated.

5.5. Discontinuing a measurement


If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the 0/I




5.6. Memory – storage and recall of the measurements


The blood-pressure monitor automatically stores each of the last 30 measurement values. By pressing the MEMORY button, the last measurement (MR1) as well as the further last 29 measurements (MR2, MR3, ..., MR30) can be displayed one after the other.











(MR1: Values of the last measurement)

(MR2-MR30: Values of the measurement before MR1)

Further information


Measurements should not occur soon after each other, since otherwise, the results will be falsified. Wait therefore for several minutes in a relaxed position, sitting or lying, before you repeat a measurement.


In order to delete all stored readings, depress the MEMORY button for at least 7 seconds, the display will show the symbol "CL" and 3 short beep sounds will be heard to indicate deletion of stored readings.



Note:

If your last stored measurement was taken in the Average Mode, you will see the "MAM" symbol in approximately 3 seconds during the deletion procedure. (The display will show these symbols during deletion of all stored measurements if the last measurement is taken from the Average Mode)







6. Printer Funtions

This unit can be used in connection with the Microlife diagnostic printer. The unit can print out measurements stored in the memory by pressing the PRINTER button. There are 2 kinds of formats available.

- b) Connect the printer with the monitor and switch the printer on (for details please refer to the manual of the printer).
- c) Click the PRINTER button of the monitor once. The display will show "Pr1" and the printer will print the last stored measurement for User 1 or 2 based on your selection.





- 6.2. How to print all measurements stored in memory and a data chart, (data print out with chart).
- a) Select correct printer objective first before you depress the PRINTER button.



b) Connect the printer with the monitor and switch the printer on Downloaded from www.Manu(forodetails:please refer to the manual of the printer).


IB BP 3BU1-5 VarB 0702 axd 13.08.2002 8:41 Uhr

6.3. Stop printing

If you would like to stop printing, click the PRINTER button during the printing process.

NOTE:

Further details regarding the printer can be found in the printer instruction manual.



7. Error messages/malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (Example: Error No. 2).



Error No.	Possible cause(s)
ERR 1	No pulse has been detected.
ERR 2	Unnatural pressure impulses influence the measurement
	result. Reason: The arm was moved during the
	Measurement (Artefact).
ERR 3	The inflation of the cuff takes too long. The cuff is not

Downloaded from www.Manualslib.com manual Free of y Seated.

Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains empty	1.Check batteries for correct polarity
when the instrument	and if necessary insert correctly.
is switched on although	2.If the display is unusual, re-insert
the batteries are in place.	batteries or exchange them.
The device frequently fails	1.Check the positioning of the cuff.
to measure the blood	2. Measure the blood-pressure again in
pressure values, or the	peace and quiet under observance of
values measured	the details made under point 5.
are too low (too high).	
Every measurement	1.Please read the following information
produces a different	and the points listed under «Common
value although the	sources of error». Repeat the
instrument functions	measurement.
normally and the values	Please note: Blood pressure
displayed are normal	fluctuates continually so s
	uccessive measurements will
	show some variability.
Blood pressure	1. Record the daily development of the
m wmeasured differ from	values and consult your doctor.

Further Information

The level of blood-pressure is subject to fluctuations even with healthy people. Important thereby is, that comparable **measure-ments always require the same conditions (Quiet conditions)!** If, in spite of observing all these factors, the fluctuations are larger than 15mmHg, and/or you hear irregular pulse tones on several occasions, please consult your doctor.

For licensing, the device has been subjected to strict clinical tests, by which the computer programme used to measure the blood-pressure values was tested by experienced specialist doctors in Germany. The same computer programme is used in every individual device, and has thus also been clinically tested.


The manufacture of the devices takes place according to the terms of the European standard for blood-pressure measuring devices (see technical data) under the supervision of the Technical Monitoring Association Essen (RWTüV-Essen).

You must consult your specialist dealer or chemist if there are technical problems with the blood-pressure instrument. Never attempt to repair the instrument yourself!

Any unauthorised opening of the instrument invalidates all guarantee claims!

8. Care and maintenance, recalivration

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use petrol, thinners or similar solvent. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff must not be washed!





d) Do not drop the instrument or Downloaded from www.Mantingathi throughlysin any way.

Periodical recalibration

Sensitive measuring devices must from time to time be checked for accuracy. We therefore recommend a periodical inspection of the static pressure display every 2 years.

Your specialist dealer would be pleased to provide more extensive information about this.

9. Guarantee

The blood-pressure monitor BP 3BU1-5 is guaranteed for 2 years from date of purchase. This guarantee includes the instrument and the cuff. The guarantee does not apply to damage caused by improper handling, accidents, not following the operating instructions or alterations made to the instrument by third parties.

The guarantee is only valid upon presentation of the guarantee card filled out by the dealer.

Name and company address of the responsible dealer:

IB BP 3BU1-5 VarB 0702 qxd 13.08.2002 8:41 Uhr

10. Reference to standards

Device standard: Device corresponds to the requirements of

the European standard for non-invasive

blood- pressure monitor

EN1060-1 / 12:95 EN1060-3 / 09:97

DIN 58130, NIBP – clinical investigation ANSI / AAMI SP10, NIBP – requirements

Electromagn.

compatibility: Device fulfils the stipulations of the

European standard EN 60601-1-2

Clinical testing: The clinical performance test was carried

out in Germany according to the DIN 58130 / 1997 procedure N6

(sequential).

The stipulations of the EU-Guidelines 93/42/EWG for Medical Products Class IIa have been fulfilled.

11. Technical specifications

Weight: 148g (with batteries + cuff)

Size: 85 x 77 x 75mm (including cuff)

Storage temperature: -5 to +50°C

Humidity: 15 to 85% relative humidity

maximum

Operation temperature: 10 to 40°C

Display: LCD-Display

(Liquid Crystal Display)

Measuring method: oscillometric

Pressure sensor: capacitive

Measuring range:

SYS/DIA: 30 to 280 mmHg

Pulse: 40 to 200 per minute

Cuff pressure display range: 0-299 mmHg

Memory: toring the last 2 x 30

measurements automatically

Measuring resolution: 1 mmHg

Accuracy: Pressure within ± 3 mmHg

Pulse \pm 5 % of the reading

Power source: 2 dry cells (Batteries) UM-4,

size AAA 1.5V

Accessories: storage case