voice **MONOMETER**[®]

BLUTZUCKERMESSSYSTEM BLOOD GLUCOSE MONITORING SYSTEM

Zur Eigenanwendung
For self-testing

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Seite 3

Owner's Manual

Page 59

BLUTZUCKERMESSSYSTEM

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Lieber Monometer® voice Inhaber:

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Owner's Manual

Dear Monometer[®] voice System Owner:

This manual contains important information you must know about the system. Please read it thoroughly and carefully.

The greatest feature of this system is its **speaking function** which is an acoustic aid for users especially with visual disability.

Another unique feature is its **"no code"** function which means no need to calibrate your meter, and makes you easily monitor blood glucose at home.

**IMPORTANT SAFETY INSTRUCTIONS
READ THIS BEFORE USING**

The following basic safety precautions should always be taken.

1. Close supervision is necessary when the device is used by, on, or near children, handicapped persons or invalids.
2. Use the device only for the intended use described in this manual.
3. Do not use strips and control solutions which are not supplied by the manufacturer.
4. Do not use the device if it is not working properly, or if it has suffered any damage.
5. Before testing your blood glucose, read all instructions thoroughly and practice the test. Do all quality control checks as directed and consult with a diabetes healthcare professional.

KEEP THESE INSTRUCTIONS

TABLE OF CONTENTS

IMPORTANT SAFETY INSTRUCTIONS	61	PROBLEM-SOLVING GUIDE	102
IMPORTANT INFORMATION	64	Result Indicator	103
ABOUT ALTERNATIVE SITE TESTING (AST)	65	Error Message	104
INTRODUCTION OF THE SYSTEM	67	Problem in Operation	106
Intended Use	67	SPECIFICATIONS	108
Principle of Measurement	67	SYMBOL INFORMATION	109
Contents of the System	68	SUMMARY OF OPERATION	110
Appearance and Key Function of the Meter	69		
Meter Display Segments	70		
Appearance of the Test Strip	71		
PREPARATION BEFORE USE	72		
Battery Replacement	72		
Setting the Meter and Deleting the Memory	73		
BEFORE TESTING	78		
Checking with Monometer [®] Control Solutions	78		
Important Control Solution Information	79		
Doing a Control Solution Test	80		
TESTING YOUR BLOOD	84		
Test Procedure	85		
Expected Test Results	91		
COMPARING METER AND LABORATORY RESULTS	92		
USING THE METER MEMORY	94		
Viewing Results on the Meter	94		
Viewing Results on a Personal Computer	98		
CARING FOR YOUR METER AND TEST STRIPS	99		
Cleaning	99		
Storage	99		

● IMPORTANT INFORMATION

- ▶ Severe dehydration and excessive water loss may cause false low results. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- ▶ Test results below 60 mg/dL*¹ (3.3 mmol/L) mean low blood glucose (hypoglycemia). Test results greater than 240 mg/dL*² (13.3 mmol/L) mean high blood glucose (hyperglycemia). If you get results below 60 mg/dL or above 240 mg/dL, and do not have symptoms, first repeat the test. If you have symptoms or continue to get results that fall below 60 mg/dL or above 240 mg/dL, follow the treatment advice of your healthcare professional.
- ▶ Apply only capillary whole blood sample to test your blood glucose. Applying other substances will cause wrong results.
- ▶ If you are experiencing symptoms that are not consistent with your blood glucose test results and you have followed all instructions described in this owner's manual, call your healthcare professional.
- ▶ Inaccurate results may occur in severely hypotensive individuals or patients in shock. Inaccurate low results may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with blood glucose meters.
- ▶ Please refer to your test strip package insert for additional important information.

*1: Kahn, R., and Weir, G.: Joslin's Diabetes Mellitus, 13th ed. Philadelphia: Lea and Febiger (1994), 489.

*2: Krall, L.P., and Beaser, R.S.: Joslin Diabetes Manual. Philadelphia: Lea and Febiger (1989), 261-263.

● ABOUT ALTERNATIVE SITE TESTING (AST)

Important: There are limitations for doing AST. Please consult your healthcare professional before you do AST.

What is AST?


Alternative site testing (AST) means that people use parts of the body other than fingertips to check their blood glucose levels. This system provides you to test on the palm, the forearm, the upper arm, the calf, and the thigh with the equivalent results to fingertip testing.

What's the advantage?

Fingertips feel pain more readily because they are full of nerve endings (receptors). At other body sites, since nerve endings were not so condensed, you will not feel as much pain as at the fingertip.

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at fingertip reflects these changes faster than capillary blood at other sites. Therefore when testing blood glucose during or immediately after meal, physical exercise, or any other events, **take blood sample from your finger only.**

We strongly recommend you do AST **ONLY** in the following intervals:

- ▶ In a pre-meal or fasting state (more than 2 hours since the last meal).
- ▶ Two hours or more after taking insulin.
- ▶ Two hours or more after exercise.

Do **NOT** use AST if:

- ▶ You think your blood glucose is low.
- ▶ You are unawareness of hypoglycemia.
- ▶ Your AST results do not match the way you feel.
- ▶ You are testing for hyperglycemia.
- ▶ Your routine glucose results are often fluctuating.
- ▶ You are pregnant.

How to increase the accuracy?

Stimulating blood perfusion by rubbing the puncture site prior to blood extraction has a significant influence on the glucose value obtained. Blood from the site without rubbing exhibits a measurably different glucose concentration than blood from the finger.

Please follow suggestions below before getting a drop of blood:

- ▶ Rub the puncture site about 20 seconds before penetration.
- ▶ Use a clear cap (included in the kit) instead while setting the lancing device.

● **INTRODUCTION OF THE SYSTEM**

■ **Intended Use**

The system is intended for use outside the body (in vitro diagnostic use). It should be used only for testing glucose (sugar) and only with fresh capillary whole blood samples taken from the finger and the alternative sites including the palm, the forearm, the upper arm, the calf and the thigh. The system is intended for use in the home and in clinical settings. It should not be used for the diagnosis of diabetes or for the testing of newborns.

AST in this system can be used only during steady-state blood glucose conditions described in the section of "About AST".

■ **Principle of Measurement**


The test is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current and displays the corresponding blood glucose level. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

■ Contents of the System

The system consists of three main products: the blood glucose meter, test strips, and a control solution. These products have been designed, tested, and proven to work together as a system to produce accurate blood glucose test results. Use only the same brand-name test strips and a control solution with the blood glucose meter.

Your system includes:


- | | |
|--|-------------------------------|
| ① A meter | ② Test strips |
| ③ Sterile lancets | ④ Owner's manual |
| ⑤ Sporty carrying case | ⑥ Your first time using guide |
| ⑦ Daily log book | ⑧ Warranty card |
| ⑨ Lancing device with clear cap for Alternative Site Testing | |
| ⑩ Monometer® control solution (optional) | |

NOTE

- ▶ Please make sure that all products listed above are contained and sealed in the package before using this system. If you find any imperfection in our products, please return the whole system to the place of purchase.
- ▶ Please note that control solutions are optional. Those are not included in the standard kit. Please ask your local agent for availability.

■ Appearance and Key Function of the Meter

1. LCD DISPLAY

guides you through the test using symbols and simple messages.

2. STRIP-EJECTION BUTTON

is where the used strip will be automatically ejected after you push up the button.

3. TEST SLOT

is where you insert the test strip. The meter will turn on automatically after insertion.

4. UP BUTTON & 5. DOWN BUTTON

is used to adjust setting and the speaking volume.

6. M BUTTON

located in front of the meter with "M" on it, is used to turn on the meter and enter the memory.

7. SET BUTTON

is used to set up the meter

8. EARPHONE JACK

9. DATA PORT

located at the side, is for cable connection.

10. REPEAT BUTTON

is used to replay the previous recording which you want to hear again while operating the meter.

11. BATTERY COMPARTMENT

■ Meter Display Segments

TEST STRIP SYMBOL

Appears when the meter is turned on.

TEST RESULT AREA

Displays glucose results. A decimal point appears when the measurement unit is mmol/L.

FACE SYMBOL LOW/HIGH SYMBOL

Both are result indicators. They appear together with the test result which exceeds the reference range.

SMILE SYMBOL

Appears when the test result is within reference range.

DAY AVERAGE

Indicates that the displayed test result is an average.

TIME

DATE

BLOOD DROP SYMBOL

Flashes when it is ready to apply the sample.

VOICE SYMBOL

Indicatea speaking function is on/off.

MEASUREMENT UNIT

Appears with the test result either in mg/dL or in mmol/L .

CTL SYMBOL

Appears when doing a control test and indicates that the result will not be stored in the memory.

KETONE WARNING

Appears when the test result is equal or higher than 240mg/dL (13.3mmol/L).

MEMORY SYMBOL

Appears when you review the memory.

LOW BATTERY SYMBOL

Appears when the battery power is low.

■ Appearance of the Test Strip

This system measures the amount of sugar (glucose) in whole blood. Blood is applied to the absorbent hole of the test strip and is automatically drawn into the reaction cell where the reaction takes place.

The test strip consists of the following parts:

Contact Bars

Insert this end of the test strip into the meter. Push it in firmly until it will go no further.

Test Strip Handle


Hold this part to insert the test strip into the slot.

Confirmation Window

This is where you confirm if enough blood has been applied to the absorbent hole of the strip.

Absorbent Hole

Apply a drop of blood here, The blood will be sucked automatically.

Please see pages 84~91, "Testing Your Blood", for complete instructions.

● PREPARATION BEFORE USE

■ Battery Replacement

Your meter comes with two 1.5V AAA size alkaline batteries that are already installed.

The meter will alert you when the power is getting low by displaying two different messages:


1. When  symbol is displayed on the screen: the meter is functional and the result remains accurate, but it is time to change the batteries.

What speaks?
Battery is low, please replace.

2. The  symbol appears with the “E-b” symbol, “Error” and “low”: the batteries can not provide enough power to do a test. You must change the batteries immediately.

What speaks?
Battery is dead, please replace.


To replace the battery, make sure that the meter is turned off.

STEP1 Press the edge of the battery cover and lift it up to remove.

STEP2 Remove the old batteries and replace with two 1.5V AAA size alkaline batteries.

STEP3 Close the battery cover. If the batteries are inserted correctly, you will hear a “beep”.

PLEASE NOTE

- ▶ Replacing the batteries does not affect the test results stored in memory.
- ▶ Batteries might leak chemicals if not used for a long time. Remove the batteries if you are not going to use the device for an extended period (i.e., 3 months or more).


WARNING

As with all small batteries, the batteries should be kept away from small children who still put things in their mouths. If they are swallowed, promptly see a doctor for help.

■ Setting the Meter and Deleting the Memory


Your meter comes with the time, date, unit of measurement, unit of temperature preset, speaking function and memory deletion. If you need to set these parameters, please follow below steps.

Start with the meter off. Then keep pressing the S button. The meter is now in the setting mode. You can start setting up the meter.


STEP 1. Set the Year

Press the “Up” or “Down” button until the desired year appears. Press the S button and then the month segment flashes.



STEP 2. Set the Month

Press the “Up” or “Down” button until the desired month appears. Press the S button and then the day segment flashes.



STEP 3. Set the Day

Press the “Up” or “Down” button until the desired day appears. Press the S button and then the hour segment flashes.


STEP 4. Set the Hour

Press the “Up” or “Down” button until the desired hour appears. Press the S button and then the minute segment flashes.



STEP 5. Set the Minutes

Press the “Up” or “Down” button until the desired minute appears. Press the S button and then the current unit of measurement flashes.

STEP 6. Select mg/dL or mmol/L

Press the “Up” or “Down” button to select the unit of measurement you want to use. Press the S button and then the current unit of temperature flashes.


WARNING

The milligram per deciliter (mg/dL) is the standard unit in the United States. The millimole per liter (mmol/L) is the standard unit in Canada. Use of the wrong unit of measure may cause you to misinterpret your blood glucose level, and may lead to incorrect treatment.

STEP 7. Select °C or °F

Press the “Up” or “Down” button to select the unit of temperature you want to use. Press the S button and the meter will display “VOL”, “” and flashing number.

Temperature unit is degrees Fahrenheit.



STEP 8. Select Speaking Function


Press the “Up” and “Down” buttons until you hear the desired speaking volume. Number 0 indicates that the speaking function is turned off, where “” does not display on LCD during testing. Number 1 to 7 indicates speaking volume from low to high, where “” displays on LCD during testing. Press the S button and the meter will display “dEL” with flashing “” symbol.

Voice volume is 7.

STEP 9. Delete the Memory

When “dEL” and blinking “” symbol is displayed, if you do not want to delete memory, press the S button again to skip this step and turn off the meter. If you would like to delete **ALL** memory, press the “Up” button twice. “” and “” are displayed on the meter, which means that all data stored are deleted. The memory of this meter is empty now. The meter will automatically turn off.

Delete memory record function. Press up key to clear all records.
↓ Are you sure you want to delete ALL records? Press up key again to clear all records.
↓ All memory records have been deleted.
↓ Have a nice day.

PLEASE NOTE

- ▶ The time, date, unit of measurement and unit of temperature can **ONLY be changed** in the setting mode. Therefore, when you perform a glucose testing, those parameters are not possible to be changed.
- ▶ Your meter displays 7-, 14-, 21-, 28-, 60- and 90-day averages which you can access from the meter memory. These averages are calculated from results obtained during the 7-, 14-, 21-, 28-, 60- and 90-days preceding the current date and time settings. When the date and time are changed, the 7-, 14-, 21-, 28-, 60- and 90-day averages may change.
- ▶ While the meter is in the setting mode, if no button is pressed for 3 minutes, the meter will turn off automatically.

● BEFORE TESTING

■ Checking with Monometer® Control Solutions

Monometer® control solutions contain a known amount of glucose that reacts with test strips. By comparing your control solution test results with the expected range printed on the test strip vial label, it is able to check that the meter and the test strips are working together as a system and that you are performing the test correctly. It is very important that you do this simple check routinely to make sure you get accurate results.

How often the control solution test should be performed?

- ▶ When you use this system to test your blood for the first time, practice the procedure using control solution. When you can do three tests in a row that are within the expected range, you are ready to test your blood.
- ▶ For routinely check the meter and test strips, perform a single test for each level of control solution at least once a week.

When the control solution test should be performed?

- ▶ When you first get your Glucose Meter.
- ▶ When you begin using a new vial of test strips.
- ▶ Whenever you suspect that the meter or test strips are not working properly.
- ▶ When your blood glucose test results are not consistent with how you feel, or when you think your results are not accurate.
- ▶ When your test strips are exposed to extreme environmental conditions. (See **Storage** section of this manual).
- ▶ When you want to practice running the test.
- ▶ If you drop the meter.

■ Important Control Solution Information

- ▶ Use only Monometer® control solutions.
- ▶ Check the expiration date on the control solution vial. Do not use if expired.
- ▶ Control solution, meter, and test strips should come to room temperature (20~25°C/68~77°F) before testing.
- ▶ Shake the vial, discard the first drop of control solution, and wipe off the dispenser tip to ensure a good sample and an accurate result.
- ▶ Use only for 90 days after first opening. Record the discard date (date opened plus 90 days) on the control solution vial. Discard after 90 days.
- ▶ Store the control solution tightly closed at temperatures below 30°C (86°F). Do not refrigerate.

PLEASE NOTE


The control solution range printed on the test strip vial is for Monometer® control solution only. It is used to test meter and test strip performance. It is not recommended range for your blood glucose level.

■ Doing a Control Solution Test

STEP 1. Insert the Test Strip

Insert a test strip with contact bars end first and facing up, into the test slot. The meter turns on automatically and displays the followings in sequence:

- “CHK” and “”
- “” and ambient temperature
- “” and flashing “”


What speaks?

Thank you for using this product.

- ↓ Test Strip Inserted.
- ↓ Today's date is Sunday, January 18th, 2008. The time is 10:08AM.
- ↓ The room temperature is 30.6 degrees Celsius.
- ↓ Please apply blood onto the test strip.

STEP 2. Press the **M** button

While the “” symbol appears on the display, press the M button and then “CTL” will appear on the display. With the “CTL” sign on the display, the meter will not store your test result in memory. If you decide not to perform a control solution test, press the M button again, and the “CTL” sign will disappear.

What speaks?


You are now in control solution mode.

WARNING

- ▶ Contact bars must be inserted all the way into the meter or you may get an inaccurate test result.
- ▶ Every time you perform a control solution test, you must enter into the “CTL” test mode so that the test result will not be stored in the meter memory. Failure to do so will confuse the blood glucose test result with the control solution test result in memory.

STEP 3. Obtain Control Solution.


Shake the control solution vial well. Remove the cap. Squeeze the vial, discard the first drop, and wipe off the remains on the tip of dispenser. Squeeze the vial again to get another drop of control solution on your fingertip.

STEP 4. Apply Control Solution.



Move your finger (with control solution on it) to meet the absorbent hole of the test strip and the drop will be automatically drawn into the test strip. Remove your finger until the confirmation window is filled. The meter begins to count down.

To avoid contaminating the control solution with the content of the test strip, you have to place a drop of control solution on a clean surface or on your fingertip. Do not directly apply control solution into a strip.

STEP 5. Read and Compare the Result

After counting to 0, the test result of control solution is shown on the screen. Compare this result with the range printed on the test strip vial. It should fall within this range.


Out-of-range results

If test results fall outside the range printed on the test strip vial, check the section of "**Problem in Operation**" in troubleshooting guide and repeat the test. If you continue to get out-of-range results, it means that the system may not be working properly. Do NOT test your blood. Contact the Customer Service for help.

● TESTING YOUR BLOOD

Be sure to read this section and the test strip package insert found in the test strip box carefully before testing. Make sure you have all items needed to test:

- A. Blood Glucose Meter
- B. Test Strip
- C. Lancing Device
- D. Sterile Lancet
- E. Clear Cap (For AST use)

WARNING

To reduce the chance of infection:


- ▶ Never share a lancet or the lancing device with others.
- ▶ Always use a new, sterile lancet. Lancets are for single use only.
- ▶ Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.

■ Testing Procedure


Wash and dry your hands first before starting.

STEP 1. Set the Lancing Device


Screw off the cap of lancing device. Insert a lancet into the lancet holder and push down firmly until it is fully seated.

Twist the protective disk until it separates from the lancet.


Replace the lancing device cap. Turn the cap until it is snug but not too tight.


Twist the adjustable tip in either direction until the number lines up with the Arrow:
The adjustable tip offers 5 levels of skin penetration.

- 1-2 for soft or thin skin
- 3 for average skin
- 4-5 for thick or calloused skin


Slide the ejection/cocking control back until it clicks. If it does not click the device may have been cocked when the lancet was inserted.

► Blood from sites other than the fingertip


A clear cap, together with the kit, makes it easier to get a drop of blood for AST. When you want to obtain blood from sites other than the finger, **replace the lancing device cap with the clear cap**. Turn the clear cap until it is snug but not too tight, and then slide the ejection/cocking control back until it clicks.


The lancing device is now ready for use. Set aside for later use.

STEP 2. Insert Test Strip

Insert a test strip with contact bars end first and facing up into the test slot. The meter turns on automatically and displays the followings in sequence:

- “CHK” and “”
- “” ambient temperature
- “”, and flashing “”

Thank you for using this product.

↓ Test Strip Inserted.

↓ Today's date is Sunday, January 18th, 2008. The time is 10:08AM.

↓ The room temperature is 30.6 degrees Celsius.


↓ Please apply blood onto the test strip.

STEP 3. Get a Drop of Blood

Select the puncture site either in finger or in other parts (AST). Clean the puncture site with 70% alcohol cotton and **let it air-dry**.

► Fingertip

Hold the lancing device firmly against the side of your finger. Press the release button. You will hear a click, indicating that the puncture is complete.



► Sites other than fingertip

Please refer to the section of “About AST” for available punctured sites.

After penetration, discard the first drop of blood with a clean tissue paper or cotton. Then gently squeeze the punctured area to obtain blood. But be attention **NOT to smear the blood sample**.

The volume of blood sample must be at least 0.7 microliter (• actual size).



PLEASE NOTE

- Choose a different spot each time you test. Repeated punctures in the same spot may cause soreness and calluses.
- Before you decide to do AST, please consult your health professional firstly.
- Since the first drop of blood usually contains tissue fluid and serum, which may affect the test result, it is recommended to be discarded.


STEP 4. Apply Blood into the Test Strip

When “” is flashing on the screen, apply your blood to the absorbent hole of the test strip until the confirmation window is fully covered with blood. The meter then begins to count down automatically.

STEP 5. Obtain an Accurate Result in 7 Seconds

The result of your blood glucose test is shown after the meter counts to 0. This reading is automatically stored in the meter.

PLEASE NOTE

- Do not push your finger (with blood on it) against the test strip or try to apply a smeared sample on the test strip.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to restart the test procedure.
- The blood should be completely filled the confirmation window before the meter begins to count down. If you find that the confirmation window is not filled with blood when the meter is counting, **NEVER** try to add more blood to the test strip. **Discard the test strip and retest with a new one.**
- If you have trouble filling the test strip, please contact the Customer Service for help.

STEP 6. Eject the Used Test Strip

After finishing the measurement, you can either take out the used strip by simply pushing up the Strip-Ejection button or remove the test strip directly by your hand. "OFF" is shown after the used strip is ejected and the meter will shut off automatically afterward.


PLEASE NOTE

Please make sure the ejection button is pointed into the disposal container.

STEP 7. Remove the Lancet

Always use caution when removing the lancet.

Take the lancet out carefully. Place the disk on a hard surface and push the exposed tip into the protective disk.

WARNING

The used lancet and the used test strip may be potentially biohazard. Please consult your healthcare professionals for how to properly dispose of the used lancets and the used test strip.

Expected Test Results

Blood glucose monitoring plays an important role in diabetes control. A long-term study showed that **keeping blood glucose levels close to normal** can reduce the risk of diabetes complications by up to 60%*³. The results you get with the Monometer[®] voice system can help you and your healthcare professional monitor and adjust your treatment plan to gain better control of your diabetes.

Time of day	Plasma glucose range for people without diabetes	Your target range
Fasting and before meal ⁴	70~110 mg/dL (3.9 ~6.1 mmol/L)	
2 hours after meals ⁵	Less than 140 mg/dL (7.8 mmol/L)	

*3: American Diabetes Association position statement on the Diabetes Control and Complications Trial (1993).

*4: Sacks, DB in "Carbohydrates", Burtis, CA, Ashwood, ER (ed), Tietz Textbook of Clinical Chemistry, Philadelphia, WB Saunders Company, 1999.

*5: ADA Clinical Practice Recommendations 2003.

Please work with your doctor to determine a target range that works best for you.

● COMPARING METER AND LABORATORY RESULTS

Test results from the meter and laboratory are both expressed in plasma-equivalent units. However, the result you obtain from your meter may differ somewhat from your laboratory result due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way. (See test strip package insert for typical accuracy and precision data, and for important information on Limitations.) To make an accurate comparison between meter and laboratory results, follow the guide-lines below.

Before you go to the lab:

- ▶ Perform a control solution test to make sure that the meter is working properly.
- ▶ It is best to fast for at least eight hours before doing comparison tests.
- ▶ Take your meter with you to the lab.

While at the lab:

Make sure that the samples for both tests (the meter test and the lab test are taken and tested within 15 minutes of each other).

- ▶ Wash your hands before obtaining a blood sample.
- ▶ Never use your meter with blood that has been collected in a gray-top test tube.
- ▶ Use fresh capillary blood only.

You may still have a variation from the result because blood glucose levels can change significantly over short periods, especially if you have recently eaten, exercised, taken medication, or experienced stress⁶. In addition, if you have eaten recently, the blood glucose level from a finger stick can be up to 70 mg/dL (3.9 mmol/L) higher than blood drawn from a vein (venous sample) used for a lab test⁷. Therefore, it is best to fast for eight hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of body fluid (severe dehydration) may also cause a meter result to be different from a laboratory result.

References

*6: Surwit, R.S., and Feinglos, M.N.: Diabetes Forecast (1988), April, 49-51.

*7: Sacks, D.B.: "Carbohydrates." Burtis, C.A., and Ashwood, E.R. (ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.


● USING THE METER MEMORY

■ Viewing Results on the Meter

Your meter stores the latest 450 results of your blood glucose with date and time in its memory. It also automatically calculates the averages of your blood glucose in the intervals of 7, 14, 21, 28, 60 and 90 days. You can review the results easily by the following steps.

1. Recall the Stored Test Results


STEP 1. When the meter is off, press and release the M button. The screen shows “M”. Press the M button again, “01” appears first and then the latest glucose result along with date and time will be shown on the screen.

What speaks?

Thank you for using this product.
↓ Today's date is Sunday, January 18th, 2008. The time is 10:08 PM.
↓ Press main key for memory recall mode.
↓ Memory Record for January 10th, 2008. 10:08 PM. Your blood glucose was 112 milligrams per deciliter.

STEP 2. Press the “Up” and “Down” button to recall the test results stored in the meter.

What speaks?

Memory Record for August 23th, 2007. 8:36AM. Your blood glucose was 80 milligrams per deciliter.

STEP 3. Exit the memory mode, keep pressing the M button and the meter will be turned off.

What speaks?


Have a nice day.

2. Read the Average of Blood Glucose Results:

STEP 1. When the meter is off, press and release the M button.

The screen shows "M". Keep pressing the M button for 2-3 seconds, until blinking "DAY AVG" appears. Release the M button and then the 7-day average result will appear on the screen.

Thank you for using this product.

↓ Today's date is Sunday, January 18st, 2008. The time is 10:08 PM.

↓ Press main key for memory recall mode.


↓ The 7 day average is 82 milligrams per deciliter for 66 records.

STEP 2. Press the "Up" and "Down" buttons to review the 14-, 21-, 28-, 60- and 90- day average. Like the 7-day average, the 14- day average and the times of performing test in the past 14 days will be shown on the screen.

The 14 day average is 83 milligrams per deciliter for 87 records.

STEP 3. Exit the memory mode, keep pressing the M button and the meter will be turned off.


PLEASE NOTE

- ▶ The control solution results are NOT stored in the memory (please also go to page 24 WARNING for information). The list of past results and the result average are for blood glucose results only.
- ▶ When using the meter for the first time, "----" is displayed when you recall the test results or review the average result. It means that there is no test result in memory.

No memory record available.

- ▶ Any time you want to exit the memory mode, keep pressing the M button for 5 seconds or leave it without any action for 3 minutes. The meter will turn off automatically.

■ Viewing Results on a Personal Computer


Results in memory can be transmitted to the personal computer. Health Care Software System and an Interface Cable are required before installation. The software can be downloaded from the home page of TaiDoc, Inc. The interface cable is an optional accessory. To learn more about Diabetes Management Software or to obtain an Interface Cable separately, please contact your local customer service.

STEP 1. Install Software

Install Health Care Software System on your Computer by following the instructions provided on TaiDoc's website:
www.taidoc.com

STEP 2. Connect to Personal Computer

Connect the interface cable to a serial port of your computer. With the meter turned off, connect the Interface Cable to the Data Port of the meter. "PC" will appear on the display, indicating that the meter is ready to transmit data.

Step 3 Transmit Data

Follow the instructions provided in the software to transmit data. Results transmitted will include date and time. Remove the cable and the meter will automatically turn off.

PLEASE NOTE

While the meter is connected to the PC, it is unable to perform a blood glucose test.

● CARING FOR YOUR METER AND TEST STRIPS


To avoid the meter and test strips getting dirt, dust or other contaminants, please wash and dry your hands thoroughly before use.

■ Cleaning


1. To clean the meter exterior, wipe with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft and dry cloth. Do not flush with water.
2. Do not use organic solvents to clean the meter.

■ Storage

1. Meter Storage

- ▶ Storage condition: $-20^{\circ}\text{C}\sim 60^{\circ}\text{C}$ ($-4^{\circ}\text{F}\sim 140^{\circ}\text{F}$) below 95% relative humidity.
- ▶ Always store or transport the meter in its original storage case.


- ▶ Avoid dropping and strong impact.
- ▶ Avoid direct sunlight and humidity.


2. Strip Storage


- ▶ Storage condition: 4°C~40°C (39.2°F~10°F), below 85% relative humidity. Do Not Freeze.
- ▶ Store your test strips in their original vial only. Do not transfer to other container.

- ▶ Store test strip packages in a cool and dry place. Keep away from direct sunlight and heat.
- ▶ After removing a test strip from the vial, immediately replace the vial cap and close it tightly.
- ▶ Touch the test strip with clean and dry hands.


- ▶ Use each test strip immediately after removing it from the vial.
- ▶ Write the discard date (the date opened plus 90 days) on the vial label when you first open it. Discard remaining test strips 90 days after first opening date.

- ▶ Do not bend, cut, or alter a test strip in any way.
- ▶ Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

3. Control solution storage

- ▶ Storage condition: Store the control solution tightly closed at temperatures below 30°C (86°F). Do NOT refrigerate.
- ▶ Record the discard date (date opened plus 90 days) on the control solution vial. Discard after 90 days.


● PROBLEM-SOLVING GUIDES

Following is a summary of some display messages. These messages help to identify certain problems but do not appear in all cases when a problem has occurred. Improper use may cause an inaccurate result without producing an error message.




In the event of a problem, refer to the information under action.

Never try to disassemble the meter in any circumstances. If you encounter any error messages not listed below or if you have followed the actions recommended below but the problem keeps unsolved, please call the customer service for support.




■ Result Indicator

Message	What it means
	Lo appears when your result is below measurement limit, which is less than 20 mg/dL (1.1 mmol/L).
	low appears when your result is between 20 and 69 mg/dL (1.1 and 3.8 mmol/L). It indicates that the result is below reference range.
	☺ appears when your result is in the reference range from 70 to 119 mg/dL (3.9 and 6.6 mmol/L).
“Lo” or “ low ☹” symbol indicates hypoglycemia (low blood glucose.) You should contact with your healthcare professionals for further treatment immediately.	
	high appears when your result is equal or greater 120 mg/dL (6.6 mmol/L). It indicates that the result is higher than reference range.
	KETONE? and high are shown when your result is equal or higher than 240 mg/dL (13.3 mmol/L). This indicates the possibility of ketone accumulation for Type 1 diabetes. Please seek medical assistance immediately.
	HI is displayed when your result is higher than the limit of the measurement, which is higher than 600 mg/dL (33.3 mmol/L).

■ Error Message

MESSAGE	WHAT DOES METER SAY?	WHAT IT MEANS	ACTION
	Battery is dead. Please replace.	Appears when the batteries can not provide enough power for a test.	Replace the batteries immediately.
	Strip has been used. Please replace.	Appears when inserting a used test strip.	Test with a new test strip.
	The room temperature is out of range. Thus unable to measure.	Appears when ambient temperature is below system operation range.	System operation range is 10~40°C (50~104°F). Repeat the test after the meter and test strip have reached the above temperature.
		Appears when ambient temperature is above system operation range.	

■ Error Message

MESSAGE	WHAT DOES METER SAY?	WHAT IT MEANS	ACTION
	Please ensure proper operation then try again.	Remove the strip after applying blood to the absorbent hole.	Re-test with a new test strip.
	Please remove strip.	Test strip is inserted after the meter turns off.	Remove the test strip.
	Meter malfunction. Please contact local distributor.	Problem with the meter.	Review the instructions and re-test with a new test strip. If the above steps do not work, please contact the Customer Service for help.
			
			

■ Problem in Operation

1. If the meter does not display a message after inserting a test strip:

PROBABLE CAUSE	WHAT TO DO
Battery exhausted.	Replace the batteries.
Battery incorrectly installed or absent.	Check that the battery is correctly installed.
Test strip inserted upside down or incompletely.	Insert the test strip correctly with the contact bars end first and facing up.
Defective meter.	Please contact the Customer Service.

2. If the test does not start after applying the sample:

PROBABLE CAUSE	WHAT TO DO
Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatically shutoff (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when "💧" flashing appears on the display.
Defective meter.	Please contact the Customer Service.

3. If the control solution test result is out of range.

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read the instruction thoroughly and repeat the test again.
Do not shake the control solution vial very well.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiration date or the discarded date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should come to room temperature (20~25°C/68~77°F) before testing.
Test strip deterioration.	Repeat the test with a new test strip.
Meter malfunction.	Contact the Customer Service.

● SPECIFICATIONS

Model No.: Monometer[®] voice

Dimension & Weight: 95mm(L) x 55mm(W) x 18mm(H) 75 g

Power source: two 1.5V AAA alkaline batteries

Display: LCD

Memory: 450 measurement results with date and time

External output: RS232 PC interface

Auto electrode inserting detection

Auto sample loading detection

Auto reaction time count-down

Auto turn-off after 3 minutes without action

Temperature Warning

Operating condition: 10°C~40°C, below 85% R.H. (noncondensing)

Storage/Transportation condition: 4°C~40°C, below 85% R.H.

Measurement Range: 20~600mg/dL (1.1~33.3mmol/L)

Measurement Units: mg/dL or mmol/L

The device has been certified to meet the electrical and safety requirements of:

IEC 60601-1, EN 60601-1, IEC 61010-1, EN 61010-1, EN 61010-2-101, EN 60601-1-2, EN 61326

● SYMBOL INFORMATION


Symbol	Bedeutung
	In vitro diagnostic medical device
	Do not reuse
	Consult instructions for use
	Keep away from sunlight
	Keep dry
	Temperature limitation
	Use by
	Date of manufacture
	Batch code
	Manufacturer
SN	Serial number
	Caution, consult accompanying documents
	Biological risks
	Authorized representative in the European Community
	Corresponds to the current European Guidelines

● SUMMARY OF OPERATION

This summary is intended only for quick reference and can not be taken as the substitute for the owner's manual. Please read the entire manual before you begin doing test.

STEP 1. Insert test strip

The meter is automatically turned on and displayed room temperature.

What speaks?

Thank you for using this product.



↓ Test Strip Inserted.

↓ Today's date is Sunday, January 18th, 2008. The time is 10:08AM.

↓ The room temperature is 30.6 degrees Celsius.

STEP 2. Apply sample.

Touch and hold the drop of blood to the absorbent hole of test strip until confirmation window is completely filled.

What speaks?

Now testing.

STEP 3. Obtain accurate results in 7 seconds

The meter starts to count down. Never try to add any blood into the absorbent hole even if you find that the confirmation window is not fulfilled. Discard the strip and retest with a new strip.

What speaks?

You blood glucose is 92 milligram per deciliter.

STEP 4. Eject the used strip automatically.

STEP 5. Discard the used strip and the lancet according to your local regulations.