## Ssmart Talking Blood Pressure Monitor Model: BPU3210S

### User Manual

### WARNING

- This monitor is for home daily use only and is not meant to replace any medical device.
- The monitor should be used by adults whose age is at least 18 or above.
- No modification of this equipment is allowed.
- Do not modify this equipment without authorization of the manufacturer.
- If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.
- Remove the batteries whenever you are planning to store the unit for a long period of time.

## CONTENTS

Introduction .....	3	Application .....	10
Package contents .....	3	Connect with your mobile device .....	10
Accessories .....	3	Start data transmission .....	10
Product view / key features .....	4	Review the records .....	11
Setting up the monitor .....	4	Unique device ID .....	12
Assemble the unit .....	4	Troubleshooting .....	12
Installing and replacing the batteries .....	5	Safety and care instructions .....	12
Monitor settings .....	5	Safety precautions .....	12
Positioning the arm cuff .....	6	Caring for your blood pressure monitor .....	13
Preparing to take blood pressure measurement .....	6	Useful reference .....	14
Taking a measurement .....	7	About blood pressure .....	14
Power off .....	7	How the blood pressure monitor works .....	15
Understanding the reading .....	7	Operation frequency .....	15
Display .....	7	Specifications .....	15
Hypertension .....	8	Electromagnetic compatibility descriptions .....	16
Irregular heartbeat / arrhythmia .....	8	About Oregon Scientific .....	21
Movement / reading error .....	8	Authorized representative .....	21
Voice announcement .....	8	EU-Declaration of conformity .....	22
View records .....	9	Symbol description(s) .....	22
Recalling measurements stored in memory .....	9	FCC statement .....	22
Clear records .....	9		
Deleting a single record .....	9		
Deleting all records .....	10		

## INTRODUCTION

Thank you for selecting the Oregon Scientific™ Ssmart Talking Blood Pressure Monitor (BPU321OS). The BPU321OS features multi-language voice announcement and playback of blood pressure measurements for added convenience and has been designed to provide you with many years of reliable service. The monitor can help you measure and track the following metrics:

- Systolic pressure
- Diastolic pressure
- Average morning / evening blood pressure
- Heart rate

Our BPU321OS product is able to store up to 100 measurements per user (maximum 2 users). The monitor carries Bluetooth Smart technology (Bluetooth Low Energy,  Bluetooth) that is compatible with your Smart Phone and tablet PC. Now you can easily view and systematically manage your blood pressure data through our App “Ssmart Blood Pressure Monitor”, which is not only available in Apple App Store, but also in Google Play Store.

Readings taken by the blood pressure monitor are equivalent to those obtained by a trained observer using the cuff and stethoscope auscultation method. Clinical performance was successfully done against ISO 81060-2: 2009 and international protocol. The B.H.S. (British Hypertension Society) has rated this product “recommended for clinical and home use”, the highest grading available for blood pressure monitors.

**NOTE** Please keep this manual handy as you use your new product. It contains practical step-by-step instructions, as well as technical specifications and warnings you should know about.

## PACKAGING CONTENTS

This package should include the following items:

- 1 x BPU321OS Blood Pressure Monitor
- 1 x Cuff wrap
- 1 x Cuff holder
- 4 x AA (1.5V) batteries

## ACCESSORIES


This product comes with below optional accessories which can be purchased separately. For more information, please contact your local retailer.\*

- Replacement Arm cuff
- Storage pouch
- 6V adaptor

\* Features and accessories will not be available in all countries.

**NOTE** If you want to use a power adaptor for this Blood Pressure Monitor (BPU321OS), use ONLY the power adaptor purchased from Oregon Scientific that approved by the strict medical grade test. To order an adaptor, visit our website at [www.oregonscientific.com](http://www.oregonscientific.com). Do NOT attempt to use other brand in the market for this product that may cause unpredictable harm to the product or the user. Oregon Scientific is not liable for any misuse of the adapter.

## PRODUCT VIEW / KEY FEATURES


- A. : start taking blood pressure for user 1; view user 1 measurement records, view additional / older records; increase value of the setting; toggle setting options
  - B. : turn monitor on/off; enter setting mode; confirm settings; start taking measurement for guest (no record stored); cancel measurement; confirm deletion of individual or all data for user (1 or 2).
  - C. : start taking blood pressure for user 2; view user 2 measurement records, view additional / older records; decrease value of the setting; toggle setting options.
1. / : indicates which user (1 or 2) has been selected (no icon appears for guest).

2. : Pairing status (**bLE / Err / donE**); Clock; calendar; date and time stamp.
3. : Average morning / evening blood pressure
4. : memory is full; memory mode selected.
5. Systolic pressure
6. Air outlet
7. Diastolic pressure
8. : monitor battery low
9. : voice activation
10. Synchronizing status indicator light
11. Blood pressure classification bar
12. Heart rate counter
13. : irregular heart beat detected

## SETTING UP THE MONITOR

### ASSEMBLE THE UNIT


1. Insert the small end of the air tube into the air outlet on the monitor.
2. Snap cuff holder into place at the back of the monitor.
3. Place cuff into holder.

## INSTALLING AND REPLACING THE BATTERIES

This monitor can work with a 6V adaptor which has to be purchased separately.

1. Slide off the battery cover.
2. Install the batteries by matching the correct polarities. Always use the correct battery
3. Align the cover flat against the battery compartment, and then slide it upwards.

Above diagram displays:

1. Air outlet
2. Battery compartment

**NOTE** Flashing  indicates monitor batteries are low.

**NOTE** Replace batteries promptly to ensure accurate readings. Replace all the batteries at the same time; it is dangerous to mix old and new batteries.

**NOTE** Contact your local waste disposal authority for instructions on how to dispose of used batteries. Used batteries can be harmful to the environment, and should not be thrown out with household trash.

## MONITOR SETTINGS

**IMPORTANT** As blood pressure measurements are recorded in the memory with date and time stamp for easy analysis, set the monitor prior to first use to ensure measurements are correctly marked.

To set the monitor:

- When monitor is off, press and hold  to enter setting mode.

<b>Settings</b>	Setting options (Press  or  to adjust the value, and press  to confirm)
<b>Date</b>	Year, month, date
<b>Time</b>	Hour format, hour, minute
<b>Audio function</b>	On / off Volume control (HI / Lo) when function enabled Language: (En) English (Fr) French (CH) Putonghua (CH Hk) Cantonese (SP) Spanish

## POSITIONING THE ARM CUFF


The cuff circumference is about 22-32cm. If your arm's circumference is bigger than this, the cuff cannot be used tightly and cause inaccurate reading.

Ensure the arm cuff is correctly placed to obtain an accurate reading.


**NOTE** If your physician has diagnosed you with poor circulation in your left arm, use your right arm.

1. Choose the position from which you wish to measure, sitting or lying down.
2. Roll or push up your sleeve to expose the skin.
3. Put arm through the cuff loop.
4. Place cuff about an inch (2-3 cm) above your elbow.

Use cushions or other supports if necessary to adjust.

The air tube should run down, between the inside of your forearm and the side of your chest, in line with your middle finger.

5. Fasten the arm cuff so it fits snugly. If the cuff is too loose, the measurement will not be accurate.

## PREPARING TO TAKE BLOOD PRESSURE MEASUREMENT

You can choose to take blood pressure while sitting or lying down. Below are some helpful tips for taking a measurement.

- Be sure to set the clock before taking your first measurement, or whenever you replace the batteries, so that the date and time are stored in the memory with your history. For instructions, refer to **Monitor Settings** section.
- It is important to relax when taking your blood pressure. Try to take a 15-minute rest before you begin.
- Do not lean backward or bend your upper arm inward while taking a measurement.
- Do not sit cross-legged and avoid talking or moving your fingers and hand while taking a measurement. Rapid movements or other activities may alter your reading.
- Wait at least an hour before taking your blood pressure if you have just eaten a large meal.
- Do not smoke or drink alcohol before taking your blood pressure.
- Do not measure your blood pressure if you are under stress.
- Remove all accessories from your arm.

- Wait at least 3 minutes between measurements. This allows your blood circulation to recover.

For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.

## TAKING A MEASUREMENT

If you want the monitor to synchronize with the application of your mobile device. Please do preparation before starting a measurement.

1. Turn on the application and ensure the Bluetooth connection on your smart device is good.
2. Relax your arm.
3. Press  to turn on the unit.
4. Select user, then a voice prompts you to sit in the correct position and to remain still :

Button	User	Memory
	Guest	No records stored
	1	Record stored under User 1
	2	Record stored under User 2

5. Arm cuff will begin inflating automatically. Digits will appear on the screen.
6. Once the pressure reaches a certain level, it will slowly deflate until the results are displayed on the LCD.
7. The monitor automatically attempts to send the record wirelessly to your mobile device. If data transmission successes, **bLE** will display on the top of the monitor screen.

8. To conserve battery life, monitor will turn off 1 minute after taking a measurement. Alternatively, press to manually turn off monitor.

**NOTE** Please sanitize the monitor and main accessories for avoid cross infection.

## POWER OFF


To stop the measurement process at any time, press .

## UNDERSTANDING THE READING

For reference material, please refer to **Useful Reference** section which can be found at the end of this manual.

## DISPLAY

Example:

1. Date / Time stamp (automatically toggles)
2. User selected
3. Systolic pressure

4. Diastolic pressure
5. Heart rate
6. Blood pressure classification segments
7. Hypertension indicator bar

## HYPERTENSION

The monitor will visually display by segments under which WHO classification the measurement falls under.

Bar	Classification	Reading
█	Grade 3 Hypertension	High
	Grade 2 Hypertension	
█	Grade 1 Hypertension	Slightly high
█	High Normal	Normal
	Normal	
	Optimal	

## IRREGULAR HEARTBEAT / ARRHYTHMIA

Apart from blood pressure and heart rate, heartbeat is also an indicator of your general health condition.

With a normal heartbeat, the beats are steady at regular intervals.

An irregular heartbeat is one where there is a variance and the intervals vary.

**NOTE** If an irregular heart beat has been detected during the measurement,  will be displayed. We highly recommend you consult your physician.

## MOVEMENT / READING ERROR

If **Err 1** is displayed, please take 2-3 minutes and then refasten the cuff before taking another measurement.

Refer again to **Preparing to Take Blood Pressure Measurement** section and **Positioning the Arm Cuff** section if required.

## VOICE ANNOUNCEMENT

If audio function is enabled, the following will be announced:

- Measurement results (Systolic, Diastolic, Heart rate) ;
- WHO classification;
- Low battery;
- Error detection
-  displays when audio function is enable

## VIEW RECORDS

This monitor can store up to 100 records per user. The most recent record (1) is shown first. Each new measurement is assigned to the first (1) entry in the user record selected. All other records are pushed back one digit (e.g., 2 becomes 3 and so on), and the last record (100) is dropped from the list.

For average morning / evening blood pressure readings to be displayed, blood pressure has to measure for three consecutive days at least twice in each time period (morning / evening).

For ease of reference, the date and time when the measurement was taken is shown with each record. Memory is kept even when the batteries become exhausted and subsequently replaced.

## RECALLING MEASUREMENTS STORED IN MEMORY

Ensure the monitor is off.

Select user:

- Press ▲ to view a history of User 1's records
- Press ▼ to view a history of User 2's records

View more records for the same user:

- Press the corresponding user button repeatedly ( ▲ or ▼ ) to see additional records (average morning / evening if available or next record).

**TIP** Press and hold corresponding user button ( ▲ or ▼ ) to quickly toggle through all the entries.

To alternate between different user's records:

Press ▲ or ▼.

## CLEAR RECORDS

When memory is full,  flashes.


Apart from deleting an entire measurement history, there is also an option to delete the latest measurement record. This is useful if measurements have not been accurately recorded and need to be recorded again.

## DELETING A SINGLE RECORD

1. When monitor is off, press ▲ or ▼ to select user.
2. Press the corresponding user button repeatedly ( ▲ or ▼ ) to navigate to the record to be deleted.

**NOTE** The most recent record (1) is shown first.


3. Press and hold . **dEL ONE** displays and flashes.
  - Press  to exit to the deletion mode.
  - Press and hold  again to confirm deletion.


4. **doNE dEL** indicates the record has been deleted.

## DELETING ALL RECORDS

1. When monitor is off, press ▲ or ▼ to select user.
2. Press and hold . **dEL ONE** will be displayed.
3. Press ▲ or ▼ to toggle to **dEL ALL**.
4. Press and hold again to confirm deletion.

This screen indicates all records have been successfully deleted and there are no records left for that user.

5. Press ▲ or ▼ to navigate to the other user's records.

## APPLICATION

You can find the instructions in the application or contact your local retailer of the product.

## CONNECT WITH YOUR MOBILE DEVICE


To help user keep track of their results, the measurements recorded by the monitor can be easily uploaded and reviewed through an application by Bluetooth Low Energy. Download the application "**Smart Blood Pressure Monitor**" in Apple Store through your mobile device.

### Turn on Bluetooth:

Please ensure to turn on the Bluetooth on your mobile device before using the application.

## START DATA TRANSMISSION

You need to build up User's Profile before starting the transmission.

### Edit user's profile:

1. Tap + to add up a user.
2. In **New User** screen, tap on the profile picture for uploading a photo; tap on each row to input **Gender**, **Date of birth**, **Weight** and **Height**.
3. Tap **Done** when you finish.

## NOTE


- If you want to edit each profile, you can also tap to modify the information.
- To remove a user, tap **Edit** in **User** screen, tap the icon of the user you want to remove. Delete displays and it means the deletion is successful.

## Start data transmission:

For added convenience, the unit is programmed to send data automatically:

- At the same time as taking measurement
- After taking measurement and unit detects there is data which was not sent successfully.

## Read a measurement on the app:

1. In User screen, select a user.
2. The screen simultaneously displays the measuring value exactly the same as your blood pressure monitor. **Connecting...** and  icon appears and flashes on the screen.
3. Once the connection is done, the  icon stops flashing and **BLE Done** displays on the screen of the monitor. You can start measurement and the record is automatically saved to the application.
4. Swipe the screen to right; you can view your measurement in graph view. Rotate the screen to be horizontal, and then you can review the data history in bar view.

### NOTE

- If the  icon disappears that means the connection is not successful.
- If there is an error during the Bluetooth connection, animated **BLE Err** displays on the top of the monitor screen.

## Manually send data:

1. Select a desired data history from user 1  or user 2 .

2. Press and hold  and  at the same time.

Successful transmission

Unsuccessful transmission

3. The screen display will turn off once synchronization attempt (whether successful or not) has been completed in a minute. Synchronizing status indicator light (located above the top right corner of the screen) varies as shown below:

Description	Meaning
Flashing green	Communicating
Continuous red	No connectivity

**IMPORTANT** During transmission, pressing does not have any function. Try to press the button in a minute.

## REVIEW THE RECORDS

You can view the records from the application by tapping .

1. You can view the records by **List View** or **Calendar View**.
2. Once you view the exact data history, you will view your measurement in graph view. Rotate the screen to be horizontal, and then you can review the data history in bar view.

**Delete a record:**

1. Swipe to left on a record,  displays on the record.
2. Tap on  to confirm the deletion.

**UNIQUE DEVICE ID**

Each unit has a unique device ID which can be found on the device. You can find it on the battery compartment.

**TROUBLESHOOTING**

This section includes a list of error messages and frequently asked questions for problems you may encounter with your blood pressure monitor. If the product is not operating as you think it should, check here before arranging for servicing.

ERROR CODE	Message	Reason/Suggestion
Err	There is an error during synchronization	Repeat the synchronization procedure
Err 1	The upper arm cuff is not appropriately secured.	Refasten the cuff and then take a measurement again
Err 4	The monitor detects motion during measurement.	Do not move or talk during measurement. Please try again.
Err 5	Pressure is over 300 mmhg	Relax for a moment and then take a measurement again

EE (3-20)	Internal software/hardware error	Please remove the batteries and then re-install them for reset. If the error code still displays on the screen, please contact our Customer Service for a solution.
BLE Err (animated and displayed on the top of screen)	There is an error during the connection.	
BLE Done	The connection is successfully done.	

frequently asked questions for problems you may encounter with your blood pressure monitor. If the product is not operating as you think it should, check here before arranging for servicing.

**SAFETY AND CARE INSTRUCTIONS****SAFETY PRECAUTIONS**

Please observe the following safety precautions when setting up and using your blood pressure monitor.

- This device is intended for adult use only.
- This device is intended for non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the upper arm or for functions other than obtaining a blood pressure measurement.
- Do not confuse self-monitoring with self-diagnosis. This unit allows you to monitor your blood pressure. Do not begin or end medical treatment based solely on the measurements of this device. Consult a physician for treatment advice.

- If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure. Never change a prescribed medication without consulting your physician.
- This unit is not suitable for continuous monitoring during medical emergencies or operations.
- If the cuff pressure exceeds 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when pressures exceeds 40 kPa (300mmHg), detach the cuff from the upper arm and press the button to stop inflation.
- To avoid measurement errors, carefully read this manual before using the product.
- The continuous cuff pressure due to connection tubing kinking may cause blood flow interference and harmful injury to the patient.
- Frequent measurements can cause injury to the patient due to blood flow interference.
- Do not apply a cuff over a wound, as this can cause further injury.
- The application of the cuff on any limb where intravascular access or therapy, or an arterio-venous (A-V) shunt may cause temporary interference to blood flow and could result in injury to the patient.
- The application of the cuff and its pressurization on the arm should be on the side of a mastectomy;
- Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring ME EQUIPMENT on the same limb.
- You need to check that operation of the blood pressure monitor does not result in prolonged impairment of patient blood circulation
- Do not use this monitor close to any magnetic resonance imaging (MRI) equipment. If MRI is to be performed, remove cuff from the arm.

## CARING FOR YOUR BLOOD PRESSURE MONITOR

To ensure you receive the maximum benefit from using this product, please observe the following care guidelines.

- When not in use, store the unit in a dry place away from direct sunlight.
- Do not immerse the unit in water. If it comes in contact with water, dry it immediately with a soft lint-free cloth.
- Carefully remove dirt from the unit or cuff with a soft, slightly moistened cloth. Do not use abrasive or corrosive cleaning detergent, such as petrol, thinners or similar solvents, as these may cause damage.
- Do not wash the arm cuff or immerse it in water. Position the arm cuff in ventilated and dry place, not under direct sunlight.
- Carefully remove spots on the cuff with a damp cloth and soapsuds.
- Remove the batteries whenever you are planning to store the unit for a long period of time.
- When replacing batteries, use new batteries as specified in this user manual. Do not mix new and old batteries.
- Do not place objects such as stickers on the upper arm cuff or unit, as these may impair the measurement.
- Do not subject the unit to excessive force, shock, dust, temperature changes, or humidity. Such treatment may result in malfunction, a shorter electronic life span, damaged batteries, or distorted parts.
- Do not tamper with the internal components. Doing so will terminate the product warranty and may cause damage.
- The unit contains no user- serviceable parts.
- If you no longer need to use this product, protect the environment by bringing it to your dealer or designated collection point for proper disposal.

## USEFUL REFERENCE

### ABOUT BLOOD PRESSURE

What is blood pressure?

Blood pressure is the force generated by the blood against the walls of arteries during cardiac contraction and relaxation (e.g., the pumping action of the heart).


What are systolic pressure and diastolic pressure? When ventricles contract and pump blood out of the heart, blood pressure reaches its maximum value. This highest pressure in the cycle is known as systolic pressure. When the heart relaxes between heartbeats, the lowest blood pressure is diastolic pressure.

Why measure your blood pressure?

Blood pressure measurement can highly reflect one's health condition. High blood pressure is potentially linked to serious illnesses such as stroke, heart disease and kidney failure. Since there is no symptom most of the time, many hypertensive people do not realize they are at risk until their health is seriously threatened.

What is the standard blood pressure classification?

The figure illustrates the blood pressure classification made by World Health Organization (WHO) and International Society of Hypertension (ISH) in 1999.


Reference material: 1999 World Health Organization - International Society of Hypertension Guidelines for the management of hypertension, Journal of Hypertension, 1999, 17(2): 151-183. This is used as the basis for the coloured hypertension indication bar to the right of the screen.

#### NOTE

- Blood pressure is considered high when either the diastolic or systolic blood pressure value exceeds the normal range. When a patient's systolic and diastolic blood pressures fall into different categories, the higher category should apply.
- Only a physician can tell you your normal blood pressure range and the point at which you are at risk. Consult your physician to obtain these values. If the measurements taken with these products fall outside the range, consult your physician. Why does my blood pressure fluctuate throughout the day? Individual blood pressure varies greatly both on a daily and a seasonal

or temperature basis. These variations may be more pronounced in hypertensive patients. Normally the blood pressure rises while at work and is at its lowest during sleep. The graph below illustrates the variations over a single day with measurement taken every 5 minutes. The thick line represents sleep time. The rises in blood pressure at 4PM (A in the graph) and 12AM (B in the graph) correspond to an attack of pain and sexual intercourse (Beven, Honour & Stott, Clin. Sci. 36:329, 1969).

## HOW THE BLOOD PRESSURE MONITOR WORKS

This product uses the Oscillometric Measuring method to detect your blood pressure. Before every measurement, the unit establishes a “zero pressure” equivalent to the air pressure.

Then it starts inflating the arm cuff to 180mmHg or higher until it senses that it has blocked your blood in the artery.

After that, the deflation process starts, during which time the unit detects pressure oscillations generated by beat-to-beat

pulsatile, which is used to determine the systolic, mean and diastolic pressure, and also your pulse rate. Any motion during this period will result in an incorrect measurement.

After all readings are determined and displayed on the LCD, the measurement is finished and the arm cuff automatically deflates.

## OPERATION FREQUENCY

The device's operation frequency is between 2.40GHz to 2.4835GHz, this is a unlicensed frequency band under CFR47, Part 15 of FCC rules, there is 40 RF channels, consists of 37 data channels for communication between connected devices and 3 advertising channels for discovering devices, the communication use GFSK modulation mode with Frequency-hopping spread Spectrum technology, the Maximum peak output power is 99.25dBuV/m by FCC test report(Report reference No.: TRE1204008401- Section 4.2: Radiated Emission), the effective connection distance is up to ten meters indoor.

## SPECIFICATIONS

Application	
Measuring method	Oscillometric / non-invasive
Application	For adult use only
Measurement location	Upper arm
Maximum number of users	2 and guest
Memory	Maximum 100 per user

<b>Dimensions</b>	
Dimensions L x W x H	148 x 103 x 75 mm (5.8 x 4.1 x 3.0 in)
Weight	450g
Cuff circumference	22-32 cm
Transmission range	6 meters (19.7 feet)
<b>Measuring range</b>	
Pressure	30 to 280 mmHg (4.0 to 37.3 kPa)
Pulse	40 to 199 pulse/min
<b>Power</b>	
Power supply	4 x UM-3 (AA) 1.5V batteries 6V adaptor (separate purchase required)
Power save	Auto power off after 1 minute of no key press or after synchronisation attempt
<b>Operating environment</b>	
Operation	10°C to 40°C (50°F to 104°F)
Storage	-20°C to 60°C (-4°F to 140°F)
Humidity	10% to 83% relative humidity

## ELECTROMAGNETIC COMPATIBILITY DESCRIPTIONS

**TABLE1: GUIDANCE AND MANUFACTURER'S  
DECLARATION – ELECTROMAGNETIC EMISSION –  
FOR ALL EQUIPMENT AND SYSTEMS**

1	<b>Guidance and manufacturer's declaration – electromagnetic emission</b>		
2	The BPU3210S Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of BPU3210S Blood Pressure Monitor should assure that it is used in such an environment.		
3	Emissions test	Compliance	Electromagnetic environment - guidance
4	RF emissions CISPR 11	Group 1	The BPU3210S Blood Pressure Monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

5	RF emissions CISPR 11	Class B	The BPU321 Blood Pressure Monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
6	Harmonic emissions IEC 61000-3-2	Class A	
7	Voltage fluctuations / flicker emissions  IEC 61000-3-3	Complies	

**TABLE 2: GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY – FOR ALL EQUIPMENT AND SYSTEMS**

<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The BPU321OS Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the BPU321OS Blood Pressure Monitor should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance

Electrostatic discharge (ESD)	± 6 kV contact  ± 8 kV air	± 6 kV contact  ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
IEC 61000-4-2			
Electrostatic transient / burst	± 2 kV for power supply lines  ± 1 kV for input/output lines	± 2 kV for power supply lines  ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
IEC 61000-4-4			
Surge	± 1 kV differential mode  ± 2 kV common mode	± 1 kV differential mode  ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
IEC 61000-4-5			

Voltage dips, short interruptions and voltage variations on power supply input lines  IEC 61000-4-11	< 5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle	< 5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the BPU321OS Blood Pressure Monitor requires continued operation during power mains interruptions, it is recommended that the BPU321OS Blood Pressure Monitor be powered from an uninterruptible power supply or a battery.
	40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles	40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles	
	70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles	70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles	
	< 5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	< 5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	

Power frequency (50/60 Hz) magnetic field  IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<b>NOTE</b> $U_T$ is the a. c. mains voltage prior to application of the test level.			

**TABLE 3: GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY – FOR EQUIPMENT AND SYSTEM THAT ARE NOT LIFE-SUPPORTING**

**Guidance and manufacturer's declaration – electromagnetic immunity**

The BPU321OS Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the BPU321OS Blood Pressure Monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	Portable and mobile RF communications equipment should be used no closer to any part of the BPU3210S Blood Pressure Monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Recommended separation distance $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz

			<p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
--	--	--	--

**NOTE** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE** These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BPU321OS Blood Pressure Monitor is used exceeds the applicable RF compliance level above, the BPU321OS Blood Pressure Monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the BPU321OS Blood Pressure Monitor .

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

**TABLE 4: RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE EQUIPMENT OR SYSTEM - FOR EQUIPMENT AND SYSTEMS THAT ARE NOT LIFE-SUPPORTING**

**Recommended separation distances between portable and mobile RF communications equipment and the BPU321OS Blood Pressure Monitor**

The BPU321OS Blood Pressure Monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BPU321OS Blood Pressure Monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BPU321OS Blood Pressure Monitor as recommended below, according to the maximum output power of the communications equipment

Rated maximum output of transmitter  W	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3.5}{V_1}\right]\sqrt{P}$	$d = \left[\frac{3.5}{E_1}\right]\sqrt{P}$	$d = \left[\frac{7}{E_1}\right]\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## ABOUT OREGON SCIENTIFIC

Visit our website [www.oregonscientific.com](http://www.oregonscientific.com) to learn more about Oregon Scientific products.

For any enquiry, please contact our Customer Services at [info@oregonscientific.com](mailto:info@oregonscientific.com).

Oregon Scientific Global Distribution Limited reserves the right to interpret and construe any contents, terms and provisions in this user manual and to amend it, at its sole discretion, at any time without prior notice. To the extent that there is any inconsistency between the English version and any other language versions, the English version shall prevail.

## AUTHORIZED REPRESENTATIVE

IDT TECHNOLOGY LIMITED,  
Block C, 9/F Kaiser Estate Phase  
1, 41Man Yue Sheet, Hung Hom,  
Kowloon, Hong Kong

Oregon Scientific Italia SPA  
Centro Dir.le Colleoni, Palazzo  
Taurus 2, Viale Colleoni 3, 20864  
Agrate Brianza (MB) Italy.  
Tel: +39-039-656-181  
Fax: +39-039-643-3111

## EU-DECLARATION OF CONFORMITY

Hereby, Oregon Scientific, declares that Smart Talking Blood Pressure Monitor (model: BPU3210S) is in compliance with the essential requirements and other relevant provisions of Directive 1999/05/EC. A copy of the signed and dated Declaration of Conformity is available on request via our Oregon Scientific Customer Service.

### COUNTRIES RTTE APPROVED COMPLIED

All EU countries, Switzerland (CH)  
and Norway (N)

## SYMBOL DESCRIPTION(S)

Type BF Equipment

CE 0120

Protected by US. Patent no. 7,211,047

Use the blood pressure monitor indoors only.

## FCC STATEMENT

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.

**WARNING** Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**NOTE** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna. Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio / TV technician for help.

**WARNING** This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 centimeters between the radiator and your body.

## DECLARATION OF CONFORMITY

The following information is not to be used as contact for support or sales. Please call our customer service number (listed on our website at [www.oregonscientific.com](http://www.oregonscientific.com)), or on the warranty card for this product) for all inquiries instead.

### We

Name: Oregon Scientific, Inc.  
Address: 19861 SW 95th Ave. Tualatin,  
Oregon 97062 USA  
Telephone No.: 1-800-853-8883

### declare that the product

Product No.: BPU321OS  
Product Name: Ssmart Talking Blood Pressure Monitor  
Manufacturer: IDT Technology Limited  
Address: Block C, 9/F, Kaiser Estate,  
Phase 1, 41 Man Yue St.,  
Hung Hom, Kowloon, Hong Kong

is in conformity with Part 15 of the FCC Rules. Operation is subject to the following two conditions: 1) This device may not cause harmful interference. 2) This device must accept any interference received, including interference that may cause undesired operation.

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