
# Panasonic®

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## Operating Instructions Instrucciones de funcionamiento

Model No. EW3006  
Modelo No. EW3006

Before operating this device, please read these instructions completely and save this manual for future use.

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Panasonic Oscillometric Diagnostec™ Wrist Blood Pressure Monitor Model EW3006 is a device intended to measure systolic and diastolic blood pressure and pulse rate of an adult individual by using a pressurized cuff on the left wrist. The device is not intended for use on infants and children. The device is designed for home use only, not for ambulatory measurement (measurement recorded continuously during the day).

Specification of this device including pulse rate (30 - 160 pulse/min. +/- 5%) are listed in page 22.

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard, Manual, Electronic or automated sphygmomanometers.

If you suffer from disorder of heart rhythm, known as arrhythmia only use this blood pressure monitor in consultation with your doctor. In certain cases oscillometric measurement method can produce incorrect readings.

Flash warning system for hypertensive readings is based on blood pressure values classified in the paper: "JNC 7 Express; The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation and Treatment of High Blood Pressure; U.S. DEPARTMENT OF HEALTH AND HUMAN SERVICES; National Institute of Health; National Heart, Lung, and Blood Institute; National High Blood Pressure Education Program; NIH Publication No. 03-5233; May 2003." The display values are generally known, but not proven, to be an indicator of your blood pressure.

The EW3006 is not intended to be a diagnostic device. Contact your physician if prehypertensive or hypertensive values are indicated.

# Table of Contents

Introduction .....	3
Basics of Blood Pressure .....	3
Important Instructions Before Use.....	4
Precautions to Ensure Safe, Reliable Operation .....	6
Easily Check Your Blood Pressure Readings Against the JNC 7 Classification .....	7
Name and Function of Each Product Part.....	8
Preparations prior to taking measurements .....	9
When to Change Batteries .....	9
Positioning the Pressure Cuff .....	10
Correct Position for Blood Pressure Measurement .....	12
Taking measurements while sitting down .....	12
Measuring Your Blood Pressure .....	15
Flash warning system for hypertensive readings .....	16
Movement detection function .....	18
Irregular pulse function.....	19
Storing blood pressure readings .....	20
Calling up readings from memory .....	21
Storing the Monitor.....	22
Specifications .....	22
Care and Maintenance.....	23
Troubleshooting .....	24

Product and part specifications are subject to change without notice.

# Introduction

Thank you for purchasing the Panasonic Wrist Blood Pressure Monitor EW3006.

Measuring your own blood pressure is an important way of monitoring your health. High blood pressure (hypertension) is a major health problem which can be treated effectively once detected. Measuring your blood pressure between doctor visits on a regular basis in the comfort of your home, and keeping a record of the measurements, will help you monitor any significant changes in your blood pressure. Keeping an accurate record of your blood pressure will help your doctor diagnose and possibly prevent any health problems in the future.

## Basics of Blood Pressure

Your heart acts like a pump, sending blood surging through your blood vessels each time it contracts. Blood pressure is the pressure exerted by blood pumped from the heart on the walls of blood vessels. Systolic pressure is the pressure exerted when the heart contracts and pumps blood into the arteries. Diastolic pressure is the pressure exerted when the heart expands, or relaxes. When you or your doctor take your blood pressure, both your systolic and diastolic pressures are measured. If your blood pressure measurement is 120 over 80 (120/80), for example, your systolic pressure is 120 while your diastolic pressure is 80.

# Important Instructions Before Use

1. Do not confuse self-monitoring with self-diagnosis. Blood pressure measurements should only be interpreted by a health professional who is familiar with your medical history.
2. If you are taking medication, consult with your physician to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your physician.
3. For persons with irregular or unstable circulation resulting from diabetes, liver disease, arteriosclerosis or other medical conditions, there may be variations in blood pressure values measured at the wrist versus at the upper arm. Monitoring the trends in your blood pressure taken at either the arm or the wrist is nevertheless useful and important.
4. Blood pressure can vary based on many factors, including age, gender, weight and physical condition. In general, a person's blood pressure is lower during sleep and higher when he or she is active. Blood pressure can change easily in response to physiological changes. The setting in which a person's blood pressure is measured can also affect the results. Having one's blood pressure measured by a healthcare professional in a hospital or clinic can cause nervousness and may result in a temporarily elevated reading. Because blood pressure measurements taken in a clinical setting can vary considerably from those taken at home, a person's blood pressure should be measured not only occasionally in the doctor's office, but also on a regular basis at home. Also, if you find that your blood pressure is lower at home, this is not unusual. To accurately compare with your physician's reading, take your Panasonic blood pressure monitor to your doctor's office and compare readings in this setting.
5. People suffering from cardiac arrhythmia, vascular constriction, liver disorders or diabetes, people with cardiac pacemakers or a weak pulse, and women who are pregnant should consult their physician before measuring their blood pressure themselves. Different values may be obtained due to their condition.
6. Try to take your blood pressure measurements at the same time and under the same conditions every day.
  - The ideal time to measure your blood pressure (to obtain your so-called "base blood pressure") is in the morning just after waking up, before having breakfast and before any major activity or exercise. If this is not possible, however, try to take measurements at a specified time prior to breakfast, and before you have become active. You should relax for

about 5 minutes before taking the measurement.

- The following situations may cause substantial variations in blood pressure readings and should therefore be avoided at least 30 minutes prior to taking your blood pressure.

Blood pressure will be higher than usual:

- when you are excited or tense
- when you are taking a bath
- during exercising or soon after exercising
- when it is cold
- within one hour after eating
- after drinking coffee, tea or other beverages containing caffeine
- after smoking tobacco
- when your bladder is full
- when in a moving vehicle

Blood pressure will be lower than usual:

- after taking a bath
- after drinking alcohol

7. Measurements may be impaired if this unit is used near a television, microwave oven, X-ray equipment or other devices with strong electrical fields. To prevent such interference, use the unit at a sufficient distance from such devices or turn the devices off.
8. This unit is designed for use by adults. Consult with your physician before using this unit on a child. Do not use on infants or toddlers.
9. This unit is not suitable for continuous monitoring during medical emergencies or operations.
10. Do not use the unit for any purpose other than measuring blood pressure. Do not use the unit together with other devices.
11. Improper handling of batteries may result in battery rupture or in corrosion from battery leakage. Please observe the following to ensure proper use of batteries.
  - a. Be sure to turn off the power after use.
  - b. Do not mix different types or sizes of batteries.
  - c. Change all batteries at the same time. Do not mix old and new batteries.
  - d. Be sure to insert batteries with correct polarity, as instructed.
  - e. Remove batteries when they are worn out, and dispose of them properly according to all applicable environmental regulations.
  - f. Do not disassemble batteries or throw them into a fire.
  - g. Do not short-circuit batteries.
  - h. Do not attempt to recharge the batteries included with the unit.

## ■ Differences in Blood Pressure Values Measured at the Wrist Versus at the Upper Arm

For persons with peripheral circulatory disorders\* resulting from diabetes, liver or kidney disease, arteriosclerosis or hypertension, etc., there may be variations in blood pressure values measured at the wrist versus at the upper arm.

Therefore, always consult with a health professional rather than attempting to interpret blood pressure measurements yourself.


\* Persons with underlying medical conditions with bad circulation in their hands and feet.

## **Precautions to Ensure Safe, Reliable Operation**

1. Do not drop the unit. Protect it from sudden jars or shocks.
2. Do not insert foreign objects into any openings.
3. Do not attempt to disassemble the unit.
4. Do not crush the pressure cuff.
5. If the unit has been stored at temperatures below 0°C (32°F), leave it in a warm place for about 15 minutes before using it. Otherwise, the cuff may not inflate properly.
6. Do not store the unit in direct sunlight, high humidity or dust.

# Easily Check Your Blood Pressure Readings Against the JNC 7\* Classification

## Blood Pressure Categories

If the two blood pressure measurements (systolic, diastolic) fall into separate categories, your level is classified in the higher of the two categories. For example, a Stage 2 systolic reading, but a diastolic pressure reading in the normal range.

\* JNC 7: The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation and Treatment of High Blood Pressure  
National Institute of Health (NIH) Publication; No. 03-5233, May 2003


## IMPORTANT:

- Do not be alarmed by temporarily high or low readings because fluctuations in a person's blood pressure are not uncommon. If possible, measure and record your blood pressure at the same time every day, and consult your physician if you have questions or concerns.
- If abnormal variations in blood pressure are observed in measurement, please consult your physician.

# Name and Function of Each Product Part

Fig. 1


## Batteries

- Panasonic alkaline batteries (LR03/AAA/Micro) can be used for approximately 300 measurements (approximately 3 months, 3 times a day), at an ambient temperature of 23°C (73.4°F) at 170 mmHg for a wrist circumference of 17 cm (6-3/4").  
(The batteries included are not guaranteed for the full life of brand new alkaline batteries as specified.)
- Battery life may be shortened if the ambient temperature is low.
- The number of measurements that can be taken with a set of batteries will be drastically reduced if batteries other than alkaline batteries are used. Use only alkaline batteries.
- Please replace both batteries at the same time and ensure that they are of the same type and make.


# Preparations prior to taking measurements

Fig. 2-a

1. Remove the unit from its storage case.  
Inserting dry cell batteries.  
(1) Remove the battery case by pushing it in the direction of the arrow.  
(2) Insert the batteries, ensuring the correct polarity (Fig. 2-a).
2. Slide the battery case back into place on the main body until it clicks into place (Fig. 2-b).


Fig. 2-b

## When to Change Batteries

Change the batteries when the battery symbol (Fig. 3) appears on the display or when no data is displayed even when the On/Off button is pressed.

Fig. 3

- NOTES:
- Batteries are included only for testing purposes and with no guarantee of battery life.
  - Do not mix old and new batteries.
  - The batteries' power output is impaired at low temperatures.

**Caution:** Do not throw batteries into a fire. Doing so may result in an accident or explosion.

# Positioning the Pressure Cuff

## CAUTIONS:

- Do not use the unit for any purpose other than measuring blood pressure. Also, do not use the unit in conjunction with other devices. Improper use may cause a malfunction or personal injury.
- Keep the unit out of the reach of children.

1. Remove your watch and any jewelry from your left hand.

2. Wrap the cuff around the wrist and secure it with hook-and-loop fastener.

The unit can be used by persons with a wrist circumference of 12.5 ~ 22.0 cm (5" ~ 8-3/4") (Fig. 4-a).

Wrap the pressure cuff around the wrist in such a way that the positioning mark is aligned with the center of the wrist with your wrist bent outwards.

3. In order to secure the blood pressure cuff in the proper position on the wrist, place your right thumb on the shorter side of the cuff (Hook fastener) pressing firmly while attaching to it the longer side of the cuff (Loop fastener). The pressure cuff should be fastened snugly with the display facing outward from your inner wrist (palm side).

Make sure that the pressure cuff is fixed firmly against your wrist, while also assuring it is centered on your wrist. If the pressure cuff is too loose, this will result in a false blood pressure measurement.

Fig. 4-a


Fig. 4-b



Fig. 4-c

- To remove the unit from your wrist, pull up on the longer side of the cuff without twisting it.

Fig. 4-d


Fig. 4-e

**When taking blood pressure measurements using the right hand, position the pressure cuff as shown in the diagram.**

The unit can be used on either the left or the right wrist.

However, as it is said that there can be a difference of up to 10mmHg between the left and the right wrists, please ensure that you always use the same wrist.

# Correct Position for Blood Pressure Measurement

When measuring your blood pressure at your wrist, the pressure cuff must be at heart level to ensure consistent, accurate readings.

Always rest for about four to five minutes before taking measurements.


Press the On/Off and Start button to start a measurement. The cuff will automatically inflate and then deflate when the measurement is complete.

## Taking measurements while sitting down

Fig. 5-a


1. Rest your elbows on a table or other level surface.
2. Hold the arm so that the cuff is at heart level.
3. With the palm of the hand facing up, release any tenseness in your muscles.

When taking blood pressure measurements while seated, please place your arm on top of the storage case once you have attached the pressure cuff.

If measurements are taken with the wrist or hand in the following positions, accurate measurements may not be possible.

Fig. 5-b


Wrist is bent inwards.


Wrist is tense and fist is being made.

If the blood pressure monitor is positioned lower than the heart, the measured blood pressure value may be high.

If the blood pressure monitor is positioned higher than the heart, the measured blood pressure value may be low.

If the posture applies pressure to the stomach (leaning forward), accurate measurements may not be possible.

Do not touch the blood pressure monitor while taking measurements.

Do not move your body, your arm, or your fingers while taking measurements.

Do not talk while taking measurements.

**Do not use close to cellular phones or other devices which emit electromagnetic waves. Doing so may result in malfunction.**

**IMPORTANT:**

- Measure your blood pressure at the same time and under the same conditions every day if possible. The ideal time to measure your blood pressure is in the morning immediately after waking up, while you are still relatively calm. If this is not possible, try to take your blood pressure prior to breakfast, before you become active.
- Relax and take measurements in a comfortable position.
- Rest for at least 5 minutes before taking a measurement.
- Keep still and do not speak during measurement. Blood pressure can fluctuate with the slightest change, such as tensing your muscles or changing your posture.

# Measuring Your Blood Pressure

Fig. 6-a

Fig. 6-b

Fig. 6-c

Fig. 6-d



Fig. 6-e

Relax during measurements.

1. Place the storage case under your arm. (Please see correct position for measurement on pages 12-14.)
2. Press the On/Off and Start button (Fig. 6-a).
  - Everything on the display will be shown for approx. 2 seconds, apart from the battery symbol (Fig. 6-b).
  - The cuff will begin to inflate automatically.

Changing inflation speeds and operating sounds occur during inflation and are a normal part of the measurements that are being performed at this time.

- Once your pulse rate is detected, the ♥ mark will flash on the display (Fig. 6-c).
3. Blood pressure values together with pulse rate are displayed when the measurement is complete (Fig. 6-d).
    - The air in the pressure cuff is released automatically.
    - Displayed readings will flash (for approx. 6 seconds) when in the high blood pressure range.
    - The **M** mark will flash on the screen to indicate that readings may be recorded. (See page 20 for information on how to record blood pressure readings.)
    - The displayed pulse rate indicates the number of heartbeats per minute as calculated from the number of heartbeats recorded during the measurement.
  4. Press the On/Off and Start button to switch off the unit (Fig. 6-e).
    - If you forget to switch the unit off, it will do so automatically after approx. 5 minutes.



# Flash warning system for hypertensive readings

If the measured values fall within the hypertension range, the reading displayed on screen flashes on and off to alert user.

According to JNC 7 Classification, values consistently in excess of 140 mmHg (Sys.) and/or 90 mmHg (Dia.) are considered to constitute high blood pressure.

Fig. 7

- Blood pressure readings will flash for approx. 6 seconds when in the high blood pressure range (Fig. 7).
  - Systolic blood pressure: 140 mmHg and over
  - Diastolic blood pressure: 90 mmHg and over
  - (Only when measurement is complete)

## When an Error Occurs During Measurement

Fig. 8

An  $E$  will appear on the display to indicate that measurement was unsuccessful and should be performed again (Fig. 8).

Before performing another measurement, always press the On/Off and Start button to turn off the unit and make sure to allow a rest period of 4 to 5 minutes before starting again.

Please consult your doctor without delay if the blood pressure indicator repeatedly draws attention to elevated blood pressure values.

**CAUTION:**

- For persons with poor circulation resulting from diabetes, liver disease, arteriosclerosis, high blood pressure or other conditions, there may be significant differences in blood pressure values measured at the wrist versus at the upper arm.
- Do not confuse self-monitoring with self-diagnosis.
- NEVER change a prescribed medication without first consulting with your physician.

# Movement detection function

If movement is detected while measurement is being carried out and more pressure is applied by the pressure cuff, the  mark will be displayed.

- Please carry out blood pressure measurement again if the  mark flashes or is lit.
- If the hand or wrist to which the pressure cuff is attached is moved while measurement is being carried out (e.g. the wrist is moved suddenly), in some cases correct measurements may not be obtained.
- When movement has been detected during measurement, the  mark flashes when the blood pressure reading obtained during that measurement is displayed.  
(If this reading is saved, the  mark will flash whenever it is recalled.)

Movement Mark	Explanation
 <p data-bbox="106 773 194 797">Flashing</p>	<p data-bbox="265 627 576 652">Movement has been detected.</p> <ul style="list-style-type: none"> <li data-bbox="265 656 712 681">· The hand or wrist has been bent suddenly.</li> <li data-bbox="265 685 824 710">· The muscles around hand or wrist have been tensed.</li> <li data-bbox="265 714 681 739">· The hand or wrist has been moved, etc.</li> </ul> <p data-bbox="265 744 909 797">With correct posture, please take the measurement again after measurement has been completed (see pages 12-14).</p>
 <p data-bbox="137 1020 168 1045">Lit</p>	<p data-bbox="265 824 971 848">Movement has been detected that will significantly affect the reading.</p> <ul style="list-style-type: none"> <li data-bbox="265 853 702 877">· Significant movement has been detected.</li> <li data-bbox="265 882 702 907">· Movement has been detected repeatedly.</li> <li data-bbox="265 911 873 936">· Significant hand or wrist movement has been detected, etc.</li> </ul> <p data-bbox="265 940 894 994">Please take the measurement again (An <b>E</b> will appear on the display).</p> <p data-bbox="265 998 883 1023">* An <b>E</b> will also appear on the display in the following cases:</p> <ul style="list-style-type: none"> <li data-bbox="265 1027 909 1081">· When the pressure cuff hook-and-loop fastener has become loose, etc.</li> </ul> <p data-bbox="291 1086 702 1110">(Please see page 25 for further details.)</p>

\* This function is designed to help you to take measurements accurately (see page 23). Even if the Movement mark is not displayed, we recommend that measurements be taken 2~3 times in a comfortable position in order to obtain an accurate reading. Remember to rest 5 minutes in between readings. (see pages 12-14 for correct position).

## Irregular pulse function

If the pulse rate varies by more than  $\pm 25\%$  from the average while a measurement is being taken, when measurement is completed the  mark will be displayed.

- The  mark will lit after approx. 6 seconds of blinking when the measurement is completed.
- If the  mark is lit, please take the measurement again.
- In some cases, if the pulse rate fluctuates greatly while a measurement is being taken, it may be impossible to obtain an accurate reading.

Fig. 9-a



Fig. 9-b

- If you save a reading where an irregular pulse was recorded, the  mark will be displayed when you recall that reading (Fig. 9-b).

If the  mark is displayed an irregular pulse rate is detected, not necessarily due to atrial fibrillation. However if this mark continues to display, consult your physician.

\* This function is designed to act as a guide to help you to obtain accurate readings (Please see pages 12-14).

Even if the  mark is not displayed, we recommend that you take measurements 2-3 times, resting 5 minutes between measurements.

\* If the  mark is frequently displayed, please consult a health professional for your health condition.

**Don't try to interpret readings or attempt to treat any condition yourself.**

**Always follow the guidance of a health professional.**

# Storing blood pressure readings

Up to 90 readings each can be stored in the memory.

(NOTE: When replacing flat batteries, no data will be lost even when the batteries have been removed.)

Fig. 10-a

Fig. 10-b


Fig. 10-c

1. The **M** mark will flash after a measurement has been completed (Fig. 10-a).  
Press the Memory/Recall button.
2. Memory storage is completed.
  - Readings of up to 90 measurements can be stored in memory (Fig. 10-b). Readings over the memory capacity will be recorded over the oldest measurement in memory. The data numbers will also be changed.
  - Readings cannot be saved when a measurement has been unsuccessful (display of the error mark **E** for blood pressure readings). As nothing will be recorded the data number will not be changed.
  - \* Readings can be saved, however, when the **E** mark is displayed for the pulse rate.

# Calling up readings from memory

Readings can be recalled after being saved without having to switch the unit on. Readings will not be displayed if they have not been saved.

Fig. 11-a

Fig. 11-b

Fig. 11-c

Fig. 11-d

1. Press the Memory/Recall button (Fig. 11-a).
  - The average reading of all recorded measurements is displayed.
  - The **Avg.** mark will be displayed (Fig. 11-b).
  - When, for example, readings from 90 measurements have been saved the average of the readings from 90 measurements will be displayed.
  - If there is only 1 measurement recorded, **Avg.** mark will not be displayed.
2. Press the Memory/Recall button to display readings starting with from the most recent measurement.
  - Each time the Memory/Recall button is pressed readings are displayed in order of the most recent measurement (Fig. 11-c, 11-d).
3. Press the On/Off and Start button to switch off the unit. If you forget to switch off the unit it will do so automatically after approx. 30 seconds.

## Deleting all data stored in memory

- (1) Recall stored data using the Memory/Recall button.
- (2) Press the Memory/Recall button again and hold it (for approx. 3 seconds) until all readings are deleted.  
NOTE: Readings from individual measurements cannot be deleted.
- (3) If everything has been deleted, the **M** mark and **00** will be displayed.

# Storing the Monitor

1. Fold the pressure cuff as indicated. (Fig. 12-a)
2. Put the blood pressure monitor in the storage case. (Fig. 12-b)

Fig. 12-a



Fig. 12-b

## Specifications

<b>Power source</b>	DC 3V (Two AAA/LR03 size alkaline batteries)
<b>Display</b>	Digital LCD
<b>Method of measurement</b>	Oscillometric System
<b>Measuring range</b>	Blood pressure: 0 to 280 mmHg Pulse rate: 30 to 160 beats per minute
<b>Accuracy</b>	Blood pressure: Within $\pm 3$ mmHg Pulse rate: Within $\pm 5\%$
<b>Operation temperature/ humidity range</b>	10°C to 40°C (50°F to 104°F) 15% to 90% RH
<b>Storage temperature/ humidity range</b>	-20°C to 60°C (-4°F to 140°F) 10% to 95% RH
<b>Wrist circumference measurement range</b>	Approx. 12.5 to 22.0 cm (5" to 8-3/4")
<b>Dimensions</b>	7.8 × 5.4 × 2.6 cm (3-1/8" × 2-1/8" × 1")
<b>Weight</b>	Approx. 104 g (3.7 oz.) (not including batteries)
<b>Protection against electric shock</b>	Internally powered equipment, Type BF applied part

# Care and Maintenance

Do not attempt to disassemble, repair or modify the unit.

- Doing so may cause fire or cause the unit to malfunction. It may also lead to injury.

Do not fold the pressure cuff back against itself.

Always use the pressure cuff on your wrist and nowhere else.

- Either of these actions may cause the device to malfunction.

Do not apply excessive force to the unit or drop the unit.

- Doing so may cause damage.

Do not insert dust or foreign objects into the unit.

- Doing so may result in damage.

If the unit is stored at temperatures below the freezing point, do not use it immediately.

Leave it in a warm place for at least an hour before use.

- If the unit is not allowed to warm up, it may not pressurize.

Do not touch the unit when measurement is in progress.

- Doing so may cause incorrect readings.

If the unit becomes dirty, clean it with a soft cloth moistened with warm water or soapy water.

(Do not use alcohol, benzene, or paint thinner.)

- Use of such chemicals may result in cracking or discoloration.

When the unit will not be used for a long time (30 days or more), be sure to remove the batteries.

- Otherwise, the batteries may leak and damage the unit.


The pressure cuff cannot be washed.

For storage, avoid high temperatures, high humidity, and direct sunlight.

- Storage under such conditions may result in damage.

# Troubleshooting

Display	Status before error	Check points
<p data-bbox="70 409 267 471">The E blinks on the display</p> 	<p data-bbox="350 212 643 404">Pressurization to above 280 mmHg was performed. (Re-pressurization was performed several times.)</p>	<p data-bbox="692 212 964 503">Check whether measurements are being taken correctly and whether you are assuming the correct posture. Make sure the unit is on the palm side of the wrist. (See pages 10-14.)</p>
	<p data-bbox="350 470 643 564">Pressure decreased suddenly and E appeared in the display.</p>	
	<p data-bbox="350 627 627 685">♥ mark only blinks a few times, or not at all.</p>	<p data-bbox="692 627 964 856">Check whether the pressure cuff is attached to your wrist correctly. Make sure the unit is on the palm side of the wrist. (See page 10.)</p>
	<p data-bbox="350 787 607 812">Cuff does not inflate.</p>	
<p data-bbox="350 946 576 1004">Your hand or wrist moved.</p>	<p data-bbox="692 946 954 1041">Has the unit detected movement? (See page 18.)</p>	

<b>Symptom</b>	<b>Possible cause</b>
The systolic pressure value or diastolic pressure value is high.	<ul style="list-style-type: none"> <li>• The hand with the cuff was held too low and not at heart level. (See pages 12-14.)</li> <li>• The cuff was not attached properly. (See pages 10-11.)</li> <li>• You moved your hand or body or you spoke during measurement. (See pages 12-14.)</li> </ul>
The systolic pressure value or diastolic pressure value is low.	<ul style="list-style-type: none"> <li>• The hand with the cuff was held too high and not at heart level. (See pages 12-14.)</li> <li>• You moved your hand or body or you spoke during measurement. (See pages 12-14.)</li> </ul>
The measured values differ significantly.	<ul style="list-style-type: none"> <li>• Your posture was not the same for all measurements, or measurements were taken without sufficient time in between. Always wait approximately 5 minutes before taking another reading, and above all, relax.</li> </ul>
The measured value differs from that measured by the doctor. The measured value is different for each measurement.	<ul style="list-style-type: none"> <li>• Blood pressure varies in response to minute changes in your mental state, such as your reaction to your doctor taking the measurement. One way to check this is to bring the unit to your doctor's office and directly compare the readings taken by your physician.</li> <li>• Relax for at least 5 minutes and then take another measurement.</li> </ul>
Varying sounds and pressurization speeds occur during inflation of the cuff.	<ul style="list-style-type: none"> <li>• Changes in operating sounds and pressurization speeds are normal and occur due to the fact that measurements are being made as the cuff is being inflated.</li> </ul>
The  mark flashes when readings are being taken.	<ul style="list-style-type: none"> <li>• There was a significant fluctuation in pulse when readings were being taken. (see page 19.)</li> </ul>

Symptom	Possible cause
The  mark flashes when readings are being taken.	<ul style="list-style-type: none"> <li>• Your hand or wrist moved. (see page 18.)</li> </ul>
The  mark flashes even though I did not move my hand or wrist.	<ul style="list-style-type: none"> <li>• In some cases this will be displayed when you have tensed the muscles in your arm. Relax the muscles in your arm and take the reading again. (see pages 12-14.)</li> </ul>
The measured value obtained is different from that measured at the upper arm.	<ul style="list-style-type: none"> <li>• For persons with certain circulation problems, there may be substantial differences in blood pressure values measured at the upper arm versus at the wrist.</li> </ul>

**If the unit still appears to provide unusual or erroneous readings, consult your physician. If the unit does not appear to be functioning properly, contact Panasonic at 1-800-338-0552.**

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Division of Panasonic Corporation of
**Panasonic Latin America, S.A.**