
EN English Introduction

blood pressure measurements determined with AV91 are equivalent to those obtained by a trained observer using a cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard, Electronic or Automated Sphygmomanometers. This unit is to be used by adult consumers in a home environment. Do not use this device on infants or neonates. AV91 is protected against manufacturing defects by an established International Warranty Program. For further information, you can contact the manufacturer, Rossmax International Ltd.

Attention: Consult the accompanying documents. Please read this manual carefully before use. For specific information on your own blood pressure, contact your physician. Please be sure to keep this manual.

Digital Technology
Easy to use manual inflation. This manual inflates blood pressure monitor offers accurate blood pressure and pulse measurements quickly and easily. Simply inflate the cuff to the desired pressure. The cuff will automatically deflate and show your readings on the easy to read display.

Preliminary Remarks
This Blood Pressure Monitor complies with the European regulations (CE mark "CE 01207". The quality of the device has been verified and conforms to the provisions of the EC council directive 93/42/EEC (Medical Device Directive). Annex I essential requirements and applied harmonized standards.

EN 1600-1: 1995/A2: 2009 Non-invasive sphygmomanometers - Part 1 - General requirements
EN 1600-2: 1997/A2: 2009 Non-invasive sphygmomanometers - Part 3 - Supplementary requirements for electro-mechanical blood pressure measuring systems
EN 1600-4: 2004 Non-invasive sphygmomanometers - Part 4: Test Procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.

This blood pressure monitor was designed for long service time. To ensure accurate measurements, this monitor is recommended to be re-calibrated every two years.

Blood Pressure Standard
The National High Blood Pressure Education Program Coordinating Committee has developed a blood pressure standard, classifying blood pressure ranges into 4 stages. (Ref. The Seventh report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure-Complete Report JNC-7, 2003). This blood pressure classification are based on historical data, and may not be directly applicable to any particular patient. It is important that you consult with your physician regularly. Your point at which you feel your normal blood pressure range as well as the point at which you should be considered at risk for reliable monitoring and reference of blood pressure, keeping long- term records is recommended. Please download the blood pressure log at our website www.rossmax.com.

Display Explanations
E1 / Measurement Error: Make sure the T-Plug is securely connected to the air socket and measure again quietly. Wrap the cuff correctly around the proper blood pressure for users diagnosed with hypertension, return the device to your local distributor or service center.

E1 / Air Circuit Abnormality: Make sure the T-Plug is securely connected to the air socket on the side of the unit and measure again quietly. If the errors still occur, return the device to your local distributor or service center for help.

E2 / Exceeding 300 mmHg: Switch the unit off and measure again quietly. If the error keeps occurring, return the device to your local distributor or service center.

E3 / Data Error: Remove the batteries, wait for 60 seconds, and reload. If the error keeps occurring, return the device to your local distributor or service center.

Irregular Heartbeat (IHB) Detector
The Irregular Heartbeat (IHB) Detector allows those who have an irregular heartbeat to obtain accurate measurements alerting the user of the presence of an irregular heart beat during the measurement.

Note: It is strongly recommended that you consult your physician if the IHB icon appears often ().

Installing Batteries
1. Press the battery cover in the battery cover in the direction of the arrow to open the battery compartment.

2. Install or replace 4 "AA" sized batteries in the compartment according to the indications inside the compartment.
3. Replace the battery cover by clicking in the bottom hooks first, then push in the top end of the battery cover.
4. Replace the batteries in pairs. Remove batteries when unit is not in use for extended periods of time.

You need to replace the batteries when:
1. Low battery icon appears on display.
2. The ON/OFF key is pressed and nothing appears on display.

Caution:
1. Batteries are hazardous waste. Do not dispose them together with the household garbage.
2. There are no user serviceable parts inside. Batteries or damage from old batteries are not covered by warranty.
3. Use exclusively brand batteries. Always replace with new batteries together. Use batteries of the same brand and same type.

Applying the Cuff
1. Unwrap the arm cuff, leaving the end of the cuff through the D-ring of the cuff.

2. Put your left arm through the cuff loop. The color strip indication should be positioned closer to you with the tube pointing in the direction of your arm (Fig. ②). Turn your left palm upward and place the edge of the arm cuff at approximately 15 to 25 cm above the inner side of the elbow joint (Fig. ②). Tighten the cuff by pulling the end of the cuff.

3. Center the tube over the middle of the arm. Press the hook and loop material together securely. Allow room for 2 fingers to fit between the cuff and your arm. Position the artery mark (②) over the main artery (on the inside of your arm) (Fig. ③, ④). Note: Locate the main artery by pressing with 3 fingers approximately 2 cm above the bend of your elbow on the inside of your left arm. Identify where the pulse can be felt the strongest. This is your main artery.

4. Plug in the cuff connecting tube into the unit (Fig. ⑤).
5. Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked (Fig. ⑥).
6. This cuff is suitable for your use if the arrow falls within the solid color line as shown on the right (Fig. ⑦). If the arrow falls outside the solid color line, you will need a cuff with other circumference. Contact your local dealer for additional size cuffs.

Measurement Procedures
Here are a few helpful tips to help you obtain more accurate readings:
1. Blood pressure changes with every heartbeat and is in constant fluctuation throughout the day.

2. Blood pressure recording can be affected by the position of the user, his or her physiological condition and other factors. For greatest accuracy, wait one hour after exercising, bathing, eating, drinking beverages with alcohol or caffeine, or smoking to measure blood pressure.
3. Before measurement, it is suggested that you sit quietly for at least 5 minutes as measurement taken during a relaxed state will have greater accuracy. You should not be physically tired or exhausted while taking a measurement.
4. Do not take measurements if you are under stress or tension.
5. During measurement, do not talk or move your arm or hand muscles.
6. Take your blood pressure at normal body temperature. If you are feel-

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