Gebrauchsanweisung
Instructions for Use
Mode d'emploi
Istruzioni per l'uso

(DE-2 - DE-19)

(EN-20 - EN-37)

(FR-38 - FR-55)

(IT-56 - IT-73)

visomat[®]
medic

www.visomat.de • Besser Messen.

UEBE
Germany est.1890

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Table of contents

A	Intended use	21
B	Safety instructions	
1.	Important patient instructions	22
2.	Important technical details	23
C	Operating the unit	
1.	Unit description	24
2.	Important instructions for use	25
3.	Attaching the cuff	26
4.	Fitting the stethoscope	27
5.	Inflating the cuff	27
6.	Measuring blood pressure	28
D	What you should know about blood pressure	
1.	Systolic and diastolic blood pressure values	29
2.	Reasons for measuring different values	29
3.	Why you should measure blood pressure regularly	29
4.	What are normal values for blood pressure?	29
E	Technical information	
1.	Measuring errors and their causes	30
2.	Customer service	31
3.	Technical data	32
4.	Original spare parts and accessories	33
5.	Metrological inspection	34
F	General provisions	
1.	Explanation of symbols	35
G	Maintaining the unit	36
H	Warranty	37

Intended use

A

Thank you for choosing the visomat® medic upper-arm blood pressure monitor (also referred to in the following as the unit).

It's determined for measuring blood pressure and pulse at the upper arm, those of people. This unit is recommended for patients with unstable blood pressure for domestic self administration blood pressure control and treatment support. It is suitable for patients with cardiac arrhythmia.

The unit is designed for auscultatory measurement of the systolic and diastolic blood pressure on the upper arm.

The process of inflating the cuff using the pump ball cuts off the blood flow in the upper-arm artery. During the deflation process thumping sounds can be heard through the stethoscope when the blood starts flowing in the artery again. Systole is determined when the noises begin, and diastole is determined with the last audible sound. The cuff can be fully deflated when both values have been determined.

These instructions are intended to help the user to use the blood pressure monitor and stethoscope safely and efficiently. They must be kept with the product and forwarded, if applicable.

The unit must be used in accordance with the procedures contained in these instructions and must not be used for other purposes.

Please read these instructions carefully before using the unit.

1. Important patient instructions

- Blood-pressure measurements on children and pregnant women require specialist knowledge! Please consult your doctor if you wish to measure the blood pressure of a child or a pregnant woman.
- The device does not have any pressure limitation at 150 mmHg and should therefore only be used for infants by medical experts.
- The supplied standard cuff is designed for an upper-arm circumference of 22-32 cm. Special cuffs are available for larger and smaller arm circumferences (Page EN-33).
- Do not under any circumstances place the cuff on or over any critical point, e.g. wound, aneurysm, etc. Risk of injuries! The supply via an intravascular access (infusion) or other medical monitoring devices might be interrupted.
- The airhoses pose a risk of strangulation. It should therefore be kept out of the reach of children at all times.
- If you have had a mastectomy, do not measure blood pressure on the side of the body where the breast has been removed together with the lymph nodes from the armpit.
- Measuring your own blood pressure does not constitute treatment. Do not modify of your own accord the dosage of drugs/medication prescribed by your doctor.
- Please refer to the chapter "Important instructions for use" (Page EN-25) before conducting any measurements yourself.


Safety instructions

B

2. Important technical details

- The unit must only be operated with original parts. The warranty will be invalidated if the unit is damaged by non-approved accessories!
- Never open or modify the unit or the cuff - this is a medical device and may only be opened by authorised specialists. If the unit has been opened, it must be subjected to a metrological inspection by an authorised institution.
- The cuff may be inflated on the arm only.
- Please comply with the specified operating conditions for measurement. See Technical data, Page EN-32.

1. Unit description

- | | |
|---------------------------------|---|
| 1 Air socket | 7 Check measurement arrow for arm circumference |
| 2 Stethoscope diaphragm | 8 Cuff |
| 3 Artery marking | 9 Ear clip |
| 4 Marking for arm circumference | 10 Pump ball |
| 5 Air connector | 11 Earpieces |
| 6 Air hose | 12 Pressure gauge with display dial |
| | 13 Deflating button |

Operating the unit

C

2. Important instructions for use

For a representative blood pressure value, please observe the following points both before and while taking the measurement:

- Refrain from drinking alcoholic or caffeinated beverages and smoking at least one hour before measuring.
- Please rest at least 5 minutes before measuring. Depending on the previous degree of stress/exertion, this can even require up to one hour.
- The cuff size has to be selected in accordance with the upper arm. Cuffs that are too small result in excessive blood pressure values. Measure the arm circumference and compare value to data on cuff!
- Expose the upper arm, the clothes must not impair the blood flow in or from the arm as this influences the blood pressure at the measurement point.
- Body posture must be relaxed:
 - For this purpose, sit at a table (if possible the height of a dining table, not a coffee table!).
 - Rest your back against the backrest of the chair.
 - Lay down your lower arm completely, with the palm of your hand facing upwards.
 - Place your feet on the floor and do not cross your legs.
- It is absolutely essential to sit still and remain quiet. Irregular movements, vibrations (shaking), talking and breathing heavily will affect and may even falsify the measurement.
- In order to rule out differences between left and right and to obtain comparable measurement results, it is important always to take blood pressure measurements on the same arm. Ask your doctor which arm he/she recommends you take your measurements on.
- Blood pressure is not a fixed value. It may fluctuate upwards or downwards by more than 20 mmHg on patients within a few minutes.

3. Attaching the cuff

- Expose the upper arm.
- Push the cuff onto the arm until the lower edge of the cuff is 2-3 cm above the crook of the arm (Fig. 1).
- When using the unit on the left arm the air hose runs to the unit through the middle of the crook of the arm so that the 4 cm long artery marking is located centrally above the pulse position (Fig. 2).
- When using the unit on the right arm the cuff must be rotated to the left until the artery marking is on the pulse position. The air hose then runs along the inner side of the upper arm.
- The cuff should be tightened to such an extent that 2 fingers can still be inserted between arm and cuff. Tip:
 - Angle the arm slightly
 - The upper arm muscle is tensioned slightly
 - This increases the circumference of the arm slightly
- Now pull the free cuff end tight and flip it back tightly above the Velcro fastener.
- Check whether the arrow on the cuff is within the white marking for the arm circumference.


Fig. 1

Fig. 2: Pulse position

Operating the unit

C

- Insert the air connector of the cuff into the air socket on the unit.
- Relax and lay your extended arm with the cuff on the dining table and do not move it during the measurement; do not talk. Make sure the palm of your hand is facing upwards.

4. Fitting the stethoscope

Fit the ear clip of the stethoscope so that the earpieces are firmly positioned in the ears. For a perfect stethoscope fit, turn the arms of the ear clip. It is recommended to turn the earpieces forwards slightly.

5. Inflating the cuff

- Close the air deflating valve by pressing the air deflating switch upwards.
- Using the pump ball, inflate the cuff until the needle on the dial indicates a value which is approx. 30 mmHg above the systolic blood-pressure value to be expected.
- If blood-pressure values are not known, inflate the cuff to approx. 250 mmHg. However, if you can clearly hear thumping sounds in the stethoscope, continue inflating until no more sounds can be heard.

6. Measuring blood pressure

- Using your index finger, gently press the air deflating switch and allow the air to escape slowly from the cuff.

- The pressure-gauge needle moves back slowly. You can vary the air deflation rate by adjusting the air deflating switch. (Recommendation of the WHO: 2-3 mmHg/second)

Now observe exactly the needle excursion on the dial and listen out for the noises which you can hear through the stethoscope.

- On hearing the first clearly thumping sound, read off the systolic blood-pressure value on the dial and make a note of the reading (Fig. 1).
- Listen out for rhythmic knocking in your stethoscope. When the sounds stop or at the moment of the last audible sound, read off the diastolic blood-pressure value on the dial and make a note of this reading too (Fig. 2).
- Press down the air deflating switch fully after taking the diastole reading. This initiates rapid deflation, the switch locking into place for this purpose. Remove the cuff after it has deflated fully.

Fig. 1


Fig. 2:

What you should know about blood pressure

D

1. Systolic and diastolic blood pressure values

The cardiovascular system has the important function of supplying all organs and tissues in the body with sufficient amounts of blood and of transporting metabolites. For this, the heart contracts and expands at a regular rate of about 60 to 80 times per minute. The pressure of the flowing blood on the artery walls caused by the heart contracting is termed systolic. The pressure in the ensuing relaxation phase, when the heart refills with blood, is termed diastolic. During daily measurement you determine both values.

2. Reasons for measuring different values

Our blood pressure responds to internal and external influences like a sensitive measuring instrument. It can be affected by even slight changes. This explains why often values measured with the doctor or pharmacist are higher than those measured at home in the environment you are used to. Changes in the weather, climate changes, or physical or psychological stress can have effects as well.

3. Why you should measure blood pressure regularly

Even the time of day has an influence on your blood pressure. During the day the values are generally higher than during the periods of rest at night. One-off and irregular measurements therefore say little about your actual blood pressure. A reliable assessment is possible only when measurements are taken regularly. Discuss the measurement values with your doctor.

4. What are normal values for blood pressure?

The World Health Organisation (WHO) established the following limits in mmHg (millimetre hydrargyrum / mercury) to classify the blood pressure values.

Normal value	Limit	High pressure
systolic pressure ≤ 119	systolic pressure 120–139	systolic pressure ≥ 140 and/or*
diastolic pressure ≤ 79	diastolic pressure 80–89	diastolic pressure ≥ 90

WHO 2003

* It is sufficient if one of the values is too high.

1. Measuring errors and their causes

If you fail to get any measurement results, please check for the following possible causes.

If necessary, repeat measurement after the circulation in the arm has been allowed to recover for a brief period; rest approx. 3-5 minutes for this purpose.

Failure encountered	Possible cause	Corrective action
Quiet pulse sounds	Stethoscope diaphragm not lying on the pulse position.	Reposition the stethoscope diaphragm; if necessary, have your specialised dealer or doctor show you how to do this.
	tube kinked.	Check tube.
	Earpieces clogged.	Clean earpieces.
Interference noise in the stethoscope	Movement of arm or body during measurement.	Do not move your arm, do not talk.
Cuff is not inflated	Deflating valve was not closed.	Reset the air deflating switch.
Incorrect measured values	Incorrect cuff size.	Measure your arm circumference in the middle of the upper arm and compare this with the details given on the cuff.
Unusual measured values	Incorrect cuff size. Movement or talking during measurement. Resting time not observed.	Check conditions and repeat measurement after 3 minutes.

2. Customer service

The unit may only be repaired by the manufacturer or by an expressly authorised organisation.

Please contact:

UEBE Medical GmbH
Zum Ottersberg 9
97877 Wertheim, Germany
Phone: +49 (0) 9342/924040
Fax: +49 (0) 9342/924080
E-mail: info@uebe.com
Internet: www.uebe.com

3. Technical data

Unit type:	Aneroid blood pressure monitor with pull-clip cuff and integrated stethoscope
Weight:	approx. 160 g without cuff
Measurement procedure:	Auscultatory measurement
Measurement and pressure-display range	0-300 mmHg
Error limits:	Blood pressure measurement: corresponds to EN ISO 81060-1 Pressure measurement: ± 3 mmHg
Serial number:	The unit bears a serial number SN which provides clear identification.
Cuff:	Available for private users from an arm circumference of 14-55 cm
	Available for professional users from an arm circumference of 8-55 cm
Operating conditions:	Ambient temperature 10 to 40 °C Relative humidity up to 85%, non-condensing
Storage and transportation conditions:	Ambient temperature -20 to +70 °C Relative humidity up to 85%, non-condensing
Deflation	Manual

4. Original spare parts and accessories

The following original spare parts and accessories are available from specialist dealers:

Available for professional users without Ring and stethoscope:

- Cuff type 1 (for arm circumference 8-15 cm)
Part no. 2700007 PZN 09786255
- Cuff type 2 (for arm circumference 14-21,5 cm)
Part no. 2700005 PZN 09786261
- Cuff type 3 (for arm circumference 18-25 cm)
Part no. 2700003 PZN 09786284
- Cuff type 4 (for arm circumference 26-33 cm)
Part no. 2700001 PZN 09786290
- Cuff type 5 (for arm circumference 34-43 cm)
Part no. 2700002 PZN 09786309
- Cuff type XXL (for arm circumference 43-55 cm)
Part no. 2700008 PZN 09786315

Available for private users with stethoscope:

- Cuff SM (for arm circumference 22-32 cm)
Part no. 2703001 PZN 05861184
- Cuff SL (for arm circumference 32-42 cm)
Part no. 2703002 PZN 05861209
- Cuff SS (for arm circumference 14-21 cm)
Part no. 2703005 PZN 05861244
- Cuff XXL (for arm circumference 43-55 cm)
Part no. 2703008 PZN 09786249

For further spare parts, please contact Customer Service (see Page EN-31).

Subject to technical modifications.

Manufacturer:

UEBE Medical GmbH
Zum Ottersberg 9
97877 Wertheim, Germany

5. Metrological inspection

(previously calibration)

Basically, a metrological inspection is recommended at intervals of 2 years. However, professional users in Germany have to comply with the aforementioned according to the "Regulation for Operators of Medical Devices". This can be implemented either by UEBE Medical GmbH, an authority responsible for metrology, or authorised maintenance services. For this, please observe your national provisions. Upon request, responsible authorities or authorised maintenance services receive a "Test instruction for metrological inspection" from the manufacturer. Important: No modifications, e.g. opening the unit, may be made to this unit without the manufacturer's prior permission.

General provisions

F


1. Explanation of symbols

This product complies with the Council Directive 93/42/EC from 5 September 2007 regarding medical devices and bears the mark CE 0123 (TÜV SÜD Product Service GmbH).

Please observe instructions for use

Storage and transportation conditions
Ambient temperature -20 to +70 °C

Protect against moisture/humidity
Rel. air humidity up to 85%

Keep dry

Manufacturer

Maintaining the unit

- The unit contains sensitive parts and must be protected against strong variations in temperature, air humidity, dust and direct sunlight.
- Make sure when storing the unit that no heavy objects are left lying on the unit or on the cuff and that the air hose is not kinked.
- The unit is neither impact-resistant nor shock-proof. We recommend that you have the intactness and accuracy of the unit checked after heavier falls or impacts.
- The unit is not water-proof.
- If the unit becomes dirty during use, use a soft, dry cloth to clean the unit. Do not use benzene, thinners or other strong solvents.
- Do not scrub or machine-wash the cuff. If the cuff gets dirty during use, use a synthetic cleaner and softly rub the surface.
- If the unit has been opened, it must be subjected to a metrological inspection conducted by an authorised organisation.
- Make sure that liquid cannot get into the air hose. Dry thoroughly.

Warranty

H

Warranty furnished

The device has been manufactured and tested with great care. However, in the unlikely event of a defect being detected after delivery, we provide warranty in accordance with the following terms and conditions:

1. During the warranty period of 3 years from the date of purchase we reserve the right either to repair any such defect at our expense or to supply a perfect replacement unit. The cost of returning the unit to our factory shall be borne by the sender. UEBE shall refuse to accept return deliveries that have not been paid for by the sender.
2. Excluded from the warranty are parts subject to normal wear and tear as well as damage caused by non-compliance with the instructions for use, improper handling (e.g. unsuitable power sources, breakages, leaking batteries) and/or disassembly of the unit by the purchaser. Furthermore, no claims for damages against us are substantiated by the warranty.
3. Warranty claims can only be advanced in the warranty period and by presenting proof of purchase. In the event of a warranty claim, the unit must be sent to the following address together with the proof of purchase and a description of the complaint:

UEBE Medical GmbH
Service-Center
Zum Ottersberg 9
97877 Wertheim
Germany

The cost of returning the unit to our factory shall be borne by the sender. UEBE shall refuse to accept return deliveries that have not been paid for by the sender.

4. In the case of defectiveness of the goods, the contractual rights of the purchaser to claim against the seller in accordance with § 437 German Civil Code are not limited by the warranty.

Please note: In the event of a warranty claim it is essential to attach the proof of purchase.

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© Copyright 2013 UEBE Medical GmbH