Model: **Checky**

Automatic Blood Pressure Monitor
Otomatik Tansiyon Aleti

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EN Instruction Manual

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EN Safety Information

The below signs might be in the user manual, labeling or other component. They are the requirement of standard and using.

| | | | |
|--|---|--|--|
| | Symbol for "THE OPERATION GUIDE MUST BE READ" | | Symbol for "TYPE BF APPLIED PARTS" |
| | Symbol for "COMPLIES WITH MDD 93/42/EEC REQUIREMENTS" | | |
| | Symbol for "MANUFACTURER" | | Symbol for "ENVIRONMENT PROTECTION - Wast electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice" |
| | Symbol for "SERIAL NUMBER" | | |
| | Symbol for "DIRECT CURRENT" | | |
| | Symbol for "MANUFACTURE DATE" | | |

CAUTION

Please do read this user manual carefully and thoroughly before use.

This device is intended for adult use only.

This device is intended for non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the upper arm or for functions other than obtaining a blood pressure measurement.

Do not confuse self-monitoring with self-diagnosis. This unit allows you to monitor your blood pressure. Please start or end medical treatment basing solely on physician's treatment advice.

If you are taking medication, consult your physician to determine the most appropriate time for your measurement. Never change a prescribed medication without your physician's consent.

This unit is not suitable for continuous monitoring during medical emergencies or operations.

If the pressure of the cuff exceeds 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when its pressure exceeds 40 kPa (300 mmHg), detach the cuff from the upper arm and press the homologous button to stop inflation.

Do not use the monitor under the conditions of strong electromagnetic field (e.g. mobile) that radiates interference signal or electrical fast transient / burst signal.

The device is not AP/APG equipment. It is not suitable for use in the presence of a flammable anesthetic mixture with air (or oxygen, nitrous oxide).


Please keep the unit out of reach of infants or children, since inhalation or swallowing of small parts is dangerous or even fatal.

Please use ACCESSORIES and detachable parts specified / authorised by MANUFACTURER. Otherwise, it may cause damage to the unit or danger to the user / patient.

The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5: 2009 and ISO 10993-10:2010. It will not cause any potential allergic reaction or contact injury.

Please make sure the unit functions safely and it is in proper working conditions before use.

EN LCD Display Symbols


| SYMBOL | DESCRIPTION | EXPLANATION |
|-----------------------------------|------------------------------|--|
| SYS | Systolic blood pressure | High blood pressure |
| DIA | Diastolic blood pressure | Low blood pressure |
| Puls min | Pulse | beat/minute |
| ▼ | Deflating | Exhaust the air in the cuff |
| 88:88 | Current Time | Month:Day (Hour:Minute) |
| 88 | Memory | When this symbol lights up, the measuring result displayed is retrieved from memory. |
| kPa | Unit | Measurement Unit of Blood Pressure |
| mmHg | Unit | Measurement Unit of Blood Pressure |
| Lo + □ | Low Battery | Low battery and please replace the batteries. |
| ♥ | Irregular Heartbeat Detector | Please refer to Page 31 for more information |
| G3 G2 G1 Norm Optimal | Grade | Classification of the blood pressure value. |

EN Monitor Components

Component List:

1. Cuff
2. Air pipe
3. PCBA
4. Pump
5. Valve

EN List

1. Blood Pressure Monitor (CHECKY)

3. 4*AA batteries

2. Cuff (Type BF applied part)
(22~32cm)

4. User manual

EN The Choice of Power Supply

1. Battery powered mode: 6VDC 4*AA batteries


CAUTION

In order to achieve the best performance and protect your monitor, please use the authorized / specified battery.

EN

Installing and Replacing the Batteries

1. Open the battery door.
2. Insert the batteries according to the polarity indications.
3. Close the battery door.

Battery Life: Approx. 70 days

(Battery capacity: 600 mAH. If measured twice per day, each measurement takes 35s, measuring result display takes 60s. The current for measurement is 350 mA and that for records display is 35 mA, while the current when shutdown is 25 μ A.)

Replace the batteries under following circumstances:

- displays on the LCD.
- The LCD display dims.
- When powering on the monitor, the LCD doesn't light up.

CAUTION

- Remove batteries if the device is not likely to be used for some time.
- Worn batteries are harmful to the environment. Do not dispose with daily garbage.
- Remove the old batteries from the device following your local recycling guidelines

[EN] Measurement Principle

This product uses the Oscillometric Measuring Method to detect blood pressure. Before every measurement, the unit establishes a “zero point” equivalent to the atmospheric pressure. Then it starts inflating the cuff. Meanwhile, the unit detects pressure oscillation generated by beat-to-beat pulsatile, which is used to determine the systolic pressure and diastolic pressure as well as pulse rate.

The device also compares the longest and the shortest intervals of detected pulse wave to with the average value, and then calculates the standard deviation. The monitor will light up a warning symbol when the calculated standard deviation is larger than or equal to 15.

[EN] Setting Date, Time and Measurement Unit

Please proceed to time setting before your initial use so as to ensure each piece of record are labeled with a time stamp. (Year Range: 2000-2050; Hour Format: 24-hour)

1. When the monitor is OFF, press and hold “SET” button for 3 seconds to enter year setting.


2. Press “MEM” button to change the numeral. Each press will increase the numeral by one in a cycling manner.


3. Press “SET” button to confirm the [YEAR]. Then the monitor diverts to [MONTH] and [DAY] setting automatically.


-
4. Repeat step 2 and 3 to confirm [MONTH] and [DAY].

5. Repeat step 2 and 3 to confirm [HOUR] and [MINUTE].

6. Repeat step 2 and 3 to confirm the measurement unit.


3. After confirming the measurement unit, the LCD will display "dONE" and the monitor will shut off automatically.


EN Tie the Cuff

1. Tie the cuff on your upper arm, the position the tube off-center toward the inner side of arm in line with the little finger.

2. The cuff should be snug but not too tight. You should be able to insert one finger between the cuff and your arm.

3. Sit comfortably with your left arm resting on a flat surface.


Resting for 5 minutes before measuring. Wait at least 3 minutes between measurements. This allows your blood circulation to recover. For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.

EN Start Measurement

1. After correctly positioning the cuff, press START/STOP button to turn on the monitor, and it will complete the measurement process automatically.

LCD Full Screen

Adjust to zero

Inflating and measuring

Display and save the measuring result


2. Press START/STOP button to turn off the monitor. Otherwise, the monitor will shut off within 1 minute..

EN Recall the Records


1. Press “MEM” button to access the memory.

The monitor will display the average value of the last three measurements.

1. Then press the “MEM/UP” button, the monitor will display the last measurement record. Press “MEM/UP” button or “SET/DOWN” button to rotate the history records. “MEM/UP” to go forward; “SET/DOWN” to go backward.

The exact time and date of the measurement will display alternatively.


The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back one digit (e.g., 2 becomes 3, and so on), and the last record (60) is dropped from the list.

EN Delete the Records

When you did not obtain the accurate measurement, you can clear all the measuring results by following below steps.

1. When the monitor is off, press and hold both the “MEM” button for 3 seconds

2. When the LCD display blinking “dEL ALL”, press “SET” button to confirm clearing. The LCD will display “dEL dONE”, indicating that memory clearing is complete.

3. If you wish to give up clearing, press the START/STOP button to turn off the monitor.

4. When there is no memory in the monitor, if you press the “MEM” button to look up history, the LCD will display as pictured to the right.

EN Tips for Measurement

It can cause inaccuracy if the measurement is taken in the following circumstances.

Within 1 hour
after dinner or drinking

Within 20 minutes
after taking a bath

In a very cold environment

In an electromagnetic field

Immediate measurement
after tea, coffee, smoking

When talking or moving your fingers

When you want to discharge urine


When you are diagnosed Arrhythmia

EN Maintenance

It can cause inaccuracy if the measurement is taken in the following circumstances.

Put in a dry place and avoid the sunshine


Avoid shaking and collision.

Use the slightly damp cloth to remove the dirt

Avoid immersing it in the water.
Clean it with a dry cloth in case.

Avoid dusty environment and unstable temperature surrounding

Avoid washing the cuff

Please make sure the unit functions safely and it is in proper working conditions before use. Please follow the instructions for correct replacement of interchangeable or detachable parts specified by SERVICE PERSONNEL of MANUFACTURER as "replaceable". Disposal: Degraded sensors may result in inaccurate measurement while loosened electrodes may cause the monitor's failure to power on. The expected life of this unit is two years. Please dispose of ACCESSORIES, detachable parts, and


ME EQUIPMENT according to local guidelines.

Before Use - Pick out the whole unit. Use the soft cloth to remove the dirt and apply some alcohol to disinfect the cuff.

After Use - Use the soft cloth to wipe the unit and apply some alcohol to disinfect the cuff before putting the whole unit back into the storage bag. Please always disinfect the cuff before applying to another patient.

EN What are systolic pressure and diastolic pressure?

When ventricles contract and pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.


EN What is the standard blood pressure classification?

The blood pressure classification published by World Health Organization (WHO) and International Society of Hypertension (ISH) in 1999 is as follows:

CAUTION

Only a physician can tell your normal BP range. Please contact a physician if your measuring result falls out of the range. Kindly note that only a physician could tell whether your blood pressure value has reached a dangerous point.

| Level Blood Pressure (mm Hg) | Optimal | Norm | | Hypertension | | |
|---------------------------------------|---------|---------|---------|--------------|---------|------|
| | | G1 | G2 | G3 | | |
| SYS | <120 | 121-130 | 131-140 | 141-160 | 161-180 | ≥180 |
| DIA | <80 | 81-85 | 86-90 | 91-100 | 101-110 | ≥110 |

EN Irregular Heartbeat Detector

This Blood Pressure Monitor is equipped with an intelligent function of Irregular Heartbeat (IHB) Detector. During each measurement, this equipment records the heartbeat intervals and works out the standard deviation. If the calculated value is larger than or equal to 15, this equipment will light up the IHB symbol on the screen when displaying the measuring result.

CAUTION

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

EN Why my blood pressure is varies even in one day?

1. Individual blood pressure varies every in one day, it also affected by the way you tie your cuff and the your measurement position, so please take the measurement at the same condition.
2. The varies of the pressure is greater if the person take medicine.
3. Waiting at least 4-5 minutes for another measurement.

EN Why the blood pressure I get from the hospital is different from home?

The blood pressure is different even during 24 hour because of the weather, emotion, exercise etc, specially the “white coat” in hospital which makes the results are higher than the ones at home.

The attention need to pay when you measure you blood pressure at home:
If the cuff is tied properly.
If you feel anxious pressured.
You had better take deep breath 2-3 times before beginning.
Advice:adjust yourself for 4-5 minutes until you calm down.

EN Why the blood pressure I get from the hospital is different from home?

It is ok for both arms, but there will be some different results for different person, so suggest you measure the same arm every time.

This section includes a list of error messages and frequently asked questions for problems you may encounter with your blood pressure monitor. If the products not operating as you think it should, check here before arranging for servicing.

| PROBLEM | SYMPTOM | CHECK THIS | REMEDY |
|----------------------|---|---|---|
| No power | Display is dim or will not light up. | Batteries are exhausted. Batteries are inserted incorrectly. | Replace with new batteries Insert the batteries correctly |
| Low batteries |  Show on the display | Batteries are low. | Replace with new batteries |
| Error message | E 1 shows | The cuff is not secure. | Refasten the cuff and then measure again. |
| | E 2 shows | The cuff is very tight. | Refasten the cuff and then measure again. |
| | E 3 shows | The pressure of the cuff is excessive. | Relax for a moment and then measure again. |
| | E 10 or E 11 shows | The monitor detected motion while measuring. | Movement can affect the accuracy of measurement. Relax for a moment and then measure again. |
| | E 20 or E 21 shows | Measure incorrectly. | Relax for a moment and then measure again. |
| | EExx,shows on the display. | System error. | Retake the measurement.If the problem persists, contact the retailer or our customer service department for further assistance.Refer to the warranty for contact information and return instructions. |

| | |
|---|---|
| Power supply | Battery powered mode: 6VDC 4*AA batteries |
| Display mode | Digital LCD V.A. = 93mm x 61mm |
| Measurement mode | Oscillographic testing mode |
| Measurement range | Pressure: 0kpa-40kpa (0mmHg-300mmHg) pulse value:(40-199)beat/minute |
| Accuracy | Pressure: 5°C-40°C within±0.4kpa(3mmHg) 0°C-45°C(out of 5°C-40°C) within±0.7kpa(5mmHg) pulse value:±5% |
| Normal working condition | Temperature:5°C to 40°C Relative humidity ≤85% Atmospheric pressure: 86kPa to 106kPa |
| Storage & transportation condition | Temperature:-20°C-60°C Relative Humidity 10%-93% Atmospheric Pressure: 50-106 kPa |
| Measurement perimeter of the upper arm | About 22cm~32cm |
| Weight | Approx. 289g(Excluding the dry cells) |
| External dimensions | Approx.140*120*70mm |
| Attachment | 4*AA batteries,user manual |
| Mode of operation | Continuous operation |
| Degree of protection | Type BF applied part |
| IP Classification | IP22 |
| Software Version | V01 |

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Table 1 Guidance and manufacturer's declaration – electromagnetic emissions-for all EQUIPMENT and SYSTEMS

| Guidance and manufacturer's declaration – electromagnetic emission | | |
|--|----------------|--|
| The CHECKY is intended for use in the electromagnetic environment specified below. The customer or the user of the CHECKY should assure that it is used in such an environment | | |
| Emission test | Compliance | Electromagnetic environment - guidance |
| RF emissions CISPR 11 | Group 1 | The CHECKY uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emission CISPR 11 | Class B | The CHECKY is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Not applicable | |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Not applicable | |

Table 2 Guidance and manufacturer's declaration – electromagnetic immunity – for all ME EQUIPMENT and ME SYSTEMS

| Guidance and manufacturer's declaration – electromagnetic immunity | | | |
|--|---|---|--|
| The CHECKY is intended for use in the electromagnetic environment specified below The customer of the user of the CHECKY should assure that it is used in such an environment | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±6 kV contact ±8 kV air | ±6 kV contact ±8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/burst IEC 61000-4-4 | ±2 kV for power supply lines | ±2kV for power supply lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ±1 kV line(s) to line(s) | ±1 kV differential mode | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5% U_T (>95% dip in U_T) for 0.5 cycle | <5% U_T (>95% dip in U_T) for 0.5 cycle | Mains power quality should be that of a typical commercial or hospital environment. If the user of CHECKY requires continued operation during power mains interruptions, it is recommended that CHECKY be powered from an interruptible power supply or a battery. |
| | 40% U_T (60% dip in U_T) for 5 cycles | 40% U_T (60% dip in U_T) for 5 cycles | |
| | 70% U_T (30% dip in U_T) for 25 cycles | 70% U_T (30% dip in U_T) for 25 cycles | |
| | <5% U_T (>95% dip in U_T) for 5 sec | <5% U_T (>95% dip in U_T) for 5 sec | |
| Power frequency (50Hz) magnetic field IEC 61000-4-8 | 3A/m | 3A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| NOTE U_T is the a.c. mains voltage prior to application of the test level. | | | |

Table 4 Guidance and manufacturer's declaration – electromagnetic immunity – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

| Guidance and manufacturer's declaration – electromagnetic immunity | | | |
|--|--------------------------------|------------------|--|
| The CHECKY is intended for use in the electromagnetic environment specified below. The customer of the user of the CHECKY should assure that it is used in such an environment . | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
| Conducted RF IEC 61000-4-6 | 3 Vrms 150 kHz to 80 MHz | 3 Vrms | <p>Portable and mobile RF communications equipment should be used no closer to any part of the CHECKY, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d = 1.167$</p> |
| Radiated RF IEC 61000-4-3 | 3 V/m 80 MHz to 2.5 GHz | 3 V/m | <p>$d = 1.167 \quad 80 \text{ MHz to } 800 \text{ MHz}$</p> <p>$d = 2.333 \quad 800 \text{ MHz to } 2.5 \text{ GHz}$</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacture and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p> |

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CHECKY is used exceeds the applicable RF compliance level above, the CHECKY should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the CHECKY.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Table 6 Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM – For ME EQUIPMENT or ME SYSTEM that are not LIFE-SUPPORTING

| Recommended separation distances between portable and mobile RF communications equipment at the CHECKY | | | |
|---|---|----------------------------------|-----------------------------------|
| Rated maximum output power of transmitter (W) | Separation distance according to frequency of transmitter (m) | | |
| | 150 kHz to 80 MHz $d = 1.167$ | 80 MHz to 800 MHz $d = 1.167$ | 800 MHz to 2.5 GHz $d = 2.333$ |
| 0.01 | 0.167 | 0.167 | 0.233 |
| 0.1 | 0.369 | 0.369 | 0.738 |
| 1 | 1.167 | 1.167 | 2.333 |
| 10 | 3.690 | 3.690 | 7.338 |
| 100 | 11.67 | 11.67 | 23.33 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 4 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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