Model: Handy

Blood Pressure Monitor
Tansiyon Aleti

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EN Instruction Manual

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EN Safety information

The below signs might be in the user manual, labeling or other components. They are the requirement of standard and using.

	Symbol for "THE OPERATION GUIDE MUST BE READ"		
	Symbol for "COMPLIES WITH MDD 93/42/EEC REQUIREMENTS"		Symbol for "ENVIRONMENT PROTECTION - Wast electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
	Symbol for "MANUFACTURER"		
	Symbol for "SERIAL NUMBER"		Symbol for "DIRECT CURRENT"
	Symbol for "TYPE BF APPLIED PARTS"		Symbol for "MANUFACTURE DATE"

CAUTION

Please do read this user manual carefully and thoroughly before use.

This device is intended for adult use only.

This device is intended for non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the wrist or for functions other than obtaining a blood pressure measurement.

Do not confuse self-monitoring with self-diagnosis. This unit allows you to monitor your blood pressure. Please start or end medical treatment basing solely on physician's treatment advice.

If you are taking medication, consult your physician to determine the most appropriate time for your measurement. Never change a prescribed medication without your physician's consent.

This unit is not suitable for continuous monitoring during medical emergencies or operations.

If the pressure of the cuff exceeds 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when its pressure exceeds 40 kPa (300 mmHg), detach the cuff from the wrist and press the START/STOP button to stop inflation.

Do not use the monitor under the conditions of strong electromagnetic field (e.g. medical RF equipment) that radiates interference signal or electrical fast transient / burst signal.

The maximum temperature that the applied part can be achieved is 42.5 °C while the environmental temperature is 40 °C.

The device is not AP/APG equipment. It is not suitable for use in the presence of a flammable anesthetic mixture with air (or oxygen, nitrous oxide).

Please keep the unit out of reach of infants or children, since inhalation or swallowing of small parts is dangerous or even fatal.

Please use ACCESSORIES and detachable parts specified / authorised by MANUFACTURER. Otherwise, it may cause damage to the unit or danger to the user / patient.

The material s of the cuff have been tested and found to comply with requirement s of ISO 10993-5:2009 and ISO 10993-10:2010. It won't cause any potential allergic reaction or contact


Manufacturer will make the circuit diagrams, component list, etc. available on request.

EN LCD Display Signal


SYMBOL	DESCRIPTION	EXPLANATION
SYS	Systolic Blood Pressure	High blood pressure
DIA	Diastolic Blood Pressure	Low blood pressure
Pul min	Pulse	beat/minute
■+Lo	Low Battery	Low battery and please replace the batteries.
mmHg	Unit	Measurement unit of blood pressure
kPa	Unit	Measurement unit of blood pressure
MEM 88	Memory	When this symbol lights up, the measuring result displayed is retrieved from memory. For instructions, please go to P27
▲	Inflating	Charge the air into the cuff.
▼	Deflating	Exhaust the air in the cuff.
AM PM MD 8888	Time	Hour:Minute (Day:Month / Year)
Avg	Average	Average value of measuring results
 	Grade	Please refer to P31 for classification standard of blood pressure.
之心	Irregular heartbeat Detector	Please refer to P31 for more information

EN Monitor Components

Component List:

1. PCBA
2. Air Pipe
3. Pump
4. Valve
5. Cuff


EN About Device

1.Blood Pressure Monitor
(Handy)

2.Two AAA-size batteries
3.User Manual

EN Installing and Replacing the Batteries

1. Open the battery door.
2. Insert the batteries according to the polarity indications.
Always select the authorized / specified battery:
(Two LR03 AAA-size batteries).
3. Close the battery door.

Battery Life: Approx. 58 days

(Battery capacity: 600 mAH. If measured 3 times per day, each measurement takes 30s, with the measuring results displayed for 10 seconds and checked memory once per day, which takes 60s. The current for measurement is 350 mA, and that for records display is 15 mA while the current when shutdown is 25 uA.)

Replace the batteries under following circumstances:

- +Lo displays on the LCD
- The LCD display dims
- When powering on the monitor, the LCD doesn't light up.


CAUTION

- Remove batteries if the device is not likely to be used for some time.
- Worn batteries are harmful to the environment. Do not dispose with daily garbage.
- Remove the old battery from the device following your local recycling guidelines.

EN Setting Date, Time and Measurement Unit

Please proceed to time setting before your initial use so as to ensure each piece of record are labeled with a time stamp. (Year Range: 2010-2050)

1. When the monitor is OFF, press the “SET” button to enter the time mode, then press and hold “SET” button for 3 seconds to enter the [YEAR] setting.


2. Press “MEM” button to change the numeral. Each press will increase the numeral by one in a cycling manner.


3. Press “SET” button to confirm [YEAR]. Then the monitor diverts to the [MONTH] and [DAY] setting automatically.


4. Repeat step 2 and 3 to confirm [MONTH] and [DAY].

5. Repeat step 2 and 3 to confirm the [HOUR] and [MINUTE].


6. Repeat step 2 and 3 to confirm the measurement unit.

7. After confirming the measurement unit, the LCD will display all the settings you have done one by one and the monitor will shut off automatically.

EN Positioning the Cuff

1. Remove all accessories (watch, bracelet, etc) from your left wrist. If your physician has diagnosed you with poor circulation in your left wrist, use your right wrist.
 2. Roll or push up your sleeve to expose the skin.
 3. Apply the cuff to your left wrist with your palm facing up.
 4. Position the edge of the cuff about 1-1.5cm.
 5. Fasten the wrist cuff around your wrist, leaving no extra room between the cuff and your skin. If the cuff is too loose, the measurement will not be accurate.
-
- Resting for 5 minutes before measuring.
 - Wait at least 3 minutes between measurements. This allows your blood circulation to recover.
 - For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same wrist, or as directed by a physician.


EN Start Measurement


1. After correctly positioning the cuff, press START button to turn on the monitor, and it will complete the measurement process automatically.


Adjust to zero

Inflating and measuring

Display and save the measuring result


2. Press START/STOP button to turn off the monitor. Otherwise, the monitor will shut off within 1 minute


EN Recalling the Records

1. Press “MEM” button to access the memory. The monitor will display the average of the last 3 measurement records.

2. Then press “MEM/UP” button again, the monitor will display the last measurement record. Press “MEM/UP” button or “SET/DOWN” button to rotate the history records. “MEM/UP” to go forward; “SET/DOWN” to go backward.

The exact time and date of the measurement will display alternatively.


CAUTION

The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back one digit (e.g., 2 becomes 3, and so on), and the last record (60) is dropped from the list.

EN Deleting the Records

When you did not obtain the accurate measurement, you can clear all the measuring results by following below steps.

1. Press the “MEM” button to access the memory. Then press and hold both the “MEM” button for 3 seconds, the LCD will display “dEL ALL”.


2. If you would like to clear the memory, press the “SET” button. The monitor will display “dEL dOnE” in response, and then turn off automatically.

3. If you wish to give up clearing, press “START/STOP” to turn off the monitor.

4. When there is no memory in the monitor, if you press the “MEM” button to look up the history, the LCD will display as pictured.

EN Tips for Measurement

It can cause inaccuracy if the measurement is taken in the following circumstances.

Within 1 hour
after dinner or drinking

Immediate measurement
after tea, coffee, smoking

Within 20 minutes
after taking a bath

When talking or moving your fingers

In a very cold environment

When you want to discharge urine

EN Maintenance

To obtain the best performance, please follow below instructions.

Put in a dry place and avoid the sunshine

Avoid immersing it in the water. Clean it with a dry cloth in case.

Avoid shaking and collision.

Avoid dusty environment and unstable temperature surrounding

Use the slightly damp cloth to remove the dirt.

Avoid washing the cuff

Please make sure the unit functions safely and it is in proper working conditions before use.

Please follow the instructions for correct replacement of interchangeable or detachable parts specified by SERVICE PERSONNEL of MANUFACTURER as "replaceable".

Disposal: Degraded sensors may result in inaccurate measurement while loosened electrodes may cause the monitor's failure to power on. The expected life of this unit is two years. Please dispose of ACCESSORIES, detachable parts, and MEDICAL EQUIPMENT according to local guidelines.

EN Cleaning


Before Use - Pick out the whole unit. Use the soft cloth to remove the dirt and apply some alcohol to disinfect the cuff before tying the cuff.

After Use - Use the soft cloth to wipe the unit and apply some alcohol to disinfect the cuff before putting the whole unit back into the storage bag.

Please always disinfect the cuff before applying to another patient.

EN What are systolic pressure and diastolic pressure?

When ventricles contract and pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.


EN What is the standard blood pressure classification?

The blood pressure classification published by World Health Organization (WHO) and International Society of Hypertension (ISH) in 1999 is as follows:

CAUTION

Only a physician can tell your normal BP range. Please contact a physician if your measuring result falls out of the range. Kindly note that only a physician could tell whether your blood pressure value has reached a dangerous point.

Blood Pressure (mm Hg)\Level	Optimal	Normal	High-normal	Mild	Moderate	Severe
SYS	<120	120-129	130-139	140-159	160-179	≥180
DIA	<80	80-84	85-89	90-99	100-109	≥110

EN Irregular Heartbeat Detector

This Blood Pressure Monitor is equipped with an intelligent function of Irregular Heartbeat (IHB) Detector. During each measurement, this equipment records the heartbeat intervals and works out the standard deviation. If the calculated value is larger than or equal to 15, this equipment will light up the IHB symbol on the screen when displaying the measuring result.

CAUTION

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heart-beat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

[EN] Why my blood pressure is varies even in one day?

1. Individual blood pressure varies every in one day, it also affected by the way you tie your cuff and the your measurement position, so please take the measurement at the same condition.
2. The varies of the pressure is greater if the person take medicine.
3. Waiting at least 4-5 minutes for another measurement.

[EN] Why the blood pressure I get from the hospital is different from home?

The blood pressure is different even during 24 hour because of the weather, emotion, exercise etc, specially the “white coat” in hospital which makes the results are higher than the ones at home.

[EN] If the result is the same if measuring on the right wrist?

It is ok for both wrists, but there will be some different results for different person, so suggest you measure the same wrist every time.

The attention need to pay when
you measure you blood
pressure at home:

If the cuff is tied properly.

If the cuff is too tight or too loose.

If the cuff is tied on the wrist.

If you feel anxious pressured.

You had better take deep breath

2-3 times before beginning.

Advice:adjust yourself for 4-5
minutes until you calm down.

This section includes a list of error messages and frequently asked questions for problems you may encounter with your blood pressure monitor. If the products not operating as you think it should, check here before arranging for service.

PROBLEM	SYMPTOM	CHECK THIS	REMEDY
No power	Display is dim or will not light up.	Batteries are exhausted. Batteries are inserted incorrectly.	Replace with new batteries Insert the batteries correctly
	 +Lo Show on the display	Batteries are low.	Replace with new batteries
Error message	E 1 shows	The cuff is not secure	Refasten the cuff and then measure again.
	E 2 shows	The cuff is very tight	Refasten the cuff and then measure again.
	E 3 shows	The pressure of the cuff is excess.	Relax for a moment and then measure again.
	E 10 or E 11 shows	The monitor detected motion while measuring.	movement can affect the measurement. Relax for a moment and then measure again.
	E 20 shows E 21 shows	Measure incorrectly	Relax for a moment and then measure again.
	EExx,shows on the display.	A calibration error occurred.	Retake the measurement. If the problem persists, contact the retailer or our customer service department for further assistance. Refer to the warranty for contact information and return instructions.

Power supply	3V 2*AAA batteries
Display mode	Digital LCD V.A. = 36mm x 41mm
Measurement mode	Oscillographic testing mode
Measurement range	Pressure: 0kpa-40kpa (0mmHg-300mmHg) pulse value:(40-199)beat/minute
Accuracy	Pressure: 5°C-40°C within±0.4kpa(3mmHg) 0°C-45°C (out of 5°C -40°C) within±0.7kpa(5mmHg) pulse value:±5%
Normal working condition	Temperature: 5°C to 40°C RH: ≤85% Atmospheric pressure: 86kPa to 106kPa
Storage & transportation condition	Temperature:-20°C to 60°C RH: 10% to 93% Atmospheric pressure: 50kPa to 106kPa
Measurement perimeter of the wrist	About 13.5cm-19.5cm
Net Weight	Approx.105g(Excluding the dry cells)
External dimensions	Approx.80×65×22mm
Attachment	2*AAA batteries,user manual
Mode of operation	Continuous operation
Degree of protection	Type BF applied part
IP Classification	IP22
Software version	V01
Device classification	Internally Powered ME Equipment

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EN EMC Guidance

Table 1 Guidance and manufacturer's declaration – electromagnetic emissions-for all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration – electromagnetic emission		
The HANDY is intended for use in the electromagnetic environment specified below. The customer or the user of the HANDY should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The HANDY must emit electro-magnetic energy in order to perform its intended function. Nearby electric equipment may be affected.
RF emission CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Table 2 Guidance and manufacturer's declaration – electromagnetic immunity – for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration – electromagnetic immunity			
The HANDY is intended for use in the electromagnetic environment specified below. The customer or the user of the HANDY should assure that it is used in such an environment			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	N/A	
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	N/A	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle	N/A	
	40% U_T (60% dip in U_T) for 5 cycles	N/A	
	70% U_T (30% dip in U_T) for 25 cycles	N/A	
	<5% U_T (>95% dip in U_T) for 5 sec	N/A	
Power frequency (50Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Table 4 Guidance and manufacturer's declaration – electromagnetic immunity – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity			
The HANDY is intended for use in the electromagnetic environment specified below. The customer or user of HANDY should assure that it's used in such an environment .			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	<p>Portable and mobile RF communications equipment should be used no closer to any part of the HANDY, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.167 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.167 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.333 \sqrt{P} \text{ 800 MHz to 2.5 GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacture and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the HANDY is used exceeds the applicable RF compliance level above, the HANDY should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the HANDY.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Table 6 Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM – for ME EQUIPMENT or ME SYSTEM that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment at the HANDY			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.167 \sqrt{P}$	80 MHz to 800 MHz $d = 1.167 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.333 \sqrt{P}$
0.01	N/A	0.117	0.233
0.1	N/A	0.369	0.738
1	N/A	1.167	2.333
10	N/A	3.690	7.377
100	N/A	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 4 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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