


OKmeter®


*Link*

OKmeter® Link  
Blood Glucose  
Monitoring System

Seeking worldwide distributors  
Open for any contract manufacturing




## STEPS



1

Insert test strip to  
turn on meter.

2



Apply small  
blood sample.

3



Obtain results in 6  
seconds.

➤ **No Coding**  
➤ **Alternative  
Site  
Testing**

- Accurate results in 6 seconds
- Very small blood sample, only 0.7  $\mu$ L
- 180 test memory with date and time
- PC download



Manufacturer

## OK Biotech Co., Ltd.

4F-1, No. 83, Sec. 2, Gongdao 5th Rd., 300  
Hsinchu City, Taiwan

TEL: +886-3-5160258 FAX: +886-3-5160028

<http://www.okbiotech.com>

<http://www.glucose-meter-manufacturer.com>

[service@okbiotech.com](mailto:service@okbiotech.com)



OK9C-EN-110817

# SYSTEM ACCURACY (ISO 15197)

**Test system:** OKmeter  
**Reference system:** YSI2300  
**Sample:** Whole blood (capillary)

### Investigative Site:

Institut für Diabetes-Technologie GmbH an der  
Universität Ulm (IDT)  
Helmholtzstrasse 20, D-89081 Ulm / Germany


### System accuracy results for glucose concentration < 4.2 mmol/L (75 mg/dL)

Within ± 0.28 mmol/L (Within ± 5 mg/dl)	Within ± 0.56 mmol/L (Within ± 10 mg/dl)	Within ± 0.83 mmol/L (Within ± 15 mg/dl)
19 / 40 (48 %)	34 / 40 (85 %)	40 / 40 (100 %)

### System accuracy results for glucose concentration > = 4.2 mmol/L (75 mg/dL)

Within ± 5 %	Within ± 10 %	Within ± 15 %	Within ± 20 %
64 / 160 (40 %)	118 / 160 (74 %)	154 / 160 (96 %)	159 / 160 (99 %)

The OKmeter blood glucose system complies with the system accuracy requirements of the ISO 15197 Standard. 199 of 200 (99.5%) results meet the requirements.



# SPECIFICATIONS

<b>Measurement Range</b>	20~600 mg/dL (1.1~33.3 mmol/L)
<b>Specimen Type</b>	Capillary Whole Blood
<b>Alternative Site Testing</b>	YES (palm, forearm, upper arm, calf, and thigh)
<b>Calibration Against</b>	Plasma-Equivalent
<b>Test Time</b>	6 seconds
<b>Sample Volume</b>	0.7 µl
<b>Sensitivity in HCT Range</b>	20-60%
<b>Operating Conditions</b>	10~40°C (50~104°F) / 10~85% R.H. Altitude ≤ 3,402 m (11,161 ft) above the sea level
<b>Storage/Delivery Conditions</b>	4~40°C (39~104°F) / ≤ 10~85% R.H. (strip included) -20~50°C (-4~122°F) / ≤ 90% R.H. (strip excluded)
<b>Memory Capacity</b>	180 test results with date and time
<b>Average Calculation</b>	7, 14, 21 and 28 days
<b>Downloadable</b>	YES (RS-232 port)
<b>Battery Type</b>	3V Li-Battery (CR 2032) X1
<b>Battery Life</b>	1,000 times of measurements
<b>Auto OFF</b>	3 minutes without action
<b>Meter Dimensions</b>	70mm(L) X 57mm(W) X 18.6mm(H)
<b>Meter Weight</b>	Approximate 43.5g (with battery)
<b>Expiration for test strip</b>	Storage: 24 months In-Use: 180 days