

Digital Blood Pressure Monitor - Wrist Type

Model No: PH 168A

The Proton PH 168A Digital Blood Pressure Monitor is an oscillometric type blood pressure monitor which produces accurate and precise blood pressure measurements.

The PH 168A is capable of reading systolic and diastolic ranges generally unattainable by other electronic blood pressure monitors. This unique design also assists in reducing average measurement time resulting in increased user comfort.

The PH 168A has a large, easy to read numeric display which shows Systolic, Diastolic and Pulse simultaneously.



Features

- Intelligent inflation by means of fuzzy logic
- Fully automatic measuring procedure
- Oscillometric measurement of blood pressure and pulse
- Measures from 0-300mmHg
- Includes 48 sets of memory
- Power source: AAA x 2 pcs Alkaline Battery
- Unit Dimension: L 70 x W 72 x H 29 mm




WRIST TYPE PH 168A

Specifications

Method of Measurement	Oscillometric
Model Number	PH 168A
Range of Measurement	Pressure 0-300mmHg, pulse 40-199 beats/minute
Accuracy Calibration	Pressure +/- 3mmHg, pulse +/- 5%
Fuzzy Pressure	Automatic Re-inflation pump system
Inflation	Automatic inflation (Air Pump)
Deflation of Pressure	Automatic air release control valve
Exhaust	Automatic exhaust valve
Display	Liquid Crystal Digital display
Unit Dimension	L70 x W72 x H29mm
Unit Weight	Approx. 140g including batteries
Sets of Memory	48 set
Cuff Size	Wrist circumference approx. 135-195mm
Operating Temperature	+10 C to +40 C, 30~85% RH
Storage Temperature	-20 C to +50 C, 10-95% RH
Power Supply	2 x AAA, (3.0V) Alkaline
Battery Life	Approx. 250 times measurement
Accessories	Carrying case, instruction manual, 2 pcs AAA Alkaline battery
Warranty	3 years

Digital Blood Pressure Monitor - Wrist Type


Model No: **PH 168A**



WRIST TYPE PH 168A

A. Wrap the pressure cuff around the wrist

- The display of the unit should be placed on the palm side of the wrist
- The wrist should be bare




B. Fasten the pressure cuff snugly

- Do not pull strongly on the pressure cuff
- Do not make the pressure cuff too tight
- Fasten the cuff leaving a 2 finger distance from the wrist joint.

Measurement while sitting down:

- Place your elbow on a table or other object (such as our carry case)
- Use the armrest to position the arm so that the pressure cuff is at the same height as the heart
- Relax your hand with the palm facing up.
- Press the on/off button, the cuff will inflate and deflate to record the systolic, diastolic and the pulse. Remove the cuff and store the unit and the cuff in its bag safely.
- If the unit is not going to be used daily, remove the battery and keep separately to avoid any battery leakage which may spoil the electric contacts of the BP Monitor.



NOTE: Please read through the Instruction Manual for detailed guidelines on use, maintenance, storage and troubleshooting.