




# TD-3140B

## UPPER-ARM BLOOD PRESSURE MONITOR


**Mb<sup>®</sup>H**  
Mobile Health




### Smartphone App for cardiovascular evaluation

-  **HI** Fast Heart Beats  
Heart rate: > 110 bpm.


---

-  **LO** Slow Heart Beats  
Heart rate: < 50 bpm.


---

-  **IHB** Irregular Heart Beats  
Heart rate: premature, delayed contraction, skipped, pause between heartbeats, extra systole.

---

-  **AF** Atrial Fibrillation  
Heart rate: > 100 bpm, quivering or irregular heartbeat.

---

-  **Cardiovascular Evaluation**  
Evaluation: cardiovascular or age assessment.



**TaiDoc**  


[www.taidoc.com](http://www.taidoc.com)



### SPECIFICATIONS

Systolic Measuring Range	50 ~ 250 mmHg
Diastolic Measuring Range	30 ~ 180 mmHg
Accuracy of Pressure	±3mmHg or ±2% of reading
Heart Rate Measuring Range	40-199 beat / minute
Accuracy of Heart Rate	±4% of reading
Operating condition	Temperature 5 °C ~ 40 °C, Humidity 15% ~ 93%
Memory Capacity	200 sets
Backlight	LCD backlight
Time to Power saving	180 Seconds
Dimension	145 (L) x 64.5 (W) x 29.6 (H) mm
Weight	240 g (without battery)
Power source	4 x AAA
Communication	Bluetooth smart
Cuff Size	M Size: 24 ~ 35 cm Wide Range: 24 ~ 43 cm



www.taidoc.com

**TaiDoc Technology Corp.**

B1-7F., No. 127, Wugong 2nd Rd., Wugu Dist.,  
New Taipei City 24888, Taiwan  
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288