

**New for
September
2013**



Introducing GDS 100 Hand-held semi-automatic blood pressure measurement and recording



A built-in pump quickly inflates the BP cuff at the touch of a button



An electronically controlled valve regulates the deflation of the cuff



A scroll wheel allows for adjustment of the rate of deflation



An LCD provides a continuous digital display of cuff pressure



Upon hearing the Korotkoff sounds the operator can activate a convenient trigger mechanism to capture the systolic and diastolic pressure



Following the capture of diastolic pressure the cuff automatically deflates and the BP result is displayed in full on the LCD

Updated SP-800 Nonin and Nellcor compatible SpO2 available September

