

# K-jump®

## Wrist Watch Blood Pressure Monitor



KP-6155



KP-6165

### Specifications

Model	KP-6155	KP-6165
Measuring Range	20 ~ 300 mmHg 40 ~ 200 pulse/minute	
Resolution	1 mmHg	
Accuracy	±3 mmHg or 2% (pressure) ±5% of readings (pulse rate)	
Measuring Method	Oscillometric method	
Operation Key	Three Keys: Power/Start, Memory, Time Setting	
Display	Liquid Crystal Display	
Inflation	Micro rolling pump	
Exhaust	Electrical Solenoid Valve	
Memory	42-set, or 85-set	
Time	Displaying date, time and year	
Operation Condition	Temperature: 10°C ~ 40°C (50°F ~ 104°F) Relative Humidity: less than 85%	
Storage Condition	Temperature: -10°C ~ +50°C (+14°F ~ +122°F) Relative Humidity: less than 85%	
Power Source	3Vdc; "AAA" type battery x 2	
Battery Life	Around 250 measurements (alkaline batteries)	
Low Battery Indication	Replace the batteries if the low battery indicator appears	
Auto Off	Automatically after 150 seconds without operation	
Size	78 x 72 x 30 mm	78 x 72 x 30 mm
Weight	130g w/ batt.	
Standard	Complies with EN 1060-1, -3 and ANSI/AAMI SP10:2002	

Specifications are subject to change for improvement without notice.

### K-jump Health Co., Ltd.

No.56,Wu Kung 5th Rd.,New Taipei Industrial Park,

New Taipei City 24890,Taiwan


Tel: 886-2-22991378 Fax: 886-2-22991386

E-mail: service@kjump.com

http://www.kjump.com



CE 0197



HealthCare at Home  
we care for Your Health