TRUST MEDICAL INSTRUMENT

TD-8201 Handheld Pulse Oximeter

-

TD-8201



- Large and Clear LED Screen
 Water-Resistant
 Shock Resistance Proof
 - Bluetooth Transmission



Product Features:

- Fast, Accurate, Flexible
- Large and Clear LED Numerical Display

MTRUST

- Water-Resistant and Shock-Resistant for Optimal Protection in All Conditions
- Ergonomic Design Provides Simple, Single-hand Operation
- Rubber Grip Offers Secure Handling
- Bluetooth Data Transmission (Option)

TRUST TD-8201 Handheld Pulse Oximeter

Product Specifications:

Pulse Oximeter -SpO₂ Display Range: SpO₂ Resolution: SpO₂ Accuracy :

SpO₂ Method : Pulse Rate Measuring Range: Pulse Rate Accuracy: Mode of Operation:

LED Display: Status Indicate: External Output:

Alarm Range Setting: Pulse Volume Adjust: LED Light Adjust:

Power Source: Battery Usage Life: Dimension: Weight:

Storage/Transportation Condition: Operating Condition: Degree of Protection: Water-Resistance: Mechanical Impact-Resistance:

EMC Standard: Safety Standard: Performance Standard:



OXM-003 SpO2 Sensor (1m, clip)

0-100% 1% 80-100%, ±2%; 70-79%, ±3%; <70%, no definition


Dual Wavelength LED 30 - 250 bpm ±1 bpm or ±1% when HR < 100 bpm Continues

SpO₂: 3 Digits, HR: 3 Digits, Pulse-Strength Bar Probe Good, Battery Status TD-8201A Non-Bluetooth TD-8201B Bluetooth smart SpO₂ High/Low , Pulse Rate High/Low 2 Level or Off 2 Level

AA Alkaline Battery x 4 40 hours 140 (L) x 70 (W) x 32 (H)mm 240g (with batteries)

-20°C-60°C; below 95% RH, non-condensing 10°C-40°C; below 95% RH, non-condensing Type BF Against Water Splash, IPX2 Against Repeat Impact from 1 M Height

EN60601-1-2 EN60601-1 ISO 80601-2-61





TaiDoc Technology Corp.



B1-7F., No.127, Wugong 2nd Rd., Wugu Dist., New Taipei City 24888, Taiwan Tel : +886-2-6625-8188 Fax : +886-2-6625-0288